

Nov 18, 2021

# Citric Acid Water Restriction

Thomas R Clarke<sup>1</sup><sup>1</sup>University of Edinburgh

2

[dx.doi.org/10.17504/protocols.io.buw4nxgw](https://dx.doi.org/10.17504/protocols.io.buw4nxgw)Thomas Clarke  
University of Edinburgh

How to prepare and carry out citric acid water restriction for motivating mice to engage in behaviour.

DOI

[dx.doi.org/10.17504/protocols.io.buw4nxgw](https://dx.doi.org/10.17504/protocols.io.buw4nxgw)

Thomas R Clarke 2021. Citric Acid Water Restriction. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.buw4nxgw>



Mouse, Behaviour, Water Restriction, Citric Acid

\_\_\_\_\_ protocol ,

May 11, 2021

Nov 18, 2021

49852

Any mice dropping to below 80% of their pre-CA baseline weight must be returned to *ad lib* water (see Duguid Lab PPL for regulations)

## Baseline Weights 3d

- 1 Handle animals and weigh for 3-5 days to establish a mean baseline weight prior to introduction of the CA water

## Preparing 2% CA Water 15m

- 2 Prepare the following:
  - 2 bottles
  - Citric acid crystals

- Sterile Gloves or weighing boats
- Scales


- 2.1 Empty one of the water bottles, place it on the scales and tare/zero
- 2.2 Pour the water from the other bottle into the one on the scales, record the weight of water
- 2.3 Calculate 2% weight of the water, this is the weight of citric acid crystals you need  
e.g. 300g of water, requires 6g of citric acid ( $300 \times 0.02 = 6$ )
- 2.4 Use the plastic side of the packet containing the sterile gloves or a weighing boat to weigh out the citric acid crystals
- 2.5 Add crystals to the bottle of water, put on the lid/spout and shake with finger covering the opening until dissolved
- 2.6 Label the bottle on 2 sides as "2% Citric Acid". Add a 'DO NOT REMOVE' label to the base of the bottle

Starting water restriction

5m

- 3 Replace normal water bottle with 2% CA bottle, perform behaviour with 10% sucrose as a reward

- 3.1 Weigh animals daily to ensure they do not drop below 80% of their baseline weight



Plot from Urai 2021 (<https://www.eneuro.org/content/8/1/ENEURO.0230-20.2020>) showing typical progression of weight loss for mice during CA water restriction

- 3.2 If any animal drops below 80% of their baseline weight, return to *ad lib* water per the PPL