



Apr 18, 2022

Immunoblotting

Will Hancock-Cerutti^{1,2,3}, Pietro De Camilli^{1,3}¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;²Interdisciplinary Neuroscience Program and MD-PhD Program, Yale University School of Medicine, New Haven, Connecticut 06510, USA;³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

1

dx.doi.org/10.17504/protocols.io.bp2l6be9zgqe/v1

William Hancock-Cerutti

This protocol describes collection of protein from cultured cells and immunoblotting, including immunoblotting of large proteins using a proprietary tris-acetate buffer system.

[dxwxbkz7.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bp2l6be9zgqe/v1Will Hancock-Cerutti, Pietro De Camilli 2022. Immunoblotting. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bp2l6be9zgqe/v1>

Immunoblot, Western blot, Tris-glycine, Tris-acetate, siRNA, cGAMP, HT-DNA, ASAPCRN

 protocol ,

Jul 08, 2021

Apr 18, 2022

Jul 08, 2021



Urmilas

Jul 14, 2021



William Hancock-Cerutti

51400

Solutions to prepare:

DMEM solution:

A	B
FBS	10%
Penicillin	100 U/ml
Streptomycin	100 mg/ml
L-glutamine	2 mM

RIPA buffer:

NaCl	150 mM
Tris	10 mM
EDTA	0.5 mM
NP40	0.50%

Supplemented immediately before use with Protease Inhibitor Cocktail (Roche) and PhosStop phosphatase inhibitor (Roche)

TBS:

A	B
Tris-Cl	50 mM
NaCl adjust pH to 7.5	150 mM

TBST: TBS with 0.1% TWEEN-20 (Sigma-Aldrich)

3x Laemmli buffer:

A	B
Tris-HCl	188 mM
SDS	3%
Glycerol	30%
Bromophenol blue	0.01%
β -mercaptoethanol	15%

Tris-glycine running buffer:

A	B
Tris base	25 mM
Glycine	192 mM
SDS in milliQ water	0.10%

Tris-glycine transfer buffer:

A	B
Tris	25 mM
Glycine	192 mM
Methanol in milliQ water	20%

Chill to 4°C

Cell culture and treatments

3d

- 1 Culture the HeLa-M cells at **37 °C** in 5% CO₂ and DMEM containing 10% FBS, **100 U/ml** penicillin, **100 mg/mL** streptomycin, and **2 Milimolar (mM)** L-glutamine (all from Gibco).
- 2 For any given experiment, plate the cells at such density so as to be approximately 90% confluent at the time of lysis.
- 3 For experiments using siRNA, transfect 60 pmols of the indicated siRNA using **6 µL** Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per well according to manufacturer protocol. Lyse the cells **72:00:00** after siRNA transfection. 3d
- 4 For experiments using cGAMP, transfect **8 µg/L** of cGAMP using **6 µL** Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per well according to manufacturer protocol.
- 5 For experiments using herring testes (HT)-DNA, transfect **1 µg** HT-DNA using **3 µL** Lipofectamine 2000 (ThermoFisher) in Opti-MEM (Gibco) per well according to manufacturer protocol.

Cell lysis and sample preparation

40m

- 6 Supplement RIPA buffer with Protease Inhibitor Cocktail (Roche) and PhosStop phosphatase inhibitor (Roche) and chill **On ice**.



Aspirate media from cells and rinse cells with PBS **On ice**. Aspirate PBS thoroughly.



Pipette RIPA lysis buffer onto cells and scrape cells using a cell lifter (Corning).



Pipette lysis buffer containing cell mass into Eppendorf tube.

10



30m

Incubate Eppendorf tube  **On ice** for  **00:30:00**.

11



Every 10 minutes, pipette lysis mixture up and down 10 times with a P-200 pipette tip (a total of 3 cycles).

NOTE: Take care not to introduce bubbles.

12



10m

Centrifuge at  **15000 x g** for  **00:10:00** at  **4 °C** and collect the post-nuclear supernatant in a new Eppendorf tube.

NOTE: Samples can be snap frozen in liquid nitrogen at this step and stored at -70°C.

13

Determine protein concentration in sample using Pierce BCA assay (ThermoFisher).

14



Prepare samples at desired concentration and add 3x Laemmli buffer.

Gel electrophoresis and immunoblotting (Tris-glycine buffer system)









4h 50m

15



5m

Incubate samples at  **95 °C** for  **00:05:00**.

- 16 During this incubation, prepare gel apparatus with Mini PROTEAN TGX 4-20% trisglycine gels (Bio-Rad) and running buffer.
- 17 Load samples into gel and run until dye front reaches bottom (90-120 V).
- 18 Remove gel and set up transfer cassette with nitrocellulose membrane.
- 19 Transfer at 120 V for  **01:00:00** at  **4 °C** in tris-glycine transfer buffer. 1h
- 20 Remove nitrocellulose membrane and stain for total protein with ponceau stain.
- 21  Wash with milliQ water.
- 22 Block membrane with 5% milk in TBST for  **01:00:00** at  **22 °C** . 1h
- 23  1h
- Incubate membrane with primary antibodies in 2.5% milk in TBST  **Overnight** at  **4 °C** .

NOTE: Optimal primary antibody incubation time and temperature can be determined empirically for a given primary antibody

24

Wash membrane with TBST. Repeat a total of 3 times.

24.1

5m

Wash membrane for  00:05:00 with TBST (1/3).

24.2

5m

Wash membrane for  00:05:00 with TBST (2/3).



24.3

5m

Wash membrane for  00:05:00 with TBST (3/3).

25

1h

Incubate membrane with secondary antibodies conjugated to IRdye 800CW or IRdye 680CW (1:10,000, Licor) in 2.5% milk in TBST for  01:00:00 at  22 °C .

26

Wash membrane with TBST. Repeat a total of 3 times.

26.1

5m

Wash membrane for  00:05:00 with TBST (1/3).

26.2

5m

Wash membrane for  00:05:00 with TBST (2/3).

26.3

5m

Wash membrane for  00:05:00 with TBST (3/3).

27 

Wash membrane with TBS. Repeat a total of 3 times.

27.1 

5m

Wash membrane for  00:05:00 with TBS (1/3).

27.2 

5m

Wash membrane for  00:05:00 with TBS (2/3).

27.3 

5m

Wash membrane for  00:05:00 with TBS (3/3).

28 

Image membranes using a Licor Odyssey Infrared Imager.



Gel electrophoresis and immunoblotting (Tris-acetate buffer system)

20h 55m

29 For VPS13C immunoblotting, lyse samples and collect the post-nuclear supernatant as above.

30  

10m

Mix post-nuclear supernatant with NuPAGE LDS Sample Buffer and Reducing Agent (ThermoFisher) and incubate for  00:10:00 at  70 °C .

31 

During this incubation, prepare gel apparatus with NuPage Tris-Acetate 3-8% gels and NuPage Running Buffer (ThermoFisher).

32 Remove gel and set up transfer cassette with nitrocellulose membrane.

33 Transfer at 0.05 mA for  **16:00:00** at  **4 °C** in NuPage transfer buffer (Thermofisher).^{16h}

34 Remove nitrocellulose membrane and stain for total protein with ponceau stain.

35 

Wash with milliQ water.

36 Block membrane with 5% milk in TBST for  **01:00:00** at  **22 °C** .^{1h}

37 

2h

Incubate membrane with primary antibodies in 2.5% milk in TBST for  **02:00:00** at  **22 °C**

NOTE: Optimal primary antibody incubation time and temperature can be determined empirically for a given primary antibody

38 

Wash membrane for with TBST. Repeat a total of 3 times.

38.1 

5m

Wash membrane for  **00:05:00** with TBST (1/3).

38.2 

5m

Wash membrane for  00:05:00 with TBST (2/3).

38.3





5m

Wash membrane for  00:05:00 with TBST (3/3).

39



1h

Incubate membrane with secondary antibodies conjugated to IRdye 800CW or IRdye 680CW (1:10,000, Licor) in 2.5% milk in TBST for  01:00:00 at  22 °C .

40



Wash membrane with TBST. Repeat a total of 3 times.

40.1



5m

Wash membrane for  00:05:00 with TBST (1/3).

40.2



5m

Wash membrane for  00:05:00 with TBST (2/3).

40.3



5m

Wash membrane for  00:05:00 with TBST (3/3).

41



Wash membrane with TBS. Repeat a total of 3 times.

41.1



5m

Wash membrane for  00:05:00 with TBS (1/3).

41.2



5m

Wash membrane for  00:05:00 with TBS (2/3).

41.3



5m

Wash membrane for  00:05:00 with TBS (3/3).

42



Image membranes using a Licor Odyssey Infrared Imager.