



DEC 07, 2023

General Microtome Sectioning of Formalin-Fixed Paraffin Embedded (FFPE) blocks

Zane

Toby J Curless^{1,2}, Jaunmuktane^{1,3,2}

¹Department of Clinical and Movement Neurosciences, UCL Institute of Neurology, Queen Square, London WC1N 3BG, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20 815, USA.;

³QSBB, 1 Wakefield Street London Greater London WC1N 1PJ, UK

OPEN ACCESS

ASAP Collaborative Research Network



DOI:

dx.doi.org/10.17504/protocols.io.bp2l69morlqe/v1

Protocol Citation: Toby J Curless, Zane Jaunmuktane 2023. General Microtome Sectioning of Formalin-Fixed Paraffin Embedded (FFPE) blocks. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l69morlqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 17, 2023



Toby J Curless

Department of Clinical and Movement Neurosciences, UCL Insti...

ABSTRACT

This protocol describes the steps for the sectioning of formalin-fixed paraffin embedded blocks in preparation for further tissue processing such as immunohistochemistry.

MATERIALS

SuperfrostTM charged microscope slides
Microtome
Fine tweezers
Fine paintbrush (optional)
Water bath
Slide holder

SAFETY WARNINGS



Follow approved local risk assessment and institutional policies on the use of sharps.

BEFORE START INSTRUCTIONS

- Check all equipment, supplies and consumables
- Set-up water bath to required temperature
- Ensure blade is fit-for-purpose and clean
- Label/barcode SuperFrost Plus slides

Last Modified: Dec 07, 2023

PROTOCOL integer ID:
77172




Keywords: ASAPCRN,
Tissue, Sectioning,
Parkinson's disease,
Microtome, Formalin-fixed
Paraffin Embedded

Funders

Acknowledgement:

The Michael J. Fox Foundation
for Parkinson's Research
(MJFF) and the Aligning
Science Across Parkinson's
(ASAP) Initiative
Grant ID: ASAP-000478

Preparation

- 1 Check all equipment, consumables and supplies are in place.
- 2 Place a clean blade into the blade holder of the microtome.
- 3 Place the FFPE block onto ice, with the face of the block on the surface of the  On ice until the block is cold through.
- 4 Heat  1 L distilled H₂O to  40 °C in a water bath.



Sectioning

12h




- 5 Once both block and water are at the correct temperature, insert the block into the chuck of the microtome.

- 6 Adjust the levers on the microtome to position the face of the block parallel with the blade.
- 7 Lower the block until it is level with the blade and move the block closer to the blade using the handle on the microtome.
- 8 Select the thickness of choice on the microtome (usually $\pm 4-7 \mu\text{m}$).
- 9 Section and trim into the block until the whole tissue cross-section is present in the sections.
- 10 Using tweezers or a fine paintbrush, gently pull sections down until a ribbon of tissue sections has formed.
- 11 Pick up the ribbon with tweezers and place into the water bath to allow wrinkles to disappear from the tissue.



A ribbon of four consecutive sections floating in the water bath to remove wrinkles.

12 Separate sections using tweezers and collect individual sections on SuperFrost™ charged microscope slides.

13 Allow slides to drain and air dry vertically at  Room temperature for >  12:00:00 prior to further  processing.

