



Mar 01, 2022

# Tissue Homogenization and GELFrEE for Top-Down Proteomics

 **GELFrEE Protocol**

Bryon Drown<sup>1</sup>, Jeannie Camarillo<sup>1</sup>, Cameron Lloyd-Jones<sup>1</sup>, Neil Kelleher<sup>1</sup>

<sup>1</sup>Northwestern University

1



[dx.doi.org/10.17504/protocols.io.b5qxq5xn](https://dx.doi.org/10.17504/protocols.io.b5qxq5xn)

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Kelleher Research Group



Kelleher KRG Research Group

Northwestern University, National Resource for Translational...

Solid human tissue is homogenized, protein is extracted, and samples prefractionated and desalted. Protein samples are ready for top-down proteomics.

DOI

[dx.doi.org/10.17504/protocols.io.b5qxq5xn](https://dx.doi.org/10.17504/protocols.io.b5qxq5xn)

Bryon Drown, Jeannie Camarillo, Cameron Lloyd-Jones, Neil Kelleher 2022. Tissue Homogenization and GELFrEE for Top-Down Proteomics. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b5qxq5xn>

Office of the Director of the National Institutes of Health

Grant ID: UH3 CA24663

National Institute of General Medical Sciences of the National Institutes of Health

Grant ID: P41 GM108569

National Institute of Cancer of the National Institute of Cancer of the National Institutes of Health

Grant ID: F32 CA246894

GELFrEE Protocol, Kelleher KRG Research Group

proteomics, top-down proteomics

protocol ,

Feb 28, 2022

Mar 01, 2022

58871

In steps of

[Overall protocol for top-down LC-MS/MS of human lung tissue](#)

[Overall protocol for top-down LC-MS/MS of human heart tissue](#)

[Overall protocol for top-down LC-MS/MS of human heart tissue](#)

[Overall protocol for top-down CZE-MS/MS of human small intestine tissue](#)

[Overall protocol for top-down CZE-MS/MS of human spleen tissue](#)

Text written in *italic* involves critical steps, important observations, relevant information and safety notes.

RIPA lysis buffer (50 mM Tris, 150 mM NaCl, 1% NP-40 (v/v), 0.5% sodium deoxycholate (w/v), 0.1% sodium dodecyl sulfate (w/v), pH 7.4, 1X Halt Protease and Phosphatase Inhibitor Cocktail (Thermo Scientific))

[HALT phosphatase and protease inhibitor cocktail \(100x\) Thermo Fisher](#)

**Scientific Catalog #78442**

10% GELFrEE cartridge

Acetone (stored in -20C freezer)

Sodium dodecyl sulfate

LCMS-grade methanol (Fisher Scientific #A456-4)

LCMS-grade water (Fisher Scientific #W64)

LCMS-grade acetonitrile (Fisher Scientific #A955-4)

HPLC-grade chloroform

LCMS-grade formic acid (Thermo Scientific #28905)

LCMS-grade acetic acid

## Cell lysis, protein precipitation, and lysate quantitation

3h

3h

### 1 Cryopulverize tissue sample with Retsch Cryomill

Cryomill

Cryogrinder

Retsch

20.749.0001



#### 1.1 Fill liquid nitrogen dewar and equilibrate to 1.3 bar

#### 1.2 Add frozen tissue sample to 25 mL grinding jar containing a 1-inch stainless steel ball.

1.3 Pulverize with three cycles of precooling with liquid nitrogen at 1 Hz for 3 min and grinding at 30 Hz for 1 min

2 Transfer pulverized tissue to a 15 mL conical tube and resuspend in 2 mL of cold RIPA cell lysis buffer (including 1x HALT protease/phosphatase inhibitor cocktail).

3 Sonicate lysate with a probe sonicator on ice (40% power, cycle 2 sec on, 3 sec off, 30 sec total).

Fisherbrand Model 120 Sonic  
Dismembrator  
Sonicator

Fisherbrand FB120110



Equipped with 1/8" probe

4 Transfer lysate to 1.5 mL Lo-bind Eppendorf tube.

[Protein LoBind Tubes, 1.5](#)

[mL Eppendorf Catalog #0030108116](#)

5 Transfer two 200 uL aliquots from each sample lysate to fresh 1.5mL Lo-bind Eppendorf tubes

6 Add 6 volumes of ice-cold acetone to the lysates, vortex at high speed (setting 8) for 10–20 s, and place the tubes in a –80 °C freezer for at least 30 min to precipitate the proteins. *The protein precipitate is stable in acetone overnight.*

Analog Vortex Mixer

Vortex mixer

Fisherbrand 02-215-365



7 Remove the tubes from the freezer and immediately centrifuge **21000 x g, 4°C, 00:30:00**

30m

- 8 Remove as much acetone as possible from the protein pellets without disturbing them. Place tubes in a biological safety cabinet for 5–10 min to evaporate off any excess acetone.
- 9 Resuspend the protein pellet in 500  $\mu\text{L}$  of 1% (wt/vol) SDS solution with pipetting and probe sonication *It may take multiple rounds of vigorous pipetting to completely dissolve the precipitates. Lysates may also be stored at  $-80^{\circ}\text{C}$  until further processing.*
- 10 Transfer 25  $\mu\text{L}$  of the solution to a new 1.5-mL Protein LoBind tube for the BCA assay. The tubes containing the remaining 200  $\mu\text{L}$  of resuspended proteins should be kept on ice.
- 11 A BCA assay should now be performed on the 25- $\mu\text{L}$  aliquots of sample proteome, using a 96-well plate and following the manufacturer's recommended protocol.
- 12 After collection of colorimetric data for BCA assay on a 96-well-plate reader, create a standard curve, and fit the average of the triplicate values for each sample dilution to the standard curve to determine the protein concentration of each sample, after correcting for dilutions.

#### Protein fractionation by GELFrEE

3h

- 13 Calculate 300  $\mu\text{g}$  protein
- 14 Transfer the calculated volume containing 300  $\mu\text{g}$  of protein to a new 1.5-mL Protein LoBind tube. *The leftover lysates in 1% (wt/vol) SDS can be stored for several months at  $-80^{\circ}\text{C}$ . From this point on, keep the sample tubes on ice until immediately before GELFrEE.*
- 15 Add 8  $\mu\text{L}$  of 1 M DTT, 30  $\mu\text{L}$  of GELFrEE sample loading buffer, and dH<sub>2</sub>O to a final volume of 150  $\mu\text{L}$ .
- 16 Vortex each tube for 10–20 s, and heat it at  $95^{\circ}\text{C}$  for 5 min. Centrifuge the samples at 1,000g for 2 min at  $4^{\circ}\text{C}$  to re-collect any condensation at the tube lids after heating
- 17 Perform fractionation by GELFrEE by following the manufacturer's protocol for 10% cartridges.

GELFrEE 8100 Fractionation System  
Fractionator  
Protein Discovery 42001

- 18 After fractions 1-5 of each sample has been collected into a new, 1.5-mL Protein LoBind tube, the fractionation process can end or more fractions can be collected, depending on the experiment. *This is a natural stopping point in the protocol, as GELFrEE fractions are stable for several months at -80 °C.*

Quality control of GELFrEE fraction by SDS-PAGE with silver staining

2h

- 19 Mix 5–10 µL of each fraction of interest with the appropriate amount of Laemmli sample buffer for the gel lane capacity in a new microcentrifuge tube. Boil these gel samples at 95 °C for 5 min and spin down (5,000g, 20–25 °C, 30–60 s).
- 20 Separate the proteins in each fraction by MW using SDS-PAGE, being careful not to run lower-MW proteins off the gel that may be of interest for TDMS or TDP characterization.
- 21 Remove the gel and submerge it in fixing solution in a clean container. Place on an orbital shaker for 20 min at low speed to fix the gel
- 22 Replace the fixing solution with enough washing solution to submerge the gel. Shake for 10 min at RT.
- 23 Replace the washing solution with sensitizing solution. Shake for 1 min at low speed to sensitize the gel.
- 24 Pour out the sensitizing solution and rinse the gel with two washes of dH<sub>2</sub>O for 1 min each.
- 25 Add enough silver nitrate solution to submerge the gel and incubate it for 20 min on an orbital shaker at low speed to stain the gel.
- 26 Pour out the silver nitrate solution and rinse the gel twice in dH<sub>2</sub>O for 1 min.
- 27 Add enough developing solution to submerge the gel, and develop under rapid shaking. *Discard the developing solution and replace it with fresh developing solution, once the liquid turns yellow, to reduce background staining.*
- 28 Pour out the developing solution and replace with terminating solution, once the desired staining intensity has been achieved, to terminate the development . The gel can now be scanned or photographed.

- 29 To each fraction, add 4 volumes of Optima-grade methanol (600  $\mu$ L), and vortex each tube vigorously for 20 s.
- 30 Add 1 volume (150  $\mu$ L) of chloroform to each fraction, and vortex vigorously for 20 s.
- 31 Add 3 volumes (450  $\mu$ L) of LC–MS-grade water to each fraction. The liquid should immediately become cloudy and white in color. Vortex vigorously for 20 s
- 32 Centrifuge the tubes at 21,000g for 10 min at 4 °C.
- 33 Carefully remove the tubes from the centrifuge. The separated aqueous and organic phases should be clearly observed, with a visible white protein pellet floating at the interface. *Be very careful when transporting tubes to avoid shaking them or otherwise disturbing the protein pellet more than is necessary.*
- 34 Remove and discard the top layer of solution, using extreme care to not remove or disturb the protein pellet. It may be necessary to leave a few millimeters of volume above the pellet to avoid accidentally pipetting it out.
- 35 Very slowly add 3 volumes (450  $\mu$ L) of Optima methanol to the tube, to avoid breaking up the protein pellet more than is necessary. Slowly mix the solution by carefully rocking the tube back and forth three to five times, while keeping an eye on the protein pellet and being careful to break it up as little as possible.
- 36 Centrifuge the tubes at 21,000g for 10 min at 4 °C
- 37 At this point, the protein pellet should now be at the bottom of the tube. Carefully remove and discard as much methanol as possible from the tube without disturbing the pellet, and repeat steps 32–34 two more times.
- 38 After the final methanol wash, remove as much methanol as possible from the tube without accidentally removing protein pellet particles. Place open tubes in a tube rack in an operating biological safety cabinet for 5–10 min to evaporate off any excess methanol.
- 39 If sample is to be analyzed by CZE-MS/MS, pellet is resuspended in 10  $\mu$ L 0.3% acetic acid. If sample is to be analyzed by LC-MS/MS, pellet is resuspended in 25  $\mu$ L of buffer A (5% acetonitrile, 94.8% water,

0.2% formic acid).

- 40 Centrifuge the samples at 21,000g for 5 min at 4 °C, and slowly remove the majority of the solution, pipetting from just below the surface to avoid disturbing any unnoticed particles. Transfer each sample to a separate, pre-chilled LC vial.
- 41 Transport the vials on ice to the nano-UHPLC auto sampler and analyze immediately to prevent sample degradation.