



Jul 20, 2022

Purification of human ATP13A2 for cryo-EM analysis

Sue Sim¹, eunyong_park¹¹University of California, Berkeley

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.261genmojg47/v1

Sue Sim

ABSTRACT

Purification of GFP-tagged human ATP13A2 expressed in Sf9 cells for cryo-EM analysis

DOI

dx.doi.org/10.17504/protocols.io.261genmojg47/v1

PROTOCOL CITATION

Sue Sim, eunyong_park 2022. Purification of human ATP13A2 for cryo-EM analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261genmojg47/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 17, 2022

LAST MODIFIED

Jul 20, 2022

PROTOCOL INTEGER ID

66878

MATERIALS TEXT

Lysis Buffer

50 mM Tris pH 7.5

200 mM NaCl

1 mM EDTA

1 mM DTT

10% glycerol

Plus protease inhibitors (5 µg/mL aprotinin, 5 µg/mL leupeptin, 1 µg/mL pepstatin A, and 2 mM PMSF)

Wash Buffer

25 mM Tris pH 7.5

100 mM NaCl

1 mM EDTA





1 mM DTT

0.03% DDM/ 0.006% CHS

Day 1: Crude membrane preparation and solubilization



7h

- 1 Thaw Sf9 cell pellets at room temperature (typical size around 10 g from 0.7L of culture)
- 2 All subsequent steps performed at **4 °C**
- 3 Resuspend each pellet in 30 mL Lysis Buffer (use 3x volume of cell pellet, 40 mL total volume)
- 4 Homogenize pellet with Dounce homogenizer, 100 plunges tight on ice
- 5 Pour lysate into pre-chilled 50 mL centrifuge tubes
- 6 Spin in centrifuge at **4000 x g, 4°C, 00:10:00** to remove unbroken cells 10m
- 7 Transfer supernatant to pre-chilled ultracentrifuge rotor tubes

- 8 Spin lysate in ultracentrifuge at  **100000 x g, 4°C, 01:30:00** to pellet membranes (Beckman Type 45 Ti rotor) 1h 30m
- 9 Resuspend membrane pellet in Lysis Buffer and final volume 1% DDM/0.2% CHS (1X pellet, 7X Lysis Buffer, 2X 5% DDM/1% CHS)
DDM: n-dodecyl-β-D-maltopyranoside (Anatrace)
CHS: cholesteryl hemisuccinate (Anatrace)
- 10 Solubilize by rotating end-over-end for  **02:30:00** at  **4 °C** 2h 30m
- 11 Clarify lysate in ultracentrifuge at  **100000 x g, 4°C, 01:00:00** 1h

Day 1: Bead binding and overnight 3C cleavage

3h

- 12 Equilibrate 1 mL Sepharose beads conjugated with anti-GFP nanobody with Wash Buffer
- 13 Add beads to gravity column and wash with 10 mL wash buffer
- 14 After ultracentrifugation is complete, transfer supernatant into 50 mL falcon tube
- 15 Add 1 mL equilibrated anti-GFP nanobody beads to tube
- 16 Incubate by rotating end-over-end for  **02:30:00** at  **4 °C** 2h 30m
- 17 Transfer to gravity column and let flow-through drain

- 18 Wash beads with 30 column volumes of Wash Buffer
- 19 Add 5 mL Wash Buffer to gravity column and 10 µg/mL HRV 3C protease
- 20 Incubate by rotating end-over-end overnight at **4 °C**

Day 2: SEC column 3h

- 21 Equilibrate Superose 6 Increase 10/300 GL column with Wash Buffer
- 22 Concentrate the protein to 0.5 mL using an Amicon Ultrafilter (cut-off 100kDa)
- 23 After concentration, spin protein at **17000 x g, 4°C, 00:10:00** 10m
- 24 Injected sample into FPLC
- 25 Collect peak fractions and concentrate to approximately 5-7 mg/mL for cryo-EM grid preparation