

Sep 06, 2022

# Subcloning of genotype-confirmed hPSCs clones

In 1 collection

Hanqin Li<sup>1</sup>, Yogendra Verma<sup>1</sup>, Dirk Hockemeyer<sup>1</sup>, Frank Soldner<sup>2</sup><sup>1</sup>University of California, Berkeley; <sup>2</sup>Albert Einstein College of Medicine

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.b4mkqu4w](https://dx.doi.org/10.17504/protocols.io.b4mkqu4w)

Devin E Snyder

## ABSTRACT

This protocol describes a standard procedure for subcloning of genotype-confirmed human pluripotent stem cells (hPSCs).

## General notes:

1. Throughout this protocol, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.

DOI

[dx.doi.org/10.17504/protocols.io.b4mkqu4w](https://dx.doi.org/10.17504/protocols.io.b4mkqu4w)

## PROTOCOL CITATION

Hanqin Li, Yogendra Verma, Dirk Hockemeyer, Frank Soldner 2022. Subcloning of genotype-confirmed hPSCs clones. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b4mkqu4w>



## FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's  
Grant ID: ASAP-000486

## COLLECTIONS



**Standard operating procedure for the isolation of genetically engineered hPSCs clones in a high-throughput way**

## KEYWORDS

ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 03, 2022

LAST MODIFIED

Sep 06, 2022

PROTOCOL INTEGER ID

57740

PARENT PROTOCOLS

Part of collection

[Standard operating procedure for the isolation of genetically engineered hPSCs clones in a high-throughput way](#)

## MATERIALS TEXT

| A  | B                    | C         |
|--|----------------------|-----------|
| Item                                     | Vendor               | Catalog # |
| DMEM/F12                                 | Thermo<br>Fisher     | 11320082  |
| DPBS w/o<br>Calcium and magnesium (DPBS) | Corning              | MT21031CV |
| Fetal Bovine<br>Serum (FBS)              | Corning              | 35-011-CV |
| Knockout Serum Replacement               | Thermo<br>Fisher     | 10828-028 |
| L-Glutamine                              | Sigma                | G8540     |
| Penicillin & Streptomycin (100X)         | Thermo<br>Fisher     | 15140163  |
| MEM Non-Essential Amino Acids<br>(100X)  | Thermo<br>Fisher     | 11140050  |
| Heat Stable Recombinant Human<br>FGF2    | Thermo<br>Fisher     | PHG0360   |
| Y-27632                                  | Chemdea              | CD0141    |
| 2-Mercaptoethanol                        | Sigma                | M3148     |
| 0.25% Trypsin with EDTA<br>(Trypsin)     | Thermo<br>Fisher     | 25200114  |
| Proteinase K                             | Sigma                | P6556     |
| DMSO                                     | Fisher<br>Scientific | BP231-100 |

## 1 Change medium to hPSCs medium + Rock inhibitor one day before subcloning

### 1.1 hPSC medium

| A   | B      |
|---|--------|
| DMEM/F12                                      | 385 ml |
| Fetal Bovine Serum (FBS)                      | 75 ml  |
| Knockout Serum Replacement                    | 25 ml  |
| L-Glutamine (100X)                            | 5 ml   |
| Penicillin & Streptomycin (100X)              | 5 ml   |
| MEM Non-Essential Amino Acids (100X)          | 5 ml   |
| 2-Mercaptoethanol (10,000X)                   | 50 µl  |
| Heat Stable Recombinant Human FGF2 (25ug/ml)* | 80 µl  |

\*While we prefer Heat Stable Recombinant Human FGF2, we also have used regular FGF2. Final volume: 500ml

### L-Glutamine (100X)

|                     |        |
|---------------------|--------|
| L-Glutamine, powder | 14.6 g |
| MilliQ H2O          | 500 ml |

### 2-Mercaptoethanol (10,000X)

|                   |         |
|-------------------|---------|
| 2-Mercaptoethanol | 0.78 ml |
| MilliQ H2O        | 9.22 ml |

### Heat Stable Recombinant Human FGF2 (25µg/ml)

| A                                  | B      |
|------------------------------------|--------|
| Heat Stable Recombinant Human FGF2 | 500 µg |
| 0.1% BSA                           | 20 ml  |

Final volume: 20ml



### Y-27632 (1,000X)

|         |         |
|---------|---------|
| Y-27632 | 5 mg    |
| DMSO    | 1.56 ml |

### hPSC medium + Rock inhibitor



| A                   | B      |
|---------------------|--------|
| hPSCs medium        | 500 ml |
| Y-27632<br>(1,000X) | 500 µl |

Final volume: 500ml

- 2 Wash the genotype-confirmed wells with DPBS
- 3 Add 25 µl Trypsin to those wells
- 4 Incubate  00:05:00  37 °C 5m
- 5 Add 75 µl hPSCs medium + Rock inhibitor, mix well by pipetting
- 6 Transfer dissociated cells to a 15 conical tube
- 7 Prepare three MEF wells of a 6-well plate. Seed 1/100, 1/1,1000 and the rest of the cells respectively to those wells
- 8 Change medium daily for the high density well starting from day 3. Change medium every 3 days for the low density wells in the first week, then every other day in the second week
- 9 Once the high density wells grow to 50-70% confluency, freeze.

For a detailed protocol on freezing hPSCs, refer to "Freezing of hPSCs grown on MEFs" in the "Thawing, Passaging and Freezing of hPSCs on MEFs" collection;  
[dx.doi.org/10.17504/protocols.io.b4msqu6e](https://dx.doi.org/10.17504/protocols.io.b4msqu6e)

For a detailed protocol on freezing feeder-free hPSCs, refer to "Freezing of feeder-free hPSCs" in the "Feeder-free culturing of hPSCs" collection; [dx.doi.org/10.17504/protocols.io.b4mcqu2w](https://dx.doi.org/10.17504/protocols.io.b4mcqu2w)

- 10 When big hPSCs colonies form in the low density wells, change medium to hPSCs medium + Rock inhibitor
- 11 The day after, proceed with manual colony picking
- 12 For each original clone, prepare six microcentrifuge tubes pre-added with 20 µl Trypsin
- 13 Aspirate medium and add 2 ml DMEM/F12
- 14 Change medium to DPBS for the well where colonies will be picked
- 15 Under the dissecting microscope, use mouth pipet or a fine 10 µl tip to pick one undifferentiated colony which is fully separated from other colonies.
- 16 Transfer the colony to one microcentrifuge tube from step 12.
- 17 Repeat step 15-16 to pick another five colonies
- 18 Incubate the microcentrifuge tubes at  37 °C  00:05:00 5m
- 19 Add 70 µl hPSCs medium + Rock inhibitor to each tube, pipet to mix
- 20 Seed all cells into one well of a 12-well MEF plate

21 Shake plates to distribute cells evenly

22 Culture for 7-10 days with medium change daily from day 3

23 When the subclones grow to 50-70% confluence, passage and prepare crude cell lysis for NGS genotyping

### 23.1 Crude lysis buffer (2x)

| A                             | B         |
|-------------------------------|-----------|
| KCl                           | 100 mM    |
| MgCl <sub>2</sub>             | 4 mM      |
| NP-40                         | 0.9%      |
| Tween-20                      | 0.9%      |
| Tris                          | 20 mM     |
| Proteinase K (add before use) | 100 µg/ml |

pH: 8

24 Once genotype confirmed, expand and freeze the one subclone, which shows the best hPSCs morphology and proliferates normally.

25 Test for mycoplasma, stain for pluripotent markers, and karyotyping