




Sep 08, 2020

# Open Isothermal Platform Protocol

Andrew Buultjens<sup>1</sup>, Koen Vandelannoote<sup>1</sup>, Tim Stinear<sup>1</sup><sup>1</sup>The University of Melbourne**1** Works for me dx.doi.org/10.17504/protocols.io.bk3nkyme

Open-Isothermal-Platforms Dx

Open-Isothermal-Platforms Dx

 buultjensa

DOI

[dx.doi.org/10.17504/protocols.io.bk3nkyme](https://dx.doi.org/10.17504/protocols.io.bk3nkyme)

PROTOCOL CITATION

Andrew Buultjens, Koen Vandelannoote, Tim Stinear 2020. Open Isothermal Platform Protocol.

**protocols.io**<https://dx.doi.org/10.17504/protocols.io.bk3nkyme>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 08, 2020

LAST MODIFIED

Sep 08, 2020

PROTOCOL INTEGER ID

41806

- 1 Extract RNA as per [dx.doi.org/10.17504/protocols.io.bkynkxve](https://dx.doi.org/10.17504/protocols.io.bkynkxve)
- 2 Transfer 5ml of extracted RNA to pre-aliquoted LAMP mix. Mix by pipetting up and down 5 times.
- 3 Add 5ml of H<sub>2</sub>O in negative control well 7.
- 4 Add 5ml of positive control in well 8.
- 5 Cap all wells

6 Place strips in FABL-8

7 Run FABL-8 as per manual

## 8 LAMP reaction (can refer to N1-STOP-LAMP paper)

Single reaction setup:

Reagent	Volume (ml)
Optigene Master Mix	15
N1 primer	5
Extracted RNA	5
<b>Total</b>	<b>25</b>

To reduce pipetting error, Master Mix and N1 primer are premixed according to the number of reactions desired (allow for pipetting error by preparing enough for  $\sim 1.1 \times$  the number of reactions, rounded to the nearest whole number); 20ml of this mixture is then aliquoted per reaction.