



Oct 24, 2020

PCR

Zhujun Wei¹¹Northeast Forest University

1

Works for me

dx.doi.org/10.17504/protocols.io.bnwqmfwdw

2020 iGEM NEFU China



Zhujun Wei

DOI

dx.doi.org/10.17504/protocols.io.bnwqmfwdw

PROTOCOL CITATION

Zhujun Wei 2020. PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bnwqmfwdw>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 24, 2020

LAST MODIFIED

Oct 24, 2020

PROTOCOL INTEGER ID

43696

MATERIALS TEXT

2×high Taq Master Mix
ddH₂O
Template
F/R Primer
Bio-rad S1000TM Thermo Cycler.

SAFETY WARNINGS

Please wear gloves for the experiment, don't try to touch the lid after PCR program initiation.

BEFORE STARTING

Set up a small box with ice, put DNA and 2×high Taq Master Mix into it before going into the Bio-rad S1000TM thermocycler. Set up a small box with ice, put DNA and 2×high Taq Master Mix into it before going into the Bio-rad S1000TM Thermo Cycler.

1 Program the thermocycler as follows:

Temperature	Time
98°C	5 min
98°C	15 s
Tm-4°C	15 s

72°C	(5-10s) × (1 kb)
72°C	5 min
4°C	∞

2 Add the following reagent to a PCR tube.(50 µl).

2×high Taq Master Mix (Enzyme)	25 µl
Template	2 µl
Forward Primer (10 µM)	2 µl
Reverse Primer (10 µM)	2 µl
ddH ₂ O	19 µl