



Oct 18, 2020

Plaque PCR

Jiaxin Li¹¹South China University of Technology

1

Works for me

dx.doi.org/10.17504/protocols.io.bni9mch6

Jiaxin Li

DOI

dx.doi.org/10.17504/protocols.io.bni9mch6

PROTOCOL CITATION

Jiaxin Li 2020. Plaque PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bni9mch6>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 18, 2020

LAST MODIFIED

Oct 18, 2020

PROTOCOL INTEGER ID

43329

- 1 Mark independent plaques on the back of the Petri dish with a black pen. Take 6 points in a plate and try to take points randomly and evenly.
- 2 Take 24 1.5 ml EP tubes.
 - A) Add 200 µl SM Buffer per tube
 - B) Add 2 µl chloroform to each tube

Sampling

- 3 Suck out the plaque just spotted with the pipette
- 4 Put into the prepared EP tube.
- 5 Open the metal bath to 95 degrees for lysis of the phage
- 6 Place the EP tube at room temperature for one hour to dissolve the sample sufficiently and mix several times during the period.

PCR system

- 7 KOD FX 20ul system
- 8 Add 18.4 μ l to each tube (1.6 μ l template is needed)
- 9 After an hour, put the EP tube in a metal bath and heat it for 10 minutes.
- 10 Centrifugation for 5 min, 12,000 rpm
- 11 Take 1.6ul supernatant from EP tube and add it to new EP tube.
- 12 Centrifuge the EP tube briefly with a centrifuge.
- 13 Place in PCR instrument.