



MAR 07, 2024

## 🌐 Plasmid construction

📁 In 1 collection

Dan Tudorica<sup>1</sup>

<sup>1</sup>University of California, Berkeley



Dan Tudorica  
Hurley Lab, QB3, UC Berkeley

### ABSTRACT

Via restriction enzyme digest and ligation

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.yxmvm3z5bl3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvm3z5bl3p/v1)

**Protocol Citation:** Dan Tudorica  
2024. Plasmid construction.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.yxmvm3z5bl3p/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Nov 09, 2023

**Last Modified:** Mar 07, 2024

- 1 Design PCR primers such that they contain a BglII restriction site upstream of the gene, and Sall downstream of the gene. Design cut sites such that gene will still be in frame post ligation.
- 2 Acquire purified pmCherryC1 stock, and cleave using BglII and Sall using the standard NEB protocol. Simultaneously cleave PCR product
- 3 Purify both vector and insert using 1% agarose gel, and compare to uncut controls to make sure that restriction digest went to completion. Gel extract vector and insert.
- 4 Ligate with T4 ligase at room temperature overnight following standard NEB protocol.
- 5 Transform into XL10 Gold cells, plate on ampicillin selective medium, grow overnight at 37 C
- 6 Pick colonies, grow 5 mL overnight culture, and purify plasmid using a commercially available plasmid purification kit.