



Sep 12, 2022

# Twist 96-Plex (Riptide) Library Prep

Katarina A Cohen<sup>1,2</sup>, Khai-Minh H Nguyen<sup>1,2</sup>, Oksana Polesskaya<sup>1,2</sup>,  
Abraham Palmer<sup>1,2</sup>

<sup>1</sup>UC San Diego; <sup>2</sup>Palmer Lab

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.j8nlkkm85l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkkm85l5r/v1)

CGORD

Palmer Wet-Lab Protocols



Khai-Minh H Nguyen

## ABSTRACT

This protocol is designed for Twist 96-PLex Library Prep. We use the EPmotion 5075 to add sample barcodes (can also easily be done manually with a multichannel pipette). This is a continuation of the "EPMotion - Normalization and Randomization" protocol.

## DOI

[dx.doi.org/10.17504/protocols.io.j8nlkkm85l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkkm85l5r/v1)

## PROTOCOL CITATION

Katarina A Cohen, Khai-Minh H Nguyen, Oksana Polesskaya, Abraham Palmer  
2022. Twist 96-Plex (Riptide) Library Prep . **protocols.io**  
<https://protocols.io/view/twist-96-plex-riptide-library-prep-b8khrut6>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 02, 2022

## LAST MODIFIED

Sep 12, 2022

## PROTOCOL INTEGER ID

61801

## GUIDELINES

It is optional to use the EPmotion in this protocol. Visit Twist Bioscience site for updated protocols.

## MATERIALS TEXT

Please see Twist 96 Plex Protocol on Twist Bioscience Site for full list of materials needed

### Equipment

- epMotion 5075
- 8 Channel Dispensing Tool (50uL) **Catalog No. 960001044**
- Qubit Fluorometer
- Bioanalyzer
- Nanodrop
- Pipettes
- Multichannel pipette (10ul) if EPmotion is not available.

### Consumables

- 1-50uL epT.I.P.S. Motion **Catalog No. 0030014413**
- Qubit HS Assay
- TapeStation D1000
- Reagent Wells
- Ultrapure DI Water
- Pipette Tips: 10ul, 200ul, 1000ul

### Adding Primer A with EPmotion


- 1 Review the Twist 96-Plex Protocol  
Note: Protocol was created for Document Version (DOC-001284 REV 2.0)  
[Twist Bioscience Protocol](#)
- 2 Fill ice pan with ice chips.
- 3 Defrost 4uL randomized/normalized sample plate from the "EPMotion - Normalization and Randomization" protocol **⚠ On ice**
- 4 Open epBlue application

5 Download and import ☐ **50% GC Content Adapters.export7** and


☐ **Low GC Content Adapters.export7**

- Depending on the GC content of the species you are working with, choose the protocol accordingly.
- Use 50% GC Content Adapters if samples have GC content between 40% - 60%
- Use Low GC Primers if samples have GC content below 40%
- Use High GC Primers if samples have GC content above 60%

6 Set up the EPmotion as shown on the worktable below



Worktable for High and Low GC Content Primers - 50%



Worktable for Only Low GC Content Primers

## Library Prep

### 7 Follow the Twist 96-Plex Protocol

Note: Protocol was created for Document Version (DOC-001284 REV 2.0)

- **IMPORTANT:** Size selection Bead Volumes used: **50ul:20ul (Option 3)** worked best for us.
- Make sure to record which pool barcode index for each library
- PRO TIP: We try to match the last digit of the library number (Ex. Riptide 53) with the pool barcode index used.
- In the case of Riptide **53**. We would use the pool barcode #3 for this particular library.

## QC

### 8 Perform on each library and record keep all QC for pooling steps

- Nanodrop
- Qubit (HS Assay)
- Tapestation (D1000)

### 9 Libraries should have an average fragment size between 420bp - 650bp.

- 260/280 should be around 1.80 - 2
- 260/230 should be around 2-2.2.
- **We have been able to get good data from libraries with relatively poor nanodrop purities.**
- Qubit concentrations can widely range. We get a range from 10ng/ul - 60ng/ul