



May 24, 2024

🌐 Labeling of microtubules using mouse anti- β -tubulin primary monoclonal antibody with secondary Fe-TAML-peg4-Cy5-goat anti-mouse IgG conjugate and oxidation of DAB with H₂O₂ for light and transmission electron microscopy

DOI

dx.doi.org/10.17504/protocols.io.261ge5dpdg47/v1

Steven Adams^{1,2}, Mason R. Mackey^{1,2}, Mark Ellisman^{1,2}

¹UC San Diego Health Sciences; ²National Center for Microscopy and Imaging Research

NCMIR@UCSD



National Center for Microscopy and Imaging Research

University of California, San Diego

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.261ge5dpdg47/v1**

External link: **<https://www.biorxiv.org/content/10.1101/2023.08.25.554352v2>**

Protocol Citation: Steven Adams, Mason R. Mackey, Mark Ellisman 2024. Labeling of microtubules using mouse anti- β -tubulin primary monoclonal antibody with secondary Fe-TAML-peg4-Cy5-goat anti-mouse IgG conjugate and oxidation of DAB with H₂O₂ for light and transmission electron microscopy. **protocols.io** **<https://dx.doi.org/10.17504/protocols.io.261ge5dpdg47/v1>**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working.

Created: May 17, 2024

Last Modified: May 24, 2024

Protocol Integer ID: 100230

Keywords: TEM, DAB oxidation, oxidation reaction, Correlated Light and Electron Microscopy, Fe-TAML, small molecule peroxidase

Funders Acknowledgement:

R01

Grant ID: GM138780 (MHE)

R01

Grant ID: GM086197

(SRA/DB),

R35

Grant ID: GM128859 (JTN)

NSF

Grant ID: 2014862 (MHE)

R01

Grant ID: AG081037 (MHE)




Abstract

This protocol details the labeling of microtubules using mouse anti- β -tubulin primary monoclonal antibody with secondary Fe-TAML-peg4-Cy5-goat anti-mouse IgG conjugate and oxidation of DAB with H_2O_2 for light and transmission electron microscopy.

Materials

CSB buffer:

A	B
Pipes buffer	10 mM
NaCl	150 mM
EGTA	5 mM
MgCl ₂	5 mM
Glucose monohydrate, pH 6.8	5 mM

-  Anti- β -Tubulin antibody, Mouse monoclonal, clone TUB 2.1 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T5201**
-  3,3'-Diaminobenzidine **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8001- 10G**
-  Durcupan™ ACM **Merck MilliporeSigma (Sigma-Aldrich) Catalog #44610**
- Fe-TAML-peg4-Cy5-goat anti-mouse IgG
- BSA **Sigma Catalog #A8022-100G**



Labeling of microtubules

2d 17h 59m

- 1 Culture HEK293T cells on MatTek plates containing poly-D-lysine coated glass bottom No.0 coverslips in DMEM supplement with 10% fetal bovine serum.
- 2 Rinse the cells (x3) with cytoskeleton stabilizing buffer, at 37 °C and fixed with 4% paraformaldehyde (19202, Electron Microscopy Sciences) and 0.05% glutaraldehyde (16220, Electron Microscopy Sciences) in CSB at 37 °C for 00:05:00 and for 00:25:00 at 4 °C .

30m

CSB buffer:

A	B
Pipes buffer	10 mM
NaCl	150 mM
EGTA	5 mM
MgCl2	5 mM
Glucose monohydrate, pH 6.8	5 mM

- 3 After fixation, first wash the cells with CSB (5 x 1 min) at 4 °C .
- 3.1 After fixation, first wash the cells with CSB for 00:01:00 at 4 °C (1/5)
- 3.2 After fixation, first wash the cells with CSB for 00:01:00 at 4 °C (2/5)
- 3.3 After fixation, first wash the cells with CSB for 00:01:00 at 4 °C (3/5)
- 3.4 After fixation, first wash the cells with CSB for 00:01:00 at 4 °C (4/5)
- 3.5 After fixation, first wash the cells with CSB for 00:01:00 at 4 °C (5/5)



1m



























1m

1m



1m

1m



- 4 Treat with 0.1% saponin and 0.05% glycine in CSB for 20 mins at  4 °C while on a rocker.
- 5 Wash the cells in CSB buffer with 0.05% glycine (3 x 1 min) at 4°C. 
- 5.1 Wash the cells in CSB buffer with 0.05% glycine  00:01:00 at  4 °C . (1/3) 1m 
- 5.2 Wash the cells in CSB buffer with 0.05% glycine  00:01:00 at  4 °C . (2/3) 1m 
- 5.3 Wash the cells in CSB buffer with 0.05% glycine  00:01:00 at  4 °C . (3/3) 1m 
- 6 Block the cells with 1% BSA (A8022-100G, Sigma), 1% normal goat serum (NGS) and 0.05% glycine in CSB for  00:20:00 at  4 °C . 20m
- 7 Incubate the cells with primary mouse monoclonal antibody to β -tubulin (300-fold dilution, clone Tub2.1, T5201, Sigma) for  03:00:00 at  4 °C in 1% BSA, 1% NGS and 0.05% glycine in CSB buffer. 3h 
- 8 Remove primary antibody and wash with 1% BSA, 1% NGS and 0.05% glycine in CSB (5 x 3 min) at 4°C. 
- 8.1 Wash with 1% BSA, 1% NGS and 0.05% glycine in CSB  00:03:00 at  4 °C . (1/5) 3m 
- 8.2 Wash with 1% BSA, 1% NGS and 0.05% glycine in CSB  00:03:00 at  4 °C . (2/5) 3m 
- 8.3 Wash with 1% BSA, 1% NGS and 0.05% glycine in CSB  00:03:00 at  4 °C . (3/5) 3m 



- 8.4 Wash with 1% BSA, 1% NGS and 0.05% glycine in CSB  00:03:00 at  4 °C . (4/5) 3m 
- 8.5 Wash with 1% BSA, 1% NGS and 0.05% glycine in CSB  00:03:00 at  4 °C . (5/5) 3m 
- 9 Incubate the cells with secondary Fe-TAML-peg4-Cy5-goat anti-mouse IgG conjugate in 1% BSA (0.15 ml diluted to 1ml) in 1% NGS and 0.05% glycine in CSB for  Overnight at  4 °C . 8h 
- 10 Then wash the cells (5 x 1 min) with CSB at 4°C. 
- 10.1 Wash the cells for  00:01:00 with CSB at  4 °C . (1/5) 1m 
- 10.2 Wash the cells for  00:01:00 with CSB at  4 °C . (2/5) 1m 
- 10.3 Wash the cells for  00:01:00 with CSB at  4 °C . (3/5) 1m 
- 10.4 Wash the cells for  00:01:00 with CSB at  4 °C . (4/5) 1m 
- 10.5 Wash the cells for  00:01:00 with CSB at  4 °C . (5/5) 1m 
- 11 Fix the cells with 2% glutaraldehyde in CSB for  00:20:00 at  4 °C . 20m
- 12 Wash the cells (5 x 1 min) with CSB at 4°C. 



- 12.1 Wash the cells for 00:01:00 with CSB at 4 °C . (1/5) 1m
- 12.2 Wash the cells for 00:01:00 with CSB at 4 °C . (2/5) 1m
- 12.3 Wash the cells for 00:01:00 with CSB at 4 °C . (3/5) 1m
- 12.4 Wash the cells for 00:01:00 with CSB at 4 °C . (4/5) 1m
- 12.5 Wash the cells for 00:01:00 with CSB at 4 °C . (5/5) 1m
- 13 Image the cells for fluorescence labeling.
- 14 Dissolve and add 5.4 mg of 3,3'- Diaminobenzidine (DAB) (D8001-10G, Sigma-Aldrich) in 1.0 mL of [M] 0.1 Mass Percent HCl and 9.0 mL of [M] 50 millimolar (mM) Bicine [M] 100 millimolar (mM) NaCl pH 8.3 with 10 µL H₂O₂ (final, [M] 40 millimolar (mM) from 30% stock) to the DAB solution.
- 15 Wash the cells (2 x 2 min) with 50 mM Bicine 100 mM NaCl pH 8.3 at 4°C.
- 15.1 Wash the cells for 00:02:00 with [M] 50 millimolar (mM) Bicine [M] 100 millimolar (mM) NaCl pH 8.3 at 4 °C . (1/2) 2m
- 15.2 Wash the cells for 00:02:00 with [M] 50 millimolar (mM) Bicine [M] 100 millimolar (mM) NaCl pH 8.3 at 4 °C . (2/2) 2m
- 16 Add the DAB/H₂O₂ solution to the cells by a 0.22µm Millex 33mm PES sterile filter (SLGSR33RS, Sigma-Aldrich) at Room temperature . Reaction time is 01:00:00 - 1h 30m



01:30:00 .

17 Remove the DAB solution, and wash the cells with 50 mM Bicine 100mM NaCl pH 8.3 (2 x 2 min) on ice.

17.1 Wash the cells with [M] 50 millimolar (mM) Bicine [M] 100 millimolar (mM) NaCl pH 8.3 for 00:02:00 On ice . (1/2) 2m

17.2 Wash the cells with [M] 50 millimolar (mM) Bicine [M] 100 millimolar (mM) NaCl pH 8.3 for 00:02:00 On ice . (2/2) 2m

18 Wash the cells with 0.1 M sodium cacodylate buffer pH 7.4 (3 x 2 min) on ice.

18.1 Wash the cells with [M] 0.1 Molarity (M) sodium cacodylate buffer pH 7.4 for 00:02:00 . (1/3) 2m

18.2 Wash the cells with [M] 0.1 Molarity (M) sodium cacodylate buffer pH 7.4 for 00:02:00 . (2/3) 2m

18.3 Wash the cells with [M] 0.1 Molarity (M) sodium cacodylate buffer pH 7.4 for 00:02:00 . (3/3) 2m

19 Do a final primary fixation with 2% glutaraldehyde in [M] 2 millimolar (mM) CaCl_2 [M] 0.1 Molarity (M) sodium cacodylate pH 7.4 for 00:30:00 at 4 °C . 30m

20 Remove the fixative and wash the cells (5 x 2 min) with 0.1 M sodium cacodylate (18851, Ted Pella) pH 7.4 at 4°C.

20.1 Wash the cells for 00:02:00 with [M] 0.1 Molarity (M) sodium cacodylate (18851, Ted Pella) pH 7.4 at 4 °C . (1/5) 2m



- 20.2 Wash the cells for 00:02:00 with 0.1 Molarity (M) sodium cacodylate (18851, Ted Pella) pH 7.4 at 4 °C . (2/5) 2m
- 20.3 Wash the cells for 00:02:00 with 0.1 Molarity (M) sodium cacodylate (18851, Ted Pella) pH 7.4 at 4 °C . (3/5) 2m
- 20.4 Wash the cells for 00:02:00 with 0.1 Molarity (M) sodium cacodylate (18851, Ted Pella) pH 7.4 at 4 °C . (4/5) 2m
- 20.5 Wash the cells for 00:02:00 with 0.1 Molarity (M) sodium cacodylate (18851, Ted Pella) pH 7.4 at 4 °C . (5/5) 2m
- 21 All cells are post-fixed with 1% osmium tetroxide (19150, Electron Microscopy Sciences) containing 0.8% potassium ferrocyanide, 2 millimolar (mM) CaCl₂ and in 0.1 Molarity (M) sodium cacodylate pH 7.4 for 00:30:00 at 4 °C . 30m
- 22 Wash the cells (5 x 2 min) with ddH₂O at 4°C.
- 22.1 Wash the cells for 00:02:00 with ddH₂O at 4 °C . (1/5) 2m
- 22.2 Wash the cells for 00:02:00 with ddH₂O at 4 °C . (2/5) 2m
- 22.3 Wash the cells for 00:02:00 with ddH₂O at 4 °C . (3/5) 2m
- 22.4 Wash the cells for 00:02:00 with ddH₂O at 4 °C . (4/5) 2m
- 22.5 Wash the cells for 00:02:00 with ddH₂O at 4 °C . (5/5) 2m



- 23 Dehydrate the cells by an ice-cold graded dehydration ethanol series of 20%, 50%, 70%, 90%, 100% (anhydrous) for 00:01:00 each and 3 x 100% (anhydrous) at Room temperature for 00:01:00 each. 2m
- 24 Infiltrate the cells with one-part Durcupan ACM epoxy resin (44610, Sigma-Aldrich) to one-part anhydrous ethanol (1:1) for 00:30:00 . 30m
- 25 Infiltrate the cells 3 times with 100% Durcupan resin for 01:00:00 each. 1h
- 26 Do a final change of Durcupan resin and immediately place cells in a vacuum oven at 60 °C for 48:00:00 to harden. 2d

Protocol references

Adams, Stephen R., et al. "Fe-TAMLs as a new class of small molecule peroxidase probes for correlated light and electron microscopy." *bioRxiv* (2023): 2023-08.