



Jan 15, 2021

Site Directed Mutagenesis 2016

Elizabeth Fozo¹¹In-house protocol

1

Works for me

This protocol is published without a DOI.



Eadewunm

ABSTRACT

Site-directed Mutagenesis

PROTOCOL CITATION

Elizabeth Fozo 2021. Site Directed Mutagenesis 2016. [protocols.io](https://protocols.io/view/site-directed-mutagenesis-2016-brgim3ue)
<https://protocols.io/view/site-directed-mutagenesis-2016-brgim3ue>

KEYWORDS

Site Directed Mutagenesis, Directed Mutagenesis, Mutagenesis

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 14, 2021

LAST MODIFIED

Jan 15, 2021

PROTOCOL INTEGER ID

46314

GUIDELINES

- Primers
- QuikChange PCR (enzyme from Stratagene)
- DpnI Digest
- Transformation

DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Site-directed Mutagenesis

BEFORE STARTING

Primers: Design 40-mer primers with mutation in center. Can mutate up to several aa's in a row (make longer

primers then). Order the forward and reverse complement.

QuikChange PCR (enzyme from Stratagene):

A	B
Template:	250 ng
10 μ M fwd primer:	2 μ L
10 μ M rev primer:	2 μ L
10x LA Taq or Pfu buffer:	5 μ L
10mM dNTPs:	1 μ L
dH ₂ O:	to 49 μ L
LA Taq or Pfu Turbo polymerase:	1 μ L (We use LA taq)
PCR program:	95°C - 3 minutes
	*95°C - 30 seconds
	*55 °C - 1 minutes
	68 °C - 2-2.5 min/kb
	4°C - hold

Note:

1. Newer QuikChange kit says 1 min/kb, but I always still use 2 min/kb and it works well.

2. If making two separate mutations with two sets of primers, can try using both sets in one reaction (as well as separately). I have had this work before and it is much faster.

* 18x

DpnI Digest:

- 1 Transfer PCR to Epp tube.
- 2 Remove 2 μ L from PCR.
- 3 Add 2 μ L DpnI enzyme to remaining 48 μ L.
- 4 Incubate 1+ hrs 37°C (1st digestion).
- 5 Add another 2 μ L DpnI enzyme to remaining 48 μ L (2nd digestion).
- 6 Use PCR purification kit to purify the 2nd digestion product.

After PCR purification, you will lose around 20% of your product. I always try to concentrate my 2nd digestion PCR product to a volume of around 40 μ L before I do the transformation in order to keep the same concentration between the undigested product and the digested product.

- 7 Store at -20 °C.

You can go ahead to do the transformation using the undigested PCR product and the 2nd digestion product. If you don't get quite good result after the transformation, such as no significant difference between the undigested one and the 2nd digestion one, you can do the 3rd digestion. But when you do the 3rd digestion, you need to add restriction enzyme buffer 4.

Transformation:

- 8 Make TOP10 electrocompetent cells as usual.
- 9 Transform 1 μ L DpnI-digested plasmid and 1 μ L non-DpnI plasmid as control.
- 10 Plate 1/10, 9/10 on LB + antibiotic plates.
- 11 Should see many more colonies on the -DpnI plate if the digest of the parental unmutated plasmid worked well.
- 12 Streak a minimum of 6.
- 13 Sequence 6.
- 14 Check that
 1. the desired mutation is present .
 2. the primer region did not get duplicated.Sometimes you will see that somehow the 60-mer of the primer is present in two or more tandem repeats.