

Dec 01, 2021

Part 2: Custom QXT



1

cecilia¹, Suzie Alarcon¹, Alessandro Sette¹

¹La Jolla Institute for Immunology

1



dx.doi.org/10.17504/protocols.io.bxr5pm86



Yaqian Xu

This protocol details the procedure of Custom QXT.

[d2bybg4xf.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.bxr5pm86

cecilia , Suzie Alarcon, Alessandro Sette 2021. Part 2: Custom QXT.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bxr5pm86>



SmartSeq

Custom QXT, Fragmentation, Adaptor-tagging, Ampure bead purification 1, Amplification, Ampure bead purification 2, ASAPCRN

protocol ,

Aug 27, 2021

Dec 01, 2021

Aug 27, 2021

Urmilas

Nov 30, 2021

Yaqian Xu

52765

Materials:

- cDNA
- 5x TD Buffer
- QXTEnzyme
- StorageSoln
- 80% EtOH
- Ampure beads/sample

A: Fragmentation and Adaptor-tagging

1m

- 1 FIRST, IF POSSIBLE:** Normalize cDNA samples to 2x the input mass (for **1 ng** cDNA input will be **2 ng / 10.13 µL**, for the **5 ng** input will be **10 ng / 10.13 µL**). Run quant (picogreen) to assess concentration. If normalization was successful, add one volume of water to the samples to bring it back to 1x concentration and continue with **10.16 µL**.

A			B		
1ng cDNA input:			5ng cDNA input:		
<input type="checkbox"/>			<input type="checkbox"/>		
INPUT: 1ng cDNA in 10.16uL water			INPUT: 5ng cDNA in 10.16uL water		
<input type="checkbox"/>			<input type="checkbox"/>		
QXT DILUTION: 1:20 dilution of QXT enzyme:Storage Soln., at least 2uL of dilution/sample			QXT DILUTION: 1:6 dilution of QXT enzyme:Storage Soln., at least 2uL of dilution/sample		
<input type="checkbox"/>			<input type="checkbox"/>		
MASTER MIX:			MASTER MIX:		
Reagent	1 rxn (uL)	____rxn	Reagent	1 rxn (uL)	____rxn
cDNA	10.16	N/A	cDNA	10.16	N/A
5x TD Buffer	8.84		5x TD Buffer	8.84	
QXTEnzyme	0.1		QXTEnzyme	0.33	
StorageSoln	1.9		StorageSoln	1.67	
Total:	21		Total:	21	

- 2** Place samples in thermocycler on "DNA Fragmentation" program:

Step	Temp(C)	Time(hh:mm:ss)
1	55	00:10:00
2	4	00:01:00
3	4	Hold

3 

Add  **5 µL** 0.2%SDS and  **24 µL** water/sample.

Reagent	1 rxn (uL)	____rxn (uL)
0.2%SDS	5	
H2O	24	
Totals:	29	

4  

1m

Seal, vortex, incubate at  **Room temperature** for  **00:01:00** .

B: Ampure bead purification 1 1m

5

A	B
Single cleanup	Double SPRI cleanup (experimental)
<input type="checkbox"/> Vortex beads and incubate @ RT for 30minutes	<input type="checkbox"/> Vortex beads and equilibrate @ RT for 30minutes
<input type="checkbox"/> Add 49uL bead mix/sample (1.0x)	<input type="checkbox"/> 0.5x cleanup (>600bp fragment removal): Add 25uL beads to each 50uL sample, vortex, briefly spin
<input type="checkbox"/> Seal samples, vortex, briefly spin down	<input type="checkbox"/> Incubate samples 5 mins @RT
<input type="checkbox"/> Incubate samples @RT 5 minutes	<input type="checkbox"/> Add samples to magnet rack, allow solution to clear, save supernatant to new clean tube (optional: save beads)
<input type="checkbox"/> Put samples on magnet rack, allow to clear, discard supernatant	<input type="checkbox"/> 0.2x cleanup (removal of <150bp fragments, rounding out total cleanup to 0.7x): Add 10uL bead mix/sample, seal, mix, incubate for 5 mins @RT
<input type="checkbox"/> Wash 1 with 200uL 80% EtOH	<input type="checkbox"/> Put samples on magnet rack, allow to clear, discard supernatant
<input type="checkbox"/> Wash 2 with 200uL 80% EtOH	<input type="checkbox"/> Wash 1 with 200uL 80% EtOH
<input type="checkbox"/> Remove last of the EtOH, allow samples to air dry for at least 10 minutes.	<input type="checkbox"/> Wash 2 with 200uL 80% EtOH
<input type="checkbox"/> Elute with 24uL H2O, incubate for 2 mins @RT	<input type="checkbox"/> Remove the last Ethanol and allow samples to air dry for at least 10 minutes
	<input type="checkbox"/> Elute with 24uL H2O, incubate for 2 mins @RT

C: Amplification

6



Add **24 µL** Kapa Hotstart Master Mix/sample.

7



Add **1 µL** i7 primer and **1 µL** i5 indexing primer, or **2 µL** from a multiplex plate.

NOTE THE USE OF CUSTOM PRIMERS.

8 

Seal, vortex, briefly spin down samples, put in thermocycler on “Pre-Capture PCR” protocol.

A	B	C	D
Step	Temp(C)	Time(hh:mm:ss)	
1	68	00:02:00	
2	98	00:02:00	
3	98	00:00:30	
4	57	00:00:30	
5	72	00:01:00	
6	Return to step 3	13x for 1ng cDNA inputs	11x for 5ng cDNA inputs
7	72	00:05:00	
8	4	hold	

Avoid S511 QXT i5 primer

D: Ampure bead purification 2



7m

9

Avoid S511 QXT i5 primer.

Make sure beads are equilibrated to  **Room temperature** .

10 

Add  **32.5 µL** -  **35 µL** Ampure beads/sample (0.65x-0.7x beads:sample ratio), seal, mix, brief spin.

11 

5m

Incubate sample/bead mix @  **Room temperature** for  **00:05:00** .

12 Move to magnet stand, allow solution to clear, discard supernatant.

13 




Wash 1 with  **200 µL** of 80% EtOH.

14 

Wash 2 with  **200 µL** of 80% EtOH.

15   

2m

Elute in  **20 µL** water. Add water to samples, take off magnet, seal, mix, spin, incubate @  **Room temperature** for  **00:02:00** .

16 Transfer samples back to magnet stand, allow to clear, save supernatant in separate, clean tubes.

17 REPEAT Ampure bead purification 2 for a total of 2x bead cleanups.

18 Evaluate on Tapestation.