

Jun 20, 2021

Full genome PCR amplification of all African Hepatitis B Virus genotypes

Forked from Full genome PCR amplification of all African Hepatitis B Virus genotypes

Louis S IV Le Clercq¹, Sheila Mary Bowyer¹, Simnikiwe Mayaphi¹

¹University of Pretoria

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bvykn7uw

Masters Project



Louis Stephane Le Clercq

University of the Free State, South African National Biodive...

ABSTRACT

This method is an adaptation of the full genome amplification protocol published by Gunther *et al.*, able to amplify most known genotypes. The primers bind near the fold in the genome where the DNA is partially double-stranded. Since genotype A1 of HBV frequently has two point mutations in the reverse primer binding region, a second primer was designed to circumvent them when amplifying A1 isolates. This protocol was used to amplify specimens from occult infections of HBV with genotypes A, D, and E.

ATTACHMENTS

Full-
genome_HBV_PCR_templa
te_UP-MedVir-Le-
Clercq.xlsx

DOI

dx.doi.org/10.17504/protocols.io.bvykn7uw

EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/bioproject/PRJNA737147>

PROTOCOL CITATION

Louis S IV Le Clercq, Sheila Mary Bowyer, Simnikiwe Mayaphi 2021. Full genome PCR amplification of all African Hepatitis B Virus genotypes. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bvykn7uw>

FORK NOTE

FORK FROM

Forked from Full genome PCR amplification of all African Hepatitis B Virus genotypes, Louis Stephane Le Clercq

KEYWORDS

HBV, Hepatitis, Full Genome, PCR, Genotypes

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 20, 2021

LAST MODIFIED

Jun 20, 2021

GUIDELINES

- Set up the mixes in a laminar flow cabinet [ⓘ On ice](#) .
- Amplicons can be stored at [ⓘ -20 °C](#)

MATERIALS TEXT

Reagents:

- [ⓘ Expand™ High Fidelity PCR](#)
- [System Roche Catalog #11732650001](#) [ⓘ High_fidelity_roche.pdf](#)
- [ⓘ Deoxynucleotide \(dNTP\) Solution Mix New England](#)
- [Biolabs Catalog #N0447S](#)
- [\[M\]25 Milimolar \(mM\) stock](#)
- [ⓘ Molecular grade water nuclease-free Contributed by users](#)
- Primers, [\[M\]15 Micromolar \(μM\) stock](#)

A	B	C	D
P1	1821--1841	5'- CTT TTT CAC CTC TGC CTA ATC A -3'	52.8
P2	1825--1806	5'- AAA AAG TTG CAT GGT GCT GG -3'	54.6
P2_A1	1825--1806	5'- AAA AAG TTG CAT GAT GAT GG -3'	49.3

Primers used to amplify full genomes of HBV. P1 is the forward primer and P2 or P2_A1 is the reverse primer, depending on the genotype to be amplified. Numbering is based on the EcoR1 site. The sequences and calculated T_m for each primer is indicated.

- [ⓘ SeaKem® LE](#)
- [Agarose Lonza Catalog #50004](#) Step 3.1
- [ⓘ UltraPure™ TBE Buffer, 10X Thermo](#)
- [Fisher Catalog #15581044](#) Step 3.1
- [ⓘ SYBR SAFE DNA stain Life](#)
- [Technologies Catalog #S33102](#) Step 3.1
- [ⓘ GeneRuler 1 kb DNA Ladder Thermo Fisher](#)
- [Scientific Catalog #SM0311](#) Step 3.2

Equipment:

- Thermal cycler
- Gel Documentation system

SAFETY WARNINGS

- Set up master mixes in a "DNA-free" room and laminar flow cabinet.
- Add DNA to reaction tubes in a "DNA-loading" laminar flow cabinet.
- Always dispose of biohazardous waste appropriately in accordance to lab regulations.
- Always wear gloves and a lab coat.
- Never directly look at the UV lamps.

BEFORE STARTING

- Thaw reagents/components **On ice**.
- Wipe workspace with **10 % volume** Bleach, followed by **70 % volume** Ethanol, and ddH₂O before (and after).
- UV the relevant laminar flow cabinets.

Master Mix preparation

1 Prepare the following two mixes **On ice** :

1.1 Master Mix 1:

For one **15 µl** reaction combine the following:

A	B	C	D
dNTPs	25 mM	200 µM	0.2 µL
P1	15 µM	300 nM	0.5 µL
P2 or P2_A1	15 µM	300 nM	0.5 µL
Expand Hi Fi Buffer	10x	1x	2 µL
MgCl ₂	15 mM	1.5 mM	-
ddH ₂ O	n.a.	n.a.	11.8 µL

Components for Master Mix 1 with their stock and final concentrations and volume needed for one 15 µL reaction.

1.2 Master Mix 2:

For one **5 µl** reaction combine the following:

A	B	C	D
Expand Hi Fi Buffer	10x	1x	0.5 µL
Expand Hi Fi Enzyme	3.5 U/µL	2.6 U	0.75 µL
ddH ₂ O	n.a.	n.a.	3.75 µL

Components for Master Mix 2 with their stock and final concentrations and volume needed for one 5 µL reaction.

1.3

- Add **15 µl Master Mix 1** to **5 µl template DNA** for a **20 µl** reaction in thin walled PCR tubes.
- At the Thermal cycler, add an additional **5 µl Master Mix 2** to each tube after initial denaturation.

Thermal cycling

2 Perform thermal cycling according to the following conditions:

- Initial denaturation at **94 °C** for 2 minutes.
- Cooling to **58 °C** before adding Mix 2.
- 40 cycles of:
 1. Denaturation at **94 °C** for 40 seconds
 2. Annealing at **55 °C** for 90 seconds
 3. Elongation at **68 °C** for 180 seconds*
- *Add 180 seconds every 10 cycles
- Cooling/hold at **4 °C**

Electrophoresis 50m

3 Perform TBE-gel electrophoresis to confirm success of amplification prior to amplicon clean-up.

3.1

Prepare a **0.7 % (v/v)** gel with [SeaKem® LE](#) and [Agarose Lonza Catalog #50004](#) and [UltraPure™ TBE Buffer, 10X Thermo Fisher Catalog #15581044](#), pre-stained with [SYBR SAFE DNA stain Life Technologies Catalog #S33102](#). [SYBR SAFE DNA stain Life Technologies Catalog #S33102](#) is usually added at a concentration of **1 µl Stock** per **10 mL TBE**.

3.2

Load **5 µl** of PCR product to the gel using **2 µl loading buffer** such as [Gel Loading Dye Orange \(6X\) - 4.0 ml New England Biolabs Catalog #B7022S](#) and molecular weight marker, e.g. [GeneRuler 1 kb DNA Ladder Thermo Fisher Scientific Catalog #SM0311](#), and run at **100 Volt** for **00:50:00**.

3.3


Visualize and capture gel on an appropriate imager and paired software, e.g. the

Gel Doc XR+ Gel Documentation System
Gel Documentation System
Bio-rad Laboratories 1708195 [Link](#)

3.4

Expected result:





Gel image showing single bands for individually amplified full genome HBV isolates, flanked by a molecular weight marker.