

Sep 17, 2024

Small Scale IPTG Induction | Protein Expression

DOI

dx.doi.org/10.17504/protocols.io.yxmvme47og3p/v1

Catherine M Gohar¹

¹University of Washington



Catherine M Gohar

University of Washington

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yxmvme47og3p/v1

Protocol Citation: Catherine M Gohar 2024. Small Scale IPTG Induction | Protein Expression. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yxmvme47og3p/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

Test

Created: September 17, 2024

Last Modified: September 17, 2024

Protocol Integer ID: 107736

Disclaimer

Test

Abstract

Test



Before start

It's best to start small-scale protein purification on a Monday or Tuesday. This protocol will take at least 3 days before the first freezing step.



Seed Culture

1d

- 1 Preheat the shaking tube incubator to 37°C
- 2 Label a 15 mL culture tube with the vector strain you will be using, and fill with 5 mL of LB broth with 1x Kanamycin antibiotic or 5ul of 1000x.
- 3 Place the labeled tube into the shaking incubator to warm, then grab an ice bucket and fill with ice.
- 4 Take the BL21 glycerol stock from the -80°C freezer and place immediately into the ice bucket.
- 5 Use a $200\text{ }\mu\text{L}$ pipette tip to scrape some of the ice from the tube, and place the tip directly into the culture tube with LB.
- 6 Let the culture shake at 37°C overnight to be completely saturated by the next morning.

12h

Culture Expansion








2h

- 7 Grab a sterilized 250 mL Erlenmeyer flask and fill with 45mL of Terrific Broth and 1x Kanamycin antibiotic.
- 8 Prewarm the flask in the shaking tube incubator at 37°C
- 9 Add all 5 mL of the seed culture to the Erlenmeyer flask. The culture should expand 1:10 in 1-2 hours, to an OD around 0.4-0.6.
- 10 Thaw a stock of IPTG on ice about 30 minutes before you suspect the culture will be at 0.4 OD.


1h 30m

30m



- 11 Quantify the OD of your  50 mL flask using  500 μ L of culture in a spectrometer cuvette. Be sure to blank with  500 μ L of Terrific Broth.
- 12 Once OD hits 0.4-0.6, preheat the fridge shaking incubator to  16 °C
- 13 Dilute  1 Molarity (M) IPTG stock to the appropriate concentration for protein induction by adding to expanded culture. This will most likely be a range of  100 micromolar (μ M) to  1 millimolar (mM)
- 13.1 To calculate IPTG concentration, use $V1 \times C1 = V2 \times C2$, where column E is equivalent to V2.





A	B	C	D	E
50mL	1mM	?	1000mM	50uL
50mL	0.8mM	?	1000mM	40uL
50mL	0.6mM	?	1000mM	30uL
50mL	0.4mM	?	1000mM	20uL
50mL	0.2mM	?	1000mM	10uL

- 14 Take your IPTG-induced culture and shake overnight at  16 °C for about 16-18 hours.

18h

Pellet Cells

5m

- 15 Take each  50 mL culture and pour into a  50 mL conical tube, labelled.
- 16 Take the cells to a large or spinning-bucket centrifuge and spin at  500 x g, 4°C, 00:05:00
- 17 Decant cells into a bleach waste bin labeled for culture, and dry on a paper towel for 5-10 seconds
- 18 Freeze the pellet at  -80 °C by storing in a freezer for storage, or continue to the protocol for cell lysis.

5m

10s

