



3 ▼

Jun 10, 2022

# 🌐 Haemolymph extraction of adult *Drosophila* V.3

Björn Brembs<sup>1</sup>, Christine Damrau<sup>2</sup>, Brembs' lab members<sup>2</sup>

<sup>1</sup>Universität Regensburg; <sup>2</sup>Björn Brembs lab, Universität Regensburg

1



[dx.doi.org/10.17504/protocols.io.eq2lyd5elx9k/v3](https://dx.doi.org/10.17504/protocols.io.eq2lyd5elx9k/v3)

**Brembs Lab**



**Björn Brembs**  
Universität Regensburg

This is a very simple protocol showing how to extract haemolymph from adult *Drosophila melanogaster*. (Based on protocols from Sigma Aldrich).

DOI

[dx.doi.org/10.17504/protocols.io.eq2lyd5elx9k/v3](https://dx.doi.org/10.17504/protocols.io.eq2lyd5elx9k/v3)

<https://doi.org/10.3389/fnsys.2017.00100>

Björn Brembs, Christine Damrau, Brembs' lab members 2022. Haemolymph extraction of adult *Drosophila*. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.eq2lyd5elx9k/v3>  
Björn Brembs



protocol

Damrau C, Toshima N, Tanimura T, Brembs B, Colomb J, Octopamine and Tyramine Contribute Separately to the Counter-Regulatory Response to Sugar Deficit in . *Frontiers in Systems Neuroscience* doi: [10.3389/fnsys.2017.00100](https://doi.org/10.3389/fnsys.2017.00100)

minor corrections and clarifications to the protocol's wording and description

\_\_\_\_\_ protocol ,

Jun 10, 2022

Jun 10, 2022


64321

## Fly preparation

- 1 Punch three holes in a 0.5ml Eppendorf cap and put it into a 1.5ml Eppendorf cap with removed lid.
- 2 Remove the flies' wings, spear the fly's thorax with a needle.
- 3 Collect 20 speared flies in the 0.5ml Eppendorf cap with holes, on ice.
- 4 Centrifuge the 0.5ml Eppendorf cap within the 1.5ml one (1 min, 5000 rpm, at 4°C).  
🕒 00:01:00
- 5 Discard the 0.5ml Eppendorf cap, collect the extracted hemolymph with a fine capillary.
- 6 Record the amount of haemolymph (to fill up 0.5µl you need around 50 flies).

## Enzymatic procedure

- 7 Add 19.5µl cold PBS to 0.5µl haemolymph.
- 8 Add 10µl of this mixture to 30µl Citrate Acid Buffer and 10µl of a 3% Trehalase-Citrate acid buffer solution.
- 9 Incubate over night at 37°C.  
🕒 18:00:00
- 10 Add 50µl Tris Buffer.

- 11 80µl of this mixture are added to 156.8µ Glucose oxidase (aliquot in the freezer) and 3.2µl o-Dianisidine (freshly added from the fridge).
- 12 Incubate for exactly 30 min at 37°C.  
 00:30:00
- 13 Stop reaction by adding 160µl Sulfuric Acid.
- 14 Measure at 540nm at the (nanoDrop) spectrometer.