

Oct 24, 2020

NLP screening (dot blot)

Andreea S¹¹University of Groningen

Other

dx.doi.org/10.17504/protocols.io.bkbfsjn

iGEM Groningen 2020

a.stan.6

ABSTRACT

In order to determine if the His-tagged NLP is expressed and secreted by the bacteria, we plan to use the dot blot procedure with anti-His tag anti-bodies. Dot blot is an immunological technique used for detecting proteins directly from the culture supernatant (without gel separation). Thus, the samples are directly spotted on the membrane, making it a high throughput procedure ideal for testing different secretion signals.

DOI

dx.doi.org/10.17504/protocols.io.bkbfsjn

PROTOCOL CITATION

Andreea S 2020. NLP screening (dot blot). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bkbfsjn>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 25, 2020

LAST MODIFIED

Oct 24, 2020

PROTOCOL INTEGER ID

41031

MATERIALS TEXT




MATERIALS

[TBS: 20 mM Tris-HCl 150 mM NaCl pH 7.5](#) Contributed by users[TBS-T: 0.05% Tween20 in TBS](#) Contributed by users[BSA/TBS-T: 0.1% BSA in TBS-T](#) Contributed by users

ABSTRACT

In order to determine if the His-tagged NLP is expressed and secreted by the bacteria, we plan to use the dot blot procedure with anti-His tag anti-bodies. Dot blot is an immunological technique used for detecting proteins directly from the culture supernatant (without gel separation). Thus, the samples are directly spotted on the membrane, making it a high throughput procedure ideal for testing different secretion signals.

- 1 On nitrocellulose membrane indicate the blotting region by drawing a grid by pencil.

- 2 Slowly apply  2 µl of sample on the nitrocellulose membrane at the center of the grid.
- 3 Use purified NLP-His and apply 5 points of increasing concentration of them on the nitrocellulose membrane in order to create a standard concentration curve.
- 4 Dry the membrane.
- 5 Block non-specific binding sites by washing the membrane in 5% BSA in TBS-T for 1 h at room temperature.
- 6 Incubate with primary antibody diluted at the concentration recommended by the producer for 30 min. This step may require optimizing of the concentration.
- 7 Wash 3 times with TBS-T for  00:05:00
- 8 Incubate with the secondary antibody conjugated with HRP for  00:30:00 (concentration recommended by the producer).
- 9 Wash 3 times with TBS-T (15 min, 5 min, 5 min)
- 10 Wash with TBS
- 11 Incubate with ECL reagent for 1 min. Develop the blot.
- 12 Compare the intensity of the sample with the standard curve in order to estimate the concentration.