

Aug 28, 2024

## Fear Conditioning

DOI

**dx.doi.org/10.17504/protocols.io.261ge56bdg47/v1**

Sabina Marciano<sup>1,2</sup>, Roberta Marongiu<sup>1,2</sup>

<sup>1</sup>Department of Neurological Surgery, Weill Cornell Medical College, New York, NY 10065;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA

ASAP Collaborative Rese...



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.261ge56bdg47/v1](https://dx.doi.org/10.17504/protocols.io.261ge56bdg47/v1)

**Protocol Citation:** Sabina Marciano, Roberta Marongiu 2024. Fear Conditioning. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261ge56bdg47/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 04, 2024

**Last Modified:** August 28, 2024

**Protocol Integer ID:** 97806

**Keywords:** ASAPCRN, motor impairment, mouse, behavior

**Funders Acknowledgement:**

**Aligning Science Across**

**Parkinson's**

**Grant ID:** 020608



## Abstract

Behavioral assay to measure fear learning and memory in mice.

## Materials

- MedAssociates fear conditioning chambers
- white floor inserts
- black ceiling inserts
- peppermint and limonene spray



## Fear Acquisition

- 1 Habituate mice to the room. Set up fear conditioning chambers as well as protocols (listed below) in VideoFreeze software.
- 2 Place individual mice in a lit fear conditioning apparatus scented with peppermint spray.
- 3 Protocol:

- 00:02:00 baseline followed by 5 tone-shock pairings where the shock follows the tone
  - shock = 0.5mA, 00:00:01 is delivered after the end of the tone (75dB, 2kHz, 00:00:19 )
  - empty trace intervals between the shocks last 00:01:00

3m 20s
- 4 After the protocol runs for each mouse, place the mouse back in the home cage and clean the fear conditioning chamber thoroughly.
- 5 Contextual test is performed 24:00:00 after the fear acquisition test. Cued Test is performed at least 02:00:00 after the contextual test.

1d 2h

## Contextual Test

- 6 Goal: Measure the subjects' fear response to the environment where fear conditioning was received
- 7 24 hours after fear acquisition: Mice are placed in a lightened FC apparatus scented with peppermint spray.
- 8 Protocol: Neither cue nor shock are delivered during this test session (Duration 00:07:30 )

7m 30s

## Cued Test


- 9 Goal: Measure the subjects fear response to the CS, and ideally no other features of the training environment. The key is to alter the chamber so the subjects are exposed to a new set of environmental cues.



- 10 Mice are placed in a lightened FC apparatus which is altered by placing a smooth floor insert on shock grid and black triangles on top.

Limonene spray is used, instead of peppermint scent.

- 11 Protocol:

- Baseline: 2 min followed by five tone-cues
- No shock is delivered after the end of the tone (75dB, 2 kHz,  00:00:19 )

19s

## Primary Analyses

- 12 Fear conditioning acquisition

- % Component freezing for baseline, each tone (if used), shock, and interstimulus interval
- Graphs for each component should be graphed separately
- Average motion

### Context

- % Component freezing per minute and total
- Per minute and total will be on separate graphs
- Average motion per minute and total

### Cued

- %Component freezing for baseline, each tone, and interstimulus interval
- Graphs for each component should be graphed separately
- Average motion for baseline, each tone, and interstimulus interval