



Apr 14, 2020

The role of multiprotein bridging factor 1 (MBF1c) in Funariaceae species under salt stress

Maria Victória Magalhães de Vargas^{1,2}, Guilherme Afonso Kessler de Andrade^{1,2}, Sara Navarrete Bohi Goulart^{1,2}, Bruna Mota Bernardes^{1,2}, Filipe Victoria^{1,2}

¹Núcleo de Estudos da Vegetação Antártica (NEVA), Universidade Federal do Pampa (UNIPAMPA), São Gabriel, RS, Brazil, ²neva@unipampa.edu.br

Other dx.doi.org/10.17504/protocols.io.be3sjgne



ABSTRACT

Bryophytes are one of the most successful plants in nature and can tolerate several environmental stresses, such as saline stress. MBF1 is a transcription factor that may be involved in reprogramming transcription, which in turn may confer tolerance to stress conditions. The objective of this study was to quantify PaMBF1c under saline stress in three species of moss: *Funaria hygrometrica*, *Physcomitrium acutifolium*, and *Physcomitrella patens*. All moss species underwent 12 h and 24 h treatment with salt at different concentrations, followed by qRT-PCR experiments.

In *F. hygrometrica*, the MBF1c gene was up-regulated when subjected to 12 h salt treatment at a medium concentration, and may be an important factor for the tolerance of plants to NaCl concentrations. In *P. acutifolium*, 24 h NaCl treatment at of 400, 600, and 800 mM resulted in a gradual up-regulation of MBF1c expression, reflective of the increasing NaCl concentration. In *Physcomitrella*, the expression profile of MBF1c showed a down-regulation when treated with 200, 400, and 800 mM NaCl, suggesting that MBF1c may not interfere with the saline stress response in this plant.

The diversity of these expression profiles demonstrate that the MBF1c response is not universal neither among plants nor between species of the same family, as shown in the mosses studied.

ATTACHMENTS

[Methods_Protocol.docx](#)

MATERIALS

NAME	CATALOG #	VENDOR
Qubit™ dsDNA HS Assay Kit	Q32851	Invitrogen - Thermo Fisher
GoTaq(R) Hot Start Polymerase, 2500u	M5006	Promega
Qubit RNA HS Assay Kit	Q32852	Thermo Fisher Scientific
PureLink™ RNA Mini Kit	12183020	Thermo Fisher Scientific
High-Capacity RNA-to-cDNA™ Kit	4387406	Thermo Fisher



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited