

Jun 24, 2024

# BAF\_Protocol\_013\_Lipidomics: Database Search MS-DIAL and Analysis using Metaboanalyst 6.0

DOI

**[dx.doi.org/10.17504/protocols.io.ewov192jklr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov192jklr2/v1)**

Nicholas Sherman<sup>1</sup>

<sup>1</sup>University of Virginia Biomolecular Analysis Facility Core

Biomolecular Analysis Fac...



Dilza Silva

UVA

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.ewov192jklr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov192jklr2/v1)**

**Protocol Citation:** Nicholas Sherman 2024. BAF\_Protocol\_013\_Lipidomics: Database Search MS-DIAL and Analysis using Metaboanalyst 6.0. protocols.io **<https://dx.doi.org/10.17504/protocols.io.ewov192jklr2/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 17, 2024

**Last Modified:** June 24, 2024

**Protocol Integer ID:** 101969

## Abstract

These steps represent a starting point for analysis of lipidomics data. Other data analysis and/or packages may be used.

## Materials

Raw files.

MS-DIAL Software: <https://systemsomicslab.github.io/compms/msdial/main.html>

Databases: <https://systemsomicslab.github.io/compms/msdial/main.html#MSP>

MetaboAnalyst 6.0 online tool: <https://www.metaboanalyst.ca/MetaboAnalyst/home.xhtml>

## MS-DIAL V 4.9.221218 starting

- 1 Thermo RAW files are organized in two folders - one for all positive mode raw files and another for all negative mode raw files. These files will be loaded into MS-DIAL in two separate searches. Download MSP library files for Positive and Negative mode and save them in a folder in the computer you are using to perform the searches.
- 2 Open MS-DIAL
- 3 Click File -> New project
- 4 Project file path: Browse the folder with Positive or Negative raw files  
Ionization type: Soft ionization (LC-MSMS)  
Separation type: Chromatography (LC)  
MS method type: Conventional LC/MS or data dependent MS/MS  
Data type: MS1 Profile, MS/MS Centroid  
Ion mode: Positive or Negative  
Target omics: Lipidomics  
Click "Next"
- 5 Analysis file path: browse to the folder with raw files, select all the apply for the search. Should include at least 3 blanks and 3 QC runs.  
Choose the file type for each sample (Blank, QC, sample)  
Class ID - add short group names for each replicate without adding numbers among replicates.

## Analysis parameter setting - keep default if not mentioned in this section

### 6 Data collection tab

Mass accuracy: MS1 tolerance 0.01 Da,

MS2 tolerance 0.025 Da

- Advanced:

Data collection parameters: retention time begin 0 min, retention time end 33 min, MS1 range begin 150 Da, MS1 range end 2000 Da, MSMS range begin 150 Da, MSMS range end 2000 Da, Isotope recognition: max charge number 2.

### 7 Peak Detection tab

Peak detection parameters: minimum peak height 10000 amplitude, mass slice 0.1 Da

- Advanced:

Smoothing method: linear weighted moving average

Smoothing level: 3 scan



Minimum peak width: 5 scan

## 8 **MS2Dec tab**

Deconvolution parameters: sigma window value 0.5,

MS2 abundance cut off 10 amplitude.

- Advanced:

Check Exclude after precursor ion Keep the isotope ions until 0.5 Da

Check Keep the isotopic ions w/o MS2Dec

## 9 **Identification tab**

MSP file and MS/MS identification setting: the library will be loaded automatically

Retention time tolerance: 100 min,

Accurate mass tolerance (MS1) 0.01 Da,

Accurate mass tolerance MS/MS 0.05,

Identification score cut off 80%.

## 10 **Adduct tab**

For positive mode acquisition:

Check [M+H]<sup>+</sup>

Check [M+NH<sub>4</sub>]<sup>+</sup>

For negative mode acquisition:

Check [M-H]<sup>-</sup>

Check [M-H<sub>2</sub>O-H]<sup>-</sup>

## 11 **Alignment tab**

Alignment parameters settings: Under result name add a name for the alignment file, reference file choose one of the QC runs, retention time tolerance 0.05 min, MS1 tolerance 0.015 Da

- Advanced:

Check remove feature based on blank information

Check keep "reference matched" metabolite features

Check keep "suggested (w/o MS2)" metabolite features

Check keep removable features and assign the tag

Check Gap filling by compulsion

## 12 **At the bottom**

Check together with Alignment

Click "Finish" --> processing will start.

Once the search is finished, click on the name of the file that is showing at the Alignment Navigator

Click Export - Alignment result, choose a folder to save the file, check remove features from blank, and export.

Data can be loaded into Excel/R or any other software environment for further data analysis.

**13 Data Analysis:**

load data into Excel, filter out: m/z match FALSE and identifications named w/o MS2.

Combine positive and negative identifications

Summarize the data by lipid name (ex: using R), save as txt file.

Load the summarized into R or the online MetaboAnalyst 6.0 tool.

**14**

**MetaboAnalyst 6.0** --> statistical analysis

Upload .txt data as peak intensities (check file format), no data filtering,

Sample normalization: normalization by median,

Data transformation: Log transformation (base 10), no Data Scaling.

At this point there are distinct analyses that can be done, our core provides the following as a first-step analysis:

Chemometrics analysis: PCA, PLS-DA with importance measures using VIP.

Univariate analysis: fold change and t-test -> volcano plot (two groups comparisons).

Clustering analysis: heatmap.

## Protocol references

Tsugawa H, Cajka T, Kind T, Ma Y, Higgins B, Ikeda K, Kanazawa M, VanderGheynst J, Fiehn O, Arita M. MS-DIAL: data-independent MS/MS deconvolution for comprehensive metabolome analysis. *Nat Methods*. 2015 Jun;12(6):523-6. doi: 10.1038/nmeth.3393. Epub 2015 May 4. PMID: 25938372; PMCID: PMC4449330.

Tsugawa, H., Ikeda, K., Takahashi, M. *et al.* A lipidome atlas in MS-DIAL 4. *Nat Biotechnol* **38**, 1159–1163 (2020). <https://doi.org/10.1038/s41587-020-0531-2>

Pang, Z., Lu, Y., Zhou, G., Hui, F., Xu, L., Viau, C., Spigelman, A., MacDonald, P., Wishart, D., Li, S., and Xia, J. (2024) **MetaboAnalyst 6.0: towards a unified platform for metabolomics data processing, analysis and interpretation** *Nucleic Acids Research* (doi: 10.1093/nar/gkae253)