

Sep 09, 2024 Version 1

## Electrolysis agarose gel preparation V.1

DOI

**[dx.doi.org/10.17504/protocols.io.kqdg32d4qv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg32d4qv25/v1)**

UM ILT team<sup>1</sup>, Mar Roca Cugat<sup>1</sup>

<sup>1</sup>Maastricht University

UM ILT team: Original authors

Mar Roca Cugat: Importing into Protocols.io

UM-BMS Student Team



**Mar Roca Cugat**

Maastricht University, Universitat Autònoma de Barcelona, In...

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.kqdg32d4qv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg32d4qv25/v1)**

**Protocol Citation:** UM ILT team, Mar Roca Cugat 2024. Electrolysis agarose gel preparation. **protocols.io**  
**<https://dx.doi.org/10.17504/protocols.io.kqdg32d4qv25/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 09, 2024

**Last Modified:** September 09, 2024

**Protocol Integer ID:** 107168



## Abstract

This is how students at Maastricht University's Biomedical Sciences bachelor are instructed to prepare agarose gel for electrolysis.



## Materials

- Micropipettes + tips
- Glass erlenmeyer
- Microwave
- Gel casting module
- UV-transparent gel casting tray
- Combs
- Electrophoresis system
- Power Source
- UV visualization unit

## Reagents

- 50x TAE buffer (Tris-EDTA-Acetic Acid buffer)
- Aragose
- GreenStar /!\
- DNA marker
- DNA gel loading dye



- 1 Prepare the casting tray by putting the sides together.
- 2 The silicon coated site should point towards the inside.
- 3 Apply the elastic band around the extruding parts in the middle of the sides to hold everything in place.
- 4 Flip the assembled tray in the correct position and press it lightly on the table to ensure proper alignment of the both sides.
- 5 Add the comb and put the tray in the box to catch any potential spills if not correctly assembled.
- 6 Prepare a 1% agarose gel in 1x TAE buffer. The total volume depends on the size of the casting tray; we will make  100 mL gels. Calculate and weigh the amount of agarose required for your volume. Deposit the agarose in a glass Erlenmeyer flask and add the indicated volume of 1x TAE buffer.
- 7 Pour the mixture in the casting tray with combs. Remove any air bubbles with a plastic filter tip.
- 8 Allow the gel to solidify (  00:10:00 ).
- 8.1 Remove the comb.
- 9 Insert the casting tray into the electrophoresis equipment; the side of the slots in the direction of the cathode (black).  
*Do not add the metal casting system into the electrophoresis system, just the gel.*
- 10 Prepare and add 1x TAE Buffer to the system until the gel is completely submerged.

10m