



Version 2

Apr 23, 2021

Initial Rapid Pathology Assessment of Kidney Tissue V.2

Jamie Allen¹, Carrie Romer¹, Elizabeth Neumann¹, Maya Brewer¹, Haichun Yang¹, Jeff Spraggins¹, Danielle Gutierrez¹

¹Vanderbilt University

1

Works for me

dx.doi.org/10.17504/protocols.io.bt8gnrtw

VU Biomolecular Multimodal Imaging Center

Tech. support email: jeff.spraggins@vanderbilt.edu

Jamie Allen

Vanderbilt University

ABSTRACT

Scope:

Prepare formalin fixed tissue with freezing fresh tissue protocol (dx.doi.org/10.17504/protocols.io.6wghfbw).

This protocol provides the steps to prepare formalin fixed human kidney tissue and perform histology assessment for normalcy.

Expected Outcome:

Pathological assessment of kidney tissue for use in LC-MS/MS and imaging.

DOI

dx.doi.org/10.17504/protocols.io.bt8gnrtw

PROTOCOL CITATION

Jamie Allen, Carrie Romer, Elizabeth Neumann, Maya Brewer, Haichun Yang, Jeff Spraggins, Danielle Gutierrez 2021. Initial Rapid Pathology Assessment of Kidney Tissue. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bt8gnrtw>

Version created by Jamie Allen

KEYWORDS

HuBMAP, Kidney, Quality Assessment, PAS

LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Apr 14, 2021

LAST MODIFIED

Apr 23, 2021

PROTOCOL INTEGER ID

49128

1 Tissue is embedded in protocol: dx.doi.org/10.17504/protocols.io.br4fm8tn

- 2 Section samples at 5 μ m on a microtome.
- 3 PAS stain tissue sections with protocol: dx.doi.org/10.17504/protocols.io.buaknscw
- 4 Scan slides with brightfield scanner (Leica) and save as .tiff or .jpg
- 5 Place saved images on QuPath for analysis.
- 6 Assess and record the following information for each tissue:
(%): cortex and medulla
(Yes/No): Pyramid Presence, Autolysis, Non-Renal disease (*i.e.* cancer)
(0=None – 3=Severe): Glomerular disease, Tubulointerstitial disease
- 7 Based on FFPE assessment, kidneys will be used accordingly:

LC-MS/MS normalcy:
low or no glomerular disease, and 75% or more cortex for normal
100% tumor for diseased

3D Imaging:
50:50 cortex:medulla
low or no glomerular disease
no renal disease