

•



Aug 23, 2020

## © ELISA for quantification of human immunoglobulin G (IgG) in serum or plasma.

Angel A Justiz-Vaillant<sup>1</sup>, Belkis Ferrer-Cosme<sup>2</sup>

<sup>1</sup>University of the West Indies St. Augustine; <sup>2</sup>"Saturnino Lora Torres' Provincial Teaching Clinical Surgical Hospital. Cuba

1 Works for me dx.doi.org/10.17504/protocols.io.bj7skrne

University of the West Indies angel.vaillant@sta.uwi.edu

Angel Justiz-Vaillant University of the West Indies St. Augustine

## ABSTRACT

IgG is a monomer with an approximate molecular weight of 146 Kd and a serum concentration of 9.0 mg/mL. It is synthesized mostly in the secondary immune response to pathogens. IgG can activate the classical pathway of the complement system, and it also is highly protective. The four subclasses of IgG include IgG1, IgG2, IgG3, and IgG4. IgG crosses the placentae, protecting the neonate from infectious diseases. [1]

DOI

dx.doi.org/10.17504/protocols.io.bj7skrne

## PROTOCOL CITATION

Angel A Justiz-Vaillant, Belkis Ferrer-Cosme 2020. ELISA for quantification of human immunoglobulin G (IgG) in serum or plasma.. **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.bj7skrne

## LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 23, 2020

LAST MODIFIED

Aug 23, 2020

PROTOCOL INTEGER ID

40914

- 1 An anti-human IgG coating antibody is adsorbed onto the microwells by incubation overnight at 4°C with carbonatebicarbonate buffer.
- 2 Add 50 μl of human serum or plasma. Human IgG present in the serum or plasma binds to antibodies adsorbed into the microwells.
- 3 The microplate is blocked with 3% non-fat milk-PBS buffer and later wash to remove unbound proteins.

4	Fifty (50) µl of biotin-conjugated anti-lgG antibody is added. The optimal dilution must be investigated.
5	The microplate is rewashed with PBS-Tween 20 buffer, pH 7.4.
6	One hundred $\mu l$ of streptavidin-HRP conjugate is added and it binds to the biotin-conjugated anti-lgG antibody. The optimal dilution of this conjugate must be investigated.
7	The plate is washed following incubation to remove the unbound Streptavidin-HRP.
8	Add 100 μl of 3,3',5,5'- tetramethylbenzidine (TMB; Sigma-Aldrich) into each well.
9	Incubate the microwells in the dark for 20 min.
10	A colored product is formed in proportion to the quantity of IgG present in the sample or standard.
11	The reaction is terminated by addition of 100 $\mu$ l 3M H2SO4 and the absorbance is measured at 450 nm.
12	A standard curve is made from 7 human IgG standard dilutions and the human IgG sample concentration is determined.
13	For better results place the microplate on a microplate shaker in every incubation.