



Version 2 ▼

Sep 13, 2022

Rat brain processing for histological analyses (update) V.2

miquel.vila¹¹Vall d'Hebron Research Institute

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.e6nvwj61wlmk/v2

joan.compte

ABSTRACT

Protocol for rat brain processing in order to perform histological analyses

DOI

dx.doi.org/10.17504/protocols.io.e6nvwj61wlmk/v2

PROTOCOL CITATION

miquel.vila 2022. Rat brain processing for histological analyses (update).
protocols.io
<https://protocols.io/view/rat-brain-processing-for-histological-analyses-upd-cgiatuae>
Version created by [marta.gonzalez.sepulveda](#)

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 13, 2022

LAST MODIFIED

Sep 13, 2022

PROTOCOL INTEGER ID

69922

Rat perfusion 25m

- 1 Deeply anesthetize animals with sodium pentobarbital (50 mg/kg, i.p.)

2 Perfuse through the left ventricle with saline [0.9% (wt/vol)] at room temperature (RT) 5m

3 Perfuse again with ice-cold formaldehyde solution 4% in PBS buffered for histology 15m

Postfixation 1d

4 Remove brains and post-fix them for 24 h in the same fixative 1d

Processing for microtome sectioning

5 Wash twice with 0.1 M PBS and process for paraffin embedding following standard procedures (performed by an external facility)

Alternatively, brains could be stored in 0.1 M PBS at 4 °C (not for over a month without changing the PBS) prior paraffin processing

Processing for cryostat sectioning 2d

6 Cryoprotect for 24-48 h (until they sink) in 30% sucrose at 4 °C 2d

7 Exchange sucrose for 0.1 M PBS

8 Immerse brains in cold (-30°C) 2-methylbutane for 30 s and store at -80 °C 30s

9 Include in OCT

Sectioning

10 Perform sectioning with a sliding microtome at 5-µm-thickness for paraffin samples or in a cryostat at 20- or 30-µm-thickness for frozen samples.

