



VERSION 4

MAY 09, 2023

OPEN  ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.3byl4b7rjvo5/v4

Protocol Citation: sandra.crameri 2023. TS Spurrs - unstained FFPE on Thermanox coverslips (TM - 013). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4b7rjvo5/v4> Version created by [sandra.crameri](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: May 09, 2023


Last Modified: May 09, 2023

PROTOCOL integer ID:
81628

TS Spurrs - unstained FFPE on Thermanox coverslips (TM - 013) V.4

sandra.crameri¹

¹CSIRO ACDP



sandra.crameri

ABSTRACT

Processing unstained FFPE sections on Thermanox coverslips to Spurrs resin block for EM.

GUIDELINES

All step times are the minimum time required. Longer steps are acceptable.

HEADER

- 1 SAN:







SPEC No:

OPERATOR & STEPS:

OPERATOR & STEPS:


Dewax FFPE sections

15m

- 2  xylene Contributed by users or X3B  00:05:00 5m
- 3  xylene Contributed by users or X3B  00:05:00 5m
- 4  xylene Contributed by users or X3B  00:05:00 5m


Rehydrate

15m

- 5 [M] 100 % volume  Ethanol Contributed by users  00:05:00 5m
- 6 [M] 95 % volume  Ethanol Contributed by users  00:05:00 5m

7

[M] 70 % volume

 Ethanol Contributed by users 00:05:00


5m

Conventional - with reduced times for 5um thick FFPE Secti...

2h 50m

8


[M] 2.5% % volume

 Glutaraldehyde aqueous solution 25% EM Grade ampoules ProSciTech Catalog #C002


in

[M] 0.1 Molarity (M)

Sorenson's phosphate buffer

 07.2


300 mosmol/kg

 00:20:00

20m

9


Wash buffer

 00:05:00


5m

10

[M] 1 Mass / % volume

 Osmium
Tetroxide ProSciTech


in buffer

 00:30:00

30m

11


[M] 70 % volume

 Ethanol Contributed by users 00:10:00

10m

12

[M] 95 % volume

 ethanol Contributed by users 00:05:00

5m

13

[M] 100 % volume

 ethanol Contributed by users 00:05:00

5m

- 14 [M] 100 % volume dry  ethanol Contributed by users  00:05:00 5m
- 15 [M] 100 % volume Spurrs: [M] 100 % volume dry  ethanol Contributed by users in 1:1 ratio  00:30:00 30m
- 16 [M] 100 % volume Spurrs  00:30:00 30m
- 17 [M] 100 % volume Spurrs  00:30:00 30m
- 18 Oven polymerise  65 °C  Overnight