



DEC 17, 2022

WORKS FOR ME

1

Phage amplification and concentration

In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.yxmvmnb86g3p/v1Adair Borges¹, Januka Athukoralage¹¹Arcadia Science

Arcadia Science



Arcadia Science

COMMENTS 0

ABSTRACT

This protocol details methods to amplify bacteriophages T4 and SP01, and concentrate using either PEG precipitation or filtration.

DOI

dx.doi.org/10.17504/protocols.io.yxmvmnb86g3p/v1

PROTOCOL CITATION

Adair Borges, Januka Athukoralage 2022. Phage amplification and concentration. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yxmvmnb86g3p/v1>

COLLECTIONS ⓘ

[Protocol collection: Phage DNA isolation and chemical analysis](#)

KEYWORDS

phage, gDNA, nucleoside, genome modification, HPLC, Mass spectrometry

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 21, 2022

LAST MODIFIED

Dec 17, 2022

OWNERSHIP HISTORY

Jun 21, 2022



Adair Borges

Sep 20, 2022



Arcadia Science

Sep 20, 2022



Arcadia Science

PROTOCOL INTEGER ID

65060

PARENT PROTOCOLS

Part of collection

[Protocol collection: Phage DNA isolation and chemical analysis](#)

Phage amplification

- 1 For each host, prepare a 30 mL subculture at 1:100 dilution (300 μ L in 30 mL). Grow for 2 h at 37 °C in a shaking incubator. Subculture media (LB or ENB) should be supplemented with 1 mM CaCl_2 and 1 mM MgCl_2 .
- 2 Grow a 3 mL overnight culture of phage host at 37 °C with shaking. *Escherichia coli* strain B (ATCC Strain 11303) should be propagated in Lysogeny Broth (ATCC Medium 1065), and *Bacillus subtilis* strain 168M (ATCC Strain 27370) is propagated in Enriched Nutrient broth (ATCC Medium 265).

Note

Note

Note

- 3 Add $\sim 10^7$ PFU of phage T4 (ATCC Strain 11303-b4) to the *E. coli* subculture, and $\sim 10^7$ PFU of phage SP01 (ATCC Strain 27370-b1) to the *B. subtilis* subculture. Return the infected subcultures to the incubator, at 37 °C with shaking. Grow until the culture is cleared (around 5 h). If the culture doesn't completely clear within 5 h, you can let the infection proceed overnight.

Phage isolation and concentration

- 4 Move the phage lysate to a 50 mL conical tube, and add ~ 1 mL of chloroform. Seal the tube tightly, and shake on a platform shaker for 10–30 min. Transfer to a new conical tube. Then, spin the lysate down in a centrifuge at maximum speed for 30 minutes.
- 5 Carefully pipette off the supernatant from the spun-down culture into a new 50 mL conical tube, avoiding any debris at the chloroform interface. Also avoid the chloroform. Move to a new 50 mL conical tube, and spin again at max speed for another 30 min or until supernatant is clear. Move to a new 50 mL conical tube.
- 6 To concentrate the phage lysate, we have used both PEG precipitation and a filter-concentrator based protocol. The PEG protocol requires an overnight incubation, but the filter-concentrator generally requires more hands-on time. Use whichever phage concentration protocol works best for your circumstances.
Step 6 includes a Step case.

PEG prep
Filter prep

step case

PEG prep

PEG prep requires an overnight incubation step at 4 °C.

- 7 To PEG-precipitate phage, first prepare 5× PEG precipitation solution containing 2.5 M NaCl and 20% w/v PEG8000.

Note

- 8 Add PEG precipitation solution to the phage lysate to obtain a 1× concentration of 0.5 M NaCl and 4% w/v PEG8000. Mix by inverting the tube several times and refrigerate overnight at 4 °C.
- 9 After the overnight incubation with PEG precipitation solution at 4 °C, pellet the PEG precipitated phage particles by centrifuging at 19,000 × g for 60 min at 4 °C. A small white pellet should form at the bottom of the 50 mL conical tube.
- 10 Remove the supernatant, being careful not to disturb the PEG-precipitated phage pellet.
- 11 Resuspend the PEG-precipitated phage pellet in 300 µL SM buffer. Store at 4 °C.