



Dec 13, 2021

## SM buffer

**Frej Larsen**<sup>1</sup><sup>1</sup>Copenhagen University

1

[dx.doi.org/10.17504/protocols.io.b2tdqei6](https://dx.doi.org/10.17504/protocols.io.b2tdqei6)**FOOD Micro UCPH****Frej Larsen**

SM buffer is used for diluting and storing bacteriophages.

SM buffer contains:

200 mM NaCl<sub>2</sub>

10 mM MgSO<sub>4</sub>

50 mM Tris-HCl, pH 7.5

For storage at -20 or -80° C, 15% glycerol can be added before autoclaving.

The buffer can keep at room temperature for several months.

DOI

[dx.doi.org/10.17504/protocols.io.b2tdqei6](https://dx.doi.org/10.17504/protocols.io.b2tdqei6)

Frej Larsen 2021. SM buffer. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b2tdqei6>



\_\_\_\_\_ protocol ,

Dec 13, 2021

Dec 13, 2021

55877

This recipe is for 1 liter of buffer. Make sure to recalculate the measurements if you are making a smaller or larger volume.

Ensure that your  $\text{MgSO}_4$  is a monohydrate before weighing. If not, recalculate amount required to reach correct concentration.






For 1 liter of SM buffer:

11.7 g  $\text{NaCl}_2$

1.4 g  $\text{MgSO}_4$  (monohydrate)

50 mL Tris-HCl, pH 7.5 (1 M)

950 mL demineralized water

- 1 Weigh off  **11.7 g**  $\text{NaCl}_2$  and  **1.4 g**  $\text{MgSO}_4$  and add to a clean 1L blue cap flask or other autoclave safe container.
- 2 Measure  **50 mL** Tris-HCl and add to the flask.
- 3 Measure  **950 mL** demineralized water and add to the flask.
- 4 Screw the cap on tight and mix thoroughly for  **00:00:30** . 30s
- 5 Loosen the cap and autoclave.
- 6 If there are still undissolved particles, mix until dissolved. Then let it cool at room temperature.