

Sep 11, 2024

## CHCHD2 T61I Mouse Genotyping

DOI

[dx.doi.org/10.17504/protocols.io.kxygxyejwl8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygxyejwl8j/v1)

Szu-Chi Liao<sup>1</sup>, Ken Nakamura<sup>1</sup>

<sup>1</sup>Gladstone Institutes

ASAP Collaborative Rese...



Haru Yamamoto

UCSF

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.kxygxyejwl8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygxyejwl8j/v1)

**Protocol Citation:** Szu-Chi Liao, Ken Nakamura 2024. CHCHD2 T61I Mouse Genotyping. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kxygxyejwl8j/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 19, 2024

**Last Modified:** September 11, 2024

**Protocol Integer ID:** 98456

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

**ASAP**

**Grant ID:** 020529

## Abstract

To genotype CHCHD2 T61I point mutant mice by qPCR.



## Mouse Tail Lysis

### 1 **Materials**

- Isoflurane
- Sanitized scissors
- DirectPCR lysis reagent (mouse tail, Viagen Biotech, 102-T)
- Proteinase K – Recombinant, PCR grade (Roche, 03115844001)
- Microcentrifuge tubes
- Incubator
- Water bath

2 Anesthetize mice with isoflurane.

3 Cut no longer than 2 mm of distal tail from mice.

4 Place each tail in a sterile, labeled tube.

5 Add 200 µl of DirectPCR lysis reagent.

6 Add 10 µl of proteinase K and mix well.

7 Incubate the tails at 55°C overnight.

8 To deactivate proteinase K, incubate tails for 45 min in an 85°C water bath.

9 Store extracted DNA at -20°C until use.

## qPCR Genotyping

### 10 **Materials**

- Kapa Probe Fast qPCR Master Mix (2X) kit (Kapa Biosystems, KK4715)
- EndPoint-F Primer: TGAAGATGGCCCAGATCTG



- EndPoint-R Primer: CTGAAGCCCCCAGTGAT
- WT-Probe: (5' JOE) TTAGATGGCTACCACCGCG (3' Iowa Black)
- MT-Probe: (5' 6-FAM) TTAGATGGCGATCACCAGCG (3' Iowa Black)
- Nuclease-free water

11 Make Reaction Buffer:

A	B	C
Reagent	Volume (μl)	Final Concentration
qPCR Master Mix (2X)	10	1X
High ROX (50X)	0.4	1X
EndPoint-F (10 μM)	0.8	0.4 μM
EndPoint-R (10 μM)	0.8	0.4 μM
WT-Probe (10 μM)	0.4	0.2 μM
MT-Probe (10 μM)	0.4	0.2 μM
DNA		100 ng
Nuclease-free water	Make up to 20 μl	
Total	20 μl	

12 Open SDS 2.4 software, create a new experiment and choose "standard curve".

13 Load the plate.

14 Set thermal cycling as follows:

Step	Temperature (°C)	Duration
1	95	3 min
2	95	10 s
3	60	30 s

4	Repeat steps 2-3 for 40 cycles	
5	4	$\infty$

15 In the Setup tab, add detectors for FAM and JOE signal and start the program.

16 After amplification is done, create a new experiment and choose “allelic discrimination”.

17 Assign the markers which detect FAM and JOE signals.

18 Start “Post Read”.

19 **Expected Results:**

When plotted WT against MT fluorescence signal, each genotype would show a distinct ratio/slope as below:

