



Aug 17, 2021

## Cell culture, transfection and imaging

Marianna Leonzino<sup>1,2</sup>, Andrés Guillén-Samander<sup>1,2</sup>, Pietro De Camilli<sup>1,2</sup><sup>1</sup>Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

1 Works for me



Share

[dx.doi.org/10.17504/protocols.io.bvgmn3u6](https://dx.doi.org/10.17504/protocols.io.bvgmn3u6)

Andrés Guillén-Samander

### ABSTRACT

This protocol details the general preparation of cells for imaging and also for imaging experiments involving cellular hypotonic shock and cytosolic Ca<sup>2+</sup> changes as they were performed in <https://doi.org/10.1083/jcb.202010004>.

### ATTACHMENTS

[dn3vbgtxx.pdf](#)

### DOI

[dx.doi.org/10.17504/protocols.io.bvgmn3u6](https://dx.doi.org/10.17504/protocols.io.bvgmn3u6)

### EXTERNAL LINK

<https://doi.org/10.1083/jcb.202010004>

### PROTOCOL CITATION

Marianna Leonzino, Andrés Guillén-Samander, Pietro De Camilli 2021. Cell culture, transfection and imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bvgmn3u6>

### KEYWORDS

Cell culture, Transfection, Imaging, ASAPCRN

### LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Jun 02, 2021

### LAST MODIFIED

Jul 02, 2021

### OWNERSHIP HISTORY

Jun 02, 2021 Urmilas

Jun 08, 2021 Andrés Guillén-Samander

### PROTOCOL INTEGER ID

50413

## MATERIALS TEXT

 CMV-R-

- [GECO1.2 addgene Catalog #45494](#)

## DMEM:

A	B
FBS	10%
Sodium pyruvate	1 mM
Penicillin	100 U/ml
Streptomycin	100 mg/mL
L-glutamine	2 mM

- Live-cell imaging solution (Life Technologies)
- Thapsigargin (Life Technologies)
- BAPTA-AM (Thermo Fisher Scientific)

## General preparation

6h

- 1 Culture the COS-7 or HeLa (ATCC) cells at  **37 °C** and 5% CO<sub>2</sub> in DMEM containing 10% FBS,  **1 Milimolar (mM)** sodium pyruvate,  **100 U/ml** penicillin,  **100 mg/ml** streptomycin and  **2 Milimolar (mM)** L-glutamine (all from Gibco).

**Note:** For general maintenance, when cells reached 80-90% confluency, they were deattached from the dish with Trypsin and diluted 1:20 in a new dish.



- 2 For imaging experiments, seed the cells on glass-bottomed dishes (MatTek) at a concentration of  $75 \times 10^3$  cells per dish<sup>6h</sup> and transfect transiently after  **06:00:00** using FuGene HD (Promega).



Image the cells for 36-48 hours after transfection.

- 4 Just before imaging, remove the growth medium and replace with pre-warmed live-cell imaging solution (Life Technologies).



Perform all live-cell imaging at  **37 °C** and 5% CO<sub>2</sub>.



Perform spinning-disk confocal microscopy using an Andor Dragonfly system equipped with a plan apochromat objective (63x, 1.4 NA, oil) and a Zyla scientific CMOS camera.

- 7 Identify the cells to be imaged by scanning the dish.

### Hypotonic shock experiments

- 8 Once a field of view with transfected cells is found, start the acquisition (generally use a rate of 12 frames/minutes).
- 9 Replace the live-cell imaging solution with pre-warmed distilled water.

**Suggestion:** Minimal volume of live cell imaging media is placed in the dish to allow removal with a pipette in one single step. On the contrary, maximal volume of distilled water is added to allow the highest dilution of leftover ions in the dish.

**Warning:** If media exchange is to be performed by hand, extreme care should be taken not to move the dish during the process. If substantial amount of live cell imaging media needs to be left in the dish to prevent movements during the aspiration procedure, a further media replacement step could be implemented to achieve maximal hypotonic shock.

### Experiments evaluating cytosolic Ca<sup>2+</sup> changes

20m

10

Cytosolic Ca<sup>2+</sup> was monitored by the intensity of the RFP genetically encoded Ca<sup>2+</sup> indicator for optical imaging (R-GECO; plasmid was a gift from R. Campbell, University of Alberta, Edmonton, AB, Canada; Addgene catalog no. 45494)

Once a field of view with transfected cells is found, start the acquisition (generally use a rate of 0.5 Hz).

11



To acutely increase cytosolic Ca<sup>2+</sup>, add Thapsigargin (Life Technologies) to a final concentration of **[M]2 Micromolar (μM)**.

**Note:** Live-cell imaging solution (Life Technologies) contains **[M]1.8 Milimolar (mM)** Ca<sup>2+</sup>.

12

Allow cells to recover for **⌚ 00:10:00**.

10m

13



To decrease cytosolic Ca<sup>2+</sup>, add EGTA and BAPTA-AM (Thermo Fisher Scientific) to the medium to a final concentration of **[M]4 Milimolar (mM)** and **[M]10 Micromolar (μM)**, respectively.

14



10m

Image the cells for  **00:10:00** , then stop the acquisition.