



MAR 26, 2024

🌐 Triparental mating with pSEVA protocol

Laura Gómez¹

¹CBM Severo Ochoa



Laura Gómez
CBM Severo Ochoa

ABSTRACT

This process involves bacterial conjugation, where a conjugative plasmid found in one bacterial strain facilitates the transfer of a mobilizable plasmid from a second bacterial strain to a third bacterial strain.

In the method from our lab, CC118λpir containing a pSEVA plasmid works as the donor, the *E. coli* 1047 pRK2013 strand is the helper, and the receiver is a *C. rodentium* strand.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.n92ldm4zol5b/v1

Protocol Citation: Laura Gómez 2024. Triparental mating with pSEVA protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldm4zol5b/v1>


License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 26, 2024

Insertion of the suicide plasmid by three partner conjugation

3d 2h 55m

1 Inoculate  Overnight cultures of strains:

1d

DONOR CC118λpir pSEVA Gm (10 µg/ml)

HELPER *E. coli* 1047 pRK2013 Km (50 µg/ml)



RECEIVER *C. rodentium* pACBSR Sm (50 µg/ml)

Note

We have a pACBSR plasmid encoding Cm resistance.

2 Place 20 µl spots of the helper and donor strains onto an LB plate (no antibiotics) and an additional spot of **20 µl of the helper on top of 20 µl of the donor (D+H)**.

10m

3 Leave the plate open at the flame until the spots get dry. Incubate at  37 °C  02:00:00 , facing up.



2h

4 Add **40 µl of the receiver strain on top of the D+H spot (D+H+R)** and an additional 20 µl spot of this strain alone. Wait for the spot to get dry and incubate for  Overnight  37 °C facing up.

1d 0h 20m


5 Collect the 4 patches using a sterile loop and resuspend each of them in 1 ml of LB in an eppendorf.

10m


- 6 Plate 100 µl of each tube in LB plates supplemented with **Gm + Sm**. 10m
- 7 Centrifuge the rest of the D+H+R tube at  2000 x g, 00:02:00 to pellet the cells, resuspend in 100 µl and plate as well. 5m
- 8 Incubate the plates  Overnight at 37°C. 1d



Second recombination for a scarless genomic modification

2d 10h 10m

- 9 Pick two colonies of the D+H+R plate and grow them in **LB+Sm + L-arabinose broth at 0.4%** for a minimum of  06:00:00 for the induction of the I-SceI endonuclease of the pACBSR plasmid. 6h

Note

In the  Overnight plates with only donor, helper or receiver no colonies should grow.

- 10 Insert the inoculation loop in the culture and streak on LB+Sm plates to obtain individual colonies. Incubate the plates  Overnight at 37°C. 1d 0h 10m
- 11 The next day, pick some colonies and patch them on a **LB+Sm** plate and on a **LB+Sm+Gm** plate. Incubate the plates  Overnight at 37°C. 1d 0h 30m

Note

Colonies that do not grow on the plate with Gm are those that have recombined after treatment with the endonuclease I-SceI.

- 12 Analyze by PCR and gel electrophoresis those colonies which have grown on LB+Sm but not on LB+Sm to differentiate the modified colonies from the ones which have reverted to the Wild-Type genotype. **NOTE:** The primers should hybridise outside of the homology regions selected. Upon analysis of 10 colonies you should get a about 50% of modified colonies. 3h 30m
- 13 **NOTE:** To remove the pACBSR plasmid, make 8-9 passes of the strain without Sm in liquid LB. Plate the last culture in LB plates and patch individual colonies the following day on LB and LB+Sm plates to select those sensitive to the antibiotic.