

Jul 24, 2024

## Ty + AS + Thymidine Protocol

This protocol is a draft, published without a DOI.

Eleanor Harman<sup>1</sup>

<sup>1</sup>UGA



Eleanor Harman

UGA

---

OPEN  ACCESS



**Protocol Citation:** Eleanor Harman 2024. Ty + AS + Thymidine Protocol. protocols.io <https://protocols.io/view/ty-as-thymidine-protocol-dhri354e>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 23, 2024

**Last Modified:** July 24, 2024

**Protocol Integer ID:** 103946

## Abstract

.



## Media

---

- 1 Yeast Extract 1.5g/500mL  
Tryptone 2.5g/500mL  
Thymidine 150mg/500mL
- 2 Add 80% final volume of DI water to beaker
- 3 Add Yeast Extract, Tryptone, and Thymidine until fully dissolved
- 4 Adjust pH to 5.5
- 5 Autoclave on Liquid 20 minute cycle