



MAR 29, 2024

# Proteomic Analysis of Human Ovarian Cortex and Medulla Secretome Using Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS) Acquisition by Data-Independent Acquisition (DIA) on an Orbitrap Eclipse Tribrid Mass Spectrometer

Joanna Bons<sup>1</sup>, J P Rose<sup>1</sup>, M A Watson<sup>1</sup>, B Schilling<sup>1</sup>

<sup>1</sup>Buck Institute for Research on Aging

Cellular Senescence Network (SenNet) Method Development Community

Schilling Lab

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.n92ld83bnv5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ld83bnv5b/v1)

**Protocol Citation:** Joanna Bons, J P Rose, M A Watson, B Schilling 2024. Proteomic Analysis of Human Ovarian Cortex and Medulla Secretome Using Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS) Acquisition by Data-Independent Acquisition (DIA) on an Orbitrap Eclipse Tribrid Mass Spectrometer. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ld83bnv5b/v1>

**MANUSCRIPT CITATION:**

[dx.doi.org/10.17504/protocols.io.e6nvw1px7lmk/v](https://dx.doi.org/10.17504/protocols.io.e6nvw1px7lmk/v)



M A Watson  
Buck Institute

## ABSTRACT

Ovaries from human donors were cut into 3-5 mm sections. These sections were further processed into 500 µm slices containing cortex and medulla. The slices were then processed into pieces (1 mm x 1 mm x 500 µm), and cortex and medulla pieces were cultured separately as explants in static cultures.

Explants were cultured and treated with either DMSO vehicle control or 0.1 µg/mL doxorubicin to induce

senescence. After 10 days, cortex and medulla explants were thoroughly washed with serum-free basal media and transferred to a clean plate with pre-equilibrated serum-free basal media and inserts. The conditioned media were collected after 24 hours for secretome proteomics profiling.

The concentrated conditioned media was subjected to tryptic digestion using S-trap Spin columns. The reconstituted peptide elution was desalted with C18 hydrophilic-lipophilic balance (HLB) cartridges. The final reconstituted peptides were diluted with 2% ACN and 0.1% FA. Proteolytic peptide measurement was completed using liquid chromatography-tandem mass spectrometry (LC-MS/MS) acquisition by Data-Independent Acquisition (DIA) on an Orbitrap Eclipse Tribrid mass spectrometer for peptide/protein identification and quantification.

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Mar 25, 2024

**Last Modified:** Mar 29, 2024

**PROTOCOL integer ID:** 97359

**Keywords:** Human, Ovary, Secretome, SASP, Senescence, Doxorubicin, LC-MS/MS, DIA, Orbitrap Eclipse, Mass Spectrometry, Conditioned Media

**Funders Acknowledgement:**  
NIH Grant  
Grant ID: 5U54AG075932-02  
NIH Grant  
Grant ID: 5U54AG075932-03

## Conditioned Media Concentration with Amicon Ultra Centrifugal Filters

1 [dx.doi.org/10.17504/protocols.io.e6nvw1px7lmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvw1px7lmk/v1)

## Protein Digestion with S-trap Spin Columns using Conditioned Concentrat...

2 [dx.doi.org/10.17504/protocols.io.x54v928eml3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v928eml3e/v1)

## Proteolytic Peptide Desalting with C18 HLB Cartridges

3 [dx.doi.org/10.17504/protocols.io.eq2lywdzpvx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lywdzpvx9/v1)

## LC-MS/MS Acquisition by DIA on an Orbitrap Eclipse Tribrid Mass Spectro...

4 [dx.doi.org/10.17504/protocols.io.36wgq3m15lk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgq3m15lk5/v1)

## DIA Data Processing using Spectronaut/directDIA (Biognosys): Secretome...

5 [dx.doi.org/10.17504/protocols.io.q26g71rz3gwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g71rz3gwz/v1)