



Cell cryopreservation V.2

PMAT0001¹

¹Nanyang Technological University

Version 2

Oct 21, 2020

1

Works for me

This protocol is published without a DOI.

PMAT0001

PROTOCOL CITATION

PMAT0001 2020. Cell cryopreservation. **protocols.io**
<https://protocols.io/view/cell-cryopreservation-bnp7mdrn>
Version created by **PMAT0001**

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 21, 2020

LAST MODIFIED

Oct 21, 2020

PROTOCOL INTEGER ID

43487

GUIDELINES

- This protocol is for ADHERENT CELLS

MATERIALS

NAME	CATALOG #	VENDOR
Cell counter		Thermo Fisher Scientific
DMSO	9224	J.T. Baker
Bright-Line Hemacytometer	02-671-5	Fisher Scientific
Trypan Blue Stain 0.4%		Invitrogen - Thermo Fisher

MATERIALS TEXT

- Use a DMSO bottle only for cell culture work; Open only in laminar flow hood

BEFORE STARTING

- Aseptic techniques
- Wipe down hood and any item introduced into the hood with 70% ethanol

- 1 Gently detach cells from the tissue culture flask in the following manner:


- 1.1 Aspirate media



- 1.2 Wash with  2 mL PBS, then aspirate.

1.3 Wash with  2 mL trypsin.

1.4 Incubate for  00:05:00 at  37 °C .

5m


2 Resuspend the cells in  10 mL of DMEM with FBS. Aliquot into a centrifuge tube. Of course, if there is a T25 flask, resuspend in lesser volume.

3 Aliquot a small amount of the mixture (about  200 µl) for cell counting. Ideally, cell viability should be in excess of 90% in order to achieve a good recovery after freezing. Keep at  -4 °C
- We need 1×10^6 - 5×10^6 cells for freezing.

4 Centrifuge the cell suspensions at 1500rpm for  00:05:00

5m

5 Discard supernatant carefully

6 Add about  3 mL of freezing media into the tube and re-suspend.

6.1 Freezing media preparation:
- 10% DMSO, 20% FBS, 70% appropriate cell culture medium (with serum + antibiotics)
(e.g 100uL of DMSO + 200uL FBS + 700uL cell culture medium with FBS and antibiotics)

7 Dispense the aliquots of cell suspension into the cryogenic storage vials.

8 Cool the vials gradually as follows:

1h

 -20 °C for  01:00:00

 -80 °C overnight

Liquid Nitrogen