



Aug 24, 2020

# Determination of IgA concentration by the Mancini test.

Angel A Justiz-Vaillant<sup>1</sup><sup>1</sup>University of the West Indies St. Augustine

1

Works for me

[dx.doi.org/10.17504/protocols.io.bj93kr8n](https://dx.doi.org/10.17504/protocols.io.bj93kr8n)

University of the West Indies

[angel.vaillant@sta.uwi.edu](mailto:angel.vaillant@sta.uwi.edu)

Angel Justiz-Vaillant

University of the West Indies St. Augustine

DOI

[dx.doi.org/10.17504/protocols.io.bj93kr8n](https://dx.doi.org/10.17504/protocols.io.bj93kr8n)

PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Determination of IgA concentration by the Mancini test.. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bj93kr8n>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 24, 2020

LAST MODIFIED

Aug 24, 2020

PROTOCOL INTEGER ID

40987

- 1 An appropriate anti-IgA antiserum (antibody) is poured in the center well of an agar-containing plate.
- 2 Carefully circular wells are cut and detached from the plates.
- 3 A series of standards containing known concentrations of IgA are placed in separate wells, while “unknown” human serum samples and control are placed in other wells.
- 4 A ring of precipitate forms in the area of optimal antigen-antibody concentration, meaning anti-IgA - IgA reaction as the antigen diffuses radially.
- 5 The diameters of the rings are measured and perceived normally in 48-72 hours.

6 Finally, a standard curve is developed using the ring diameters of the standards versus the concentrations.

7 A curve is then used to plot the concentration of the control and unknown IgA samples.