



Feb 13, 2022

Gentamicin Protection assay (Intracellular Survival Assay) for *Salmonella* Typhimurium/Typhi

Kasturi Chandra¹¹Department of Microbiology and Cell Biology, Indian Institute of Science, India

1

dx.doi.org/10.17504/protocols.io.b4zmqx46

IISc_DC Lab_MCBL

Kasturi Chandra

Salmonella invades the epithelial cells by inducing bacterial uptake, whereas phagocytes readily take up the bacteria by phagocytosis. Intracellular *Salmonella* is resistant to gentamicin, but extracellular bacteria are susceptible to the antibiotic, therefore, adhesion, invasion, and bacterial replication can be measured using this protocol. Inside the cells, *Salmonella* proliferates and the rate of proliferation can be determined by enumerating the bacterial CFU from the infected cells.

DOI

dx.doi.org/10.17504/protocols.io.b4zmqx46

Kasturi Chandra 2022. Gentamicin Protection assay (Intracellular Survival Assay) for *Salmonella* Typhimurium/Typhi. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4zmqx46>



protocol ,

Feb 13, 2022

Feb 13, 2022

58125

1 Epithelial cells were infected with mid-log phase culture of bacteria grown in LB (OD600 0.3)

OR

phagocytes were infected with overnight culture (OD₆₀₀ 0.3). The multiplicity of infection (MOI) of 10 was used in each case.

- 2 Bacterial attachment to host cells was enhanced by centrifuging at 600 rpm for 10 min.
- 3 After 25 min of infection, cells were washed three times with 1x PBS and treated with gentamicin (100 µg/ml in complete media) for 1 hour and then maintained with 25 µg/ml gentamicin for the rest of the experiment.
- 4 0.1% Triton-X 100 (v/v in 1x PBS) was used to lyse the cells and the lysate was plated on Salmonella-Shigella (SS) agar for *S. Typhimurium* strains and Wilson Blair (WB) agar for *S. Typhi* strains.
- 5 For adhesion assay:
The cells were washed three times with 1x PBS and lysed immediately after the initial 25 min incubation. The percent adhesion was calculated with respect to the pre-inoculum used for infection.
- 6 For invasion assay:
The cells were washed three times with 1x PBS and lysed after incubation in 100 µg/ml gentamicin treatment (i.e. 1 hour post infection) and percent invasion was calculated with respect to the pre-inoculum used for infection.
- 7 For intracellular survival assay (ICSA):
The cells were washed three times with 1x PBS and lysed at 2 hours and 18 hours post infection. CFU at 18 hours was divided by CFU at 2 hours to obtain fold replication of the intracellular bacteria.