



Mar 04, 2022

RNA Extraction and RT-qPCR

Haley Geertsma¹¹University of Ottawa

1

dx.doi.org/10.17504/protocols.io.b5wmq7c6 Haley Geertsma
University of Ottawa

This protocol is used to extract RNA from mouse brain tissue and quantify with RT-qPCR.

DOI

dx.doi.org/10.17504/protocols.io.b5wmq7c6

Haley Geertsma 2022. RNA Extraction and RT-qPCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5wmq7c6>



protocol ,

Mar 04, 2022

Mar 04, 2022

59053

- 1 Homogenize hemibrain in 3mL 1X PEPI buffer using a dounce homogenizer. 2m
1X PEPI Buffer: 5mM EDTA + 1X protease inhibitor in 1X phosphate buffered saline (PBS)
- 2 Take 3% of homogenate and add 1mL TRIzol Reagent and follow their user guide. 45m
- 3 Synthesize cDNA with 5X All-In-One RT Master Mix Kit then perform qPCR. 2h
qPCR settings: 95°C for 5 minutes, 40x (95°C for 15 seconds then 60°C for 60 seconds),
melting curve