



APR 07, 2023

🌐 QIAGEN DNeasy PowerMax Soil Kit

QIAGEN¹

¹QIAGEN



Yuanyu Cheng
McGill University

ABSTRACT

For the isolation of microbial DNA from large quantities of soil - great for samples with low microbial load

The DNeasy PowerMax Soil Kit comprises a novel and proprietary method for isolating genomic DNA from environmental samples using Inhibitor Removal Technology® (IRT). With this kit, it is possible to process samples that have proven difficult in the past due to high levels of humic-like substances. The isolated DNA has a high level of purity, which allows for successful PCR amplification from samples. Total DNA isolated from various soil types has been successfully amplified using PCR with primers specific for bacteria (*Bacillus subtilis*, *Bacillus anthracis*), fungi (yeast, mold), and actinomycetes (*Streptomyces*).

Using the DNeasy PowerMax Soil Kit, environmental samples are added to a bead-beating tube with a kit-supplied proprietary buffer for rapid and thorough homogenization. Cell lysis and DNA exposure occur by mechanical and chemical methods. Extracted genomic DNA is captured on a silica membrane in a spin column format. The DNA is washed and eluted from the membrane and is ready for PCR and other downstream applications.

Protocol successfully used by Sales et al., 2021 to characterize fish species diversity and richness in water and sediment samples

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.6qpvr4ke3gmk/v1

Protocol Citation: QIAGEN 2023. QIAGEN DNeasy PowerMax Soil Kit. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr4ke3gmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working Protocol used successfully by Sales et al., 2021 to detect fish sedDNA

Created: Dec 21, 2022

Last Modified: Apr 07, 2023

PROTOCOL integer ID:
74334

Keywords: QIAGEN, DNeasy, PowerMax, SedDNA, Sedimentary DNA

MATERIALS

DNeasy PowerMax Soil Kit (10 preparations) already includes:

MB Maxi Spin Columns * 10

PowerMax Bead Pro Tubes * 10

PowerBead Solution  20 mL

Solution C1  6.6 mL * 2

Solution C2  28 mL * 2

Solution C3.  44 mL

Solution C4  330 mL

Solution C5  30 mL * 4

Solution C6  66 mL

Collection Tubes (50ml) 5 * 8

Quick Start Protocol * 1

Additional equipments and reagents to be supplied by user:

Centrifuge capable of spinning 50 ml tubes at 2500 x g using swing-out rotor

Pipettes (1ml and 10 ml)

Vortex-Genie 2 Vortex

Vortex Adapter for 2 (50 ml) tubes (cat. no. 13000-V1-50)

SAFETY WARNINGS




Solution C5 contains ethanol and is flammable.


DO NOT add bleach or acidic solutions directly to the sample preparation waste.

Sample preparation & cell lysis

1h 27m 30s

1 **ADD**  15 mL of PowerBead Solution to a PowerMax Bead Pro Tube









1m

ADD up to  10 g of soil sample to the PowerMax Bead Pro Tube containing PowerBead Solution

VORTEX vigorously for  00:01:00















Note

Refer to manufacturer's Troubleshooting Guide before deciding on the amount of soil to process. However, higher volumes of sediment (10 g) have shown better detection rates for fish sedDNA.

- 2 **ADD**  1.2 mL of Solution C1 to the PowerMax Bead Pro Tube 30s
- VORTEX** vigorously for  00:00:30
- 3 **PLACE** PowerMax Bead Pro Tube on a vortex adapter 40m
- VORTEX** for  00:10:00 at the highest speed
- ALTERNATIVELY**, place the tube in a shaking water bath set at  65 °C and shake at maximum speed for  00:30:00
- 4 **CENTRIFUGE** at  2500 x g for  00:03:00 at  Room temperature 3m
- TRANSFER** supernatant to a clean collection tube

Inhibitor removal


1h 27m 30s

- 5 **ADD**  5 mL of Solution C2 10m
- INVERT** twice to mix
- INCUBATE** at  2 °C to  8 °C for  00:10:00
- 6 **CENTRIFUGE** at  2500 x g for  00:04:00 at  Room temperature 4m
- AVOIDING** the pellet, transfer the supernatant to a clean collection tube
- 7 **ADD**  4 mL of Solution C3 10m
- INVERT** twice to mix
- INCUBATE** at  2 °C to  8 °C for  00:10:00
- 8 **CENTRIFUGE** at  2500 x g for  00:04:00 at  Room temperature 4m
- AVOIDING** the pellet, transfer the supernatant to a clean collection tube

Bind DNA

1h 27m 30s

9 **SHAKE** to mix Solution C4

ADD  30 mL of Solution C4 to supernatant

INVERT twice to mix

10 **FILL** an MB Maxi Spin Column with the solution from step 9

2m

CENTRIFUGE at  2500 x g for  00:02:00 at  Room temperature

DISCARD the liquid flow-through

11 **ADD** a second volume of supernatant to the same MB Maxi Spin Column

2m

CENTRIFUGE at  2500 x g for  00:02:00 at  Room temperature

DISCARD the liquid flow-through

REPEAT until the entire volume has been processed. This will take up to **4 total spins**

Wash spin column

1h 27m 30s

12 **ADD**  10 mL of Solution C5

3m

CENTRIFUGE at  2500 x g for  00:03:00 at  Room temperature

DISCARD the liquid flow-through


13 **CENTRIFUGE** again at  2500 x g for  00:05:00 at  Room temperature to remove residual Solution C5

5m

CAREFULLY place the MB Maxi Spin Column in a new collection tube. Avoid splashing any residual Solution C5 onto the column

Elute the DNA

1h 27m 30s

14 **ADD**  5 mL of sterile Solution C6 to the center of MB Maxi Spin Column membrane

3m

CENTRIFUGE at  2500 x g for  00:03:00 at  Room temperature

DISCARD the MB Maxi Spin Column

DNA is now ready for downstream applications