



Jul 01, 2020

Phenzine Oxidizer Enrichment and Isolation

Lev Tsypin¹, Yinon Bar-On², Scott Saunders¹, Jared R Leadbetter¹, Dianne K Newman¹¹California Institute of Technology; ²Weizmann Institute of Science

1 Works for me dx.doi.org/10.17504/protocols.io.bh4tj8wn

 Lev Tsypin

ABSTRACT

This is a protocol for the enrichment of phenazine-1-carboxylic acid oxidizing microbes from a soil sample.

DOI

dx.doi.org/10.17504/protocols.io.bh4tj8wn

PROTOCOL CITATION

Lev Tsypin, Yinon Bar-On, Scott Saunders, Jared R Leadbetter, Dianne K Newman 2020. Phenzine Oxidizer Enrichment and Isolation. **protocols.io**
dx.doi.org/10.17504/protocols.io.bh4tj8wn

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 01, 2020

LAST MODIFIED

Jul 01, 2020

PROTOCOL INTEGER ID

38771

Medium composition

1 Stock solutions

100x Freshwater salts:

- 1.71 M sodium chloride
- 197 mM magnesium chloride
- 68 mM calcium chloride
- 671 mM potassium chloride

1000x Trace Elements:

- 20 mM HCl
- 7.5 mM FeCl₃ 6H₂O
- 480 uM H₃BO₃
- 500 uM MnCl₂ 4H₂O
- 6.8 mM CoCl₂ 6H₂O
- 1 mM NiCl₂ 6H₂O
- 12 uM CuCl₂ 2H₂O
- 500 uM ZnCl₂
- 23 uM Na₂SeO₃

- 150 µM Na₂MoO₄

1000x 13-vitamin solution

- 1 mM MOPS pH 7.2
- 100 µg/mL riboflavin
- 30 µg/mL biotin
- 100 µg/mL thiamine HCl
- 100 µg/mL L-ascorbic acid
- 100 µg/mL d-Ca-pantothenate
- 100 µg/mL folic acid
- 100 µg/mL nicotinic acid
- 100 µg/mL 4-aminobenzoic acid
- 100 µg/mL pyridoxine HCl
- 100 µg/mL lipoic acid
- 100 µg/mL Nicotinamide adenine dinucleotide
- 100 µg/mL thiamine pyrophosphate
- 10 µg/mL cyanocobalamin

2 mM reduced phenazine-1-carboxylic acid (PCA)

- In a sealed Balch tube or serum vial, sparge 2 mM oxidized PCA with 80:20 (vol:vol) H₂/CO₂ over Pd (II) until yellow-green in color
- Stir with a magnetic stirbar or otherwise agitate during sparging and overnight afterwards.

Enrichment medium (anoxic)

- 45 mM sodium bicarbonate buffer (buffer)
- 10 mM ammonium chloride (nitrogen source)
- 1 mM potassium phosphate (phosphorus source)
- 1 mM reduced phenazine-1-carboxylic acid (PCA, electron donor)
- 100 µM sodium acetate (carbon source)
- 50 µM sodium sulfide (sulfur source)
- 1x trace elements solution
- 1x freshwater salts solution
- 1x 13-vitamin solution

Isolation medium (oxic or anoxic, agar plates)

- 1 mM sodium phosphate buffer pH 7
- 10 mM ammonium chloride
- 100 µM sodium sulfate
- 1 mM sodium acetate
- 10 mM sodium nitrate
- 1x freshwater salts solution
- 1x trace elements solution
- 1x 13-vitamin solution
- 1.5 % agar

Sampling

- 2 Collect topsoil by scooping directly with 50 mL conical tube.
- 3 Grind with mortar and pestle

Pass through 2 mm sieve to remove large rocks and organic matter.

4

5 Transfer to anaerobic chamber

Inoculation of enrichment cultures

6 Prepare 96-well plate with 190 μ L of anoxic enrichment medium per well.

7 Add 20 mg soil into each well.

8 Prepare Balch tubes with 9 mL of anoxic enrichment medium in each.

9 Add 1 g soil into each tube and seal.

10 For both the 96-well plates and Balch tubes, watch for a color change over the next days. The yellow-green reduced PCA will turn clear if oxidized.

11 Passage candidate PCA-oxidizing enrichments via 1:10 dilutions into fresh medium, maintaining anoxic conditions.

Isolation

12 Plate enrichment cultures onto LB/Agar and the solid isolation medium.

13 Grow both oxically and anoxically in an anaerobic chamber.

14 Pick colonies of different morphologies and assay PCA-oxidizing capability in the enrichment medium as a pure culture.

15 Select isolates that oxidize PCA in pure culture for further study.