



2 ▼

Feb 10, 2022

Making your own electrocompetent cells

V.2

[New England Biolabs¹](#)¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.bd2zi8f6**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**
New England Biolabs

This protocol explains methods for making two 250 ml cultures of electrocompetent cells.

An Electroporation Protocol can be found [here](#).

DOI

dx.doi.org/10.17504/protocols.io.bd2zi8f6<https://www.neb.com/protocols/2012/06/21/making-your-own-electrocompetent-cells>

New England Biolabs 2022. Making your own electrocompetent cells.

protocols.io<https://dx.doi.org/10.17504/protocols.io.bd2zi8f6>

New England Biolabs



Electrocompetent cells

 protocol ,

Mar 21, 2020

Feb 10, 2022

34617

MATERIALS

 [Magnesium sulfate heptahydrate](#) **Sigma**

Aldrich Catalog #M2773

 [NaCl](#) **Sigma**

Aldrich Catalog #53014

 [Tryptone](#) **Fisher**

Scientific Catalog #BP1421-500

 [Glucose](#) **Sigma**

Aldrich Catalog #G8270

 [Magnesium chloride hexahydrate](#) **Sigma**

Aldrich Catalog #M2670

 [Potassium chloride](#) **Sigma**

Aldrich Catalog #P9333

 [Glycerol](#) **Thermo**

Fisher Catalog #17904

 [Yeast Extract](#) **Thermo**

Fisher Catalog #211930

Media

SOB:

2% tryptone

0.5% yeast extract

10 mM NaCl

2.5 mM KCl

10 mM MgCl₂

10 mM MgSO₄

SOC (for Electroporation Protocol):

SOB + 20 mM glucose

Appropriate Antibiotics for Your Application

Antibiotics for Plasmid selection

Antibiotic	Working Concentration
Ampicillin	100 µg/ml
Carbenicillin	100 µg/ml
Chloramphenicol	33 µg/ml
Kanamycin	30 µg/ml
Streptomycin	25 µg/ml
Tetracycline	15 µg/ml

Sterile 10% glycerol (can be autoclaved) is needed for the washes. The volume of 10% glycerol needed is 2X the culture volume (for example, a 500 ml culture requires 1L of 10% glycerol).

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Sterile **10 % glycerol** (can be autoclaved) is needed for the washes. The volume of **10 % glycerol** needed is **2X the culture volume** (for example, a 500 ml culture requires 1L of **10 % glycerol**).

- 1 Innoculate 1 colony from a fresh plate of the strain to be made electrocompetent into **10 mL SOB** in a 125 ml flask.



Incubate for 16-18 hours at **37 °C** and **250 rpm** .

- 3 Pre-warm two 1 L flasks containing **250 mL each of SOB** pre-warmed to **37 °C** .

- 4 Add two drops of the overnight culture to each of the pre-warmed 250 ml flasks.

- 5 Shake at **37 °C** and **250 rpm** until the cultures reach an OD₆₀₀ of 0.5-0.7.

IMPORTANT:

Be sure to turn on centrifuge and cool rotor to 4°C well in advance of harvesting cells.

Be sure to place **1 L 10% glycerol** **On ice** well in advance of harvesting cells.



Place cultures **On ice** for **00:15:00** . **From this point on the cultures must be kept**

ice cold.

7 Pour each 250 ml culture into chilled 500 ml (or 1000 ml) centrifuge bottles.

8




Centrifuge at  **5000 rpm, 00:10:00** .

9 Pour off the supernatant and aspirate any residual broth.

10



Add  **250 mL glycerol** to each of the centrifuge bottles and completely suspend the cells by pipetting up and down.

11



Centrifuge at  **5000 rpm, 00:10:00** .

12 Pour off the supernatant.

It is not necessary to aspirate any residual broth.

13 Completely suspend the cells in  **250 mL glycerol** and re-centrifuge.

14




Pour off the supernatant and suspend the cells in the **residual** glycerol by pipetting up and down.

15



At this point you can [electroporate](#) or freeze the cells away.

To freeze, add  **100 µL culture** to microcentrifuge tubes  **On ice** .

16 Once you have used all of the culture, transfer the tubes to dry ice for  **00:10:00** .

17 

Once the cultures are frozen, transfer them to a -80°C freezer.

The cultures should be good for >6 months.