

Aug 27, 2024

## DAB Staining for Tyrosine Hydroxylase (TH) on Free-floating Fixed NHP Brain Tissue

 Forked from [Standard DAB Staining for Free-floating Fixed NHP Brain Tissue](#)

DOI

[dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1)

Andreea Bostan<sup>1</sup>

<sup>1</sup>Department of Neurobiology, University of Pittsburgh

ASAP Collaborative Rese...



Andreea Bostan

University of Pittsburgh

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1)

**Protocol Citation:** Andreea Bostan 2024. DAB Staining for Tyrosine Hydroxylase (TH) on Free-floating Fixed NHP Brain Tissue . protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 27, 2024

**Last Modified:** August 27, 2024


**Protocol Integer ID:** 106519

**Keywords:** ASAPCRN, Immunostaining, DAB, NHP Brain Tissue, TH

**Funders Acknowledgement:****Aligning Science Across  
Parkinson's****Grant ID: ASAP-020519**

## Abstract

This protocol details the procedure for immunohistochemical 3,3'-Diaminobenzidine (DAB) staining for Tyrosine Hydroxylase (TH) of free-floating fixed brain tissue sections using the avidin/biotin ABC complex and the Rabbit Polyclonal Anti-TH NB300-109 (Novus).

This protocol has been tested with free-floating non-human primate (NHP) brain tissue that has been fixed (10% formalin or 4% paraformaldehyde), cryoprotected (sucrose or glycerol gradients), and cryo-sectioned  50  $\mu\text{m}$  .

## Guidelines

This protocol has been optimized for using 6 well tissue culture plates [Falcon, 353046] to react individual sections. You will need **2+ mL** solutions for *each* well plate.

It is recommended to use new culture plates for antibody incubation and discard plates used for DAB.

Plates used for other steps can be washed and re-used.



## Materials

### Tissue:

Brain tissue sections (20 - 50  $\mu\text{m}$ ).

### Materials/Equipment:

- Tissue culture plates or circular staining nets
- Orbital shaker
- Fume hood
- Nitrile Gloves
- Glass slides (charged or subbed)

### Reagents:

- Phosphate-buffered saline (PBS)
- Phosphate-buffered saline (PBS) with 0.2% Triton-X
- Akoya Blocking Reagent Powder (SKU FP1012) to make TNB Blocking Buffer/
- Hydrogen Peroxide:  $\text{H}_2\text{O}_2$  (3% or 30%)
- Distilled water:  $\text{dH}_2\text{O}$
- Primary Antibody: Novus Biologicals Tyrosine Hydroxylase (TH) Polyclonal Rabbit.
- Secondary Antibody (to match the host of the primary antibody): biotinylated Goat-Anti Rabbit.
- Normal Serum Blocking Solution (e.g., Normal Goat Serum, S-1000-20)
- Vectastain Elite ABC Peroxidase Kit (Standard) (PK-6100) (Vector Laboratories)
- Vectastain ABC-HRP Kit, Peroxidase (Rabbit IgG) (PK-4001, Vector Laboratories)
- DAB Substrate Kit

Examples:

Peroxidase (HRP) with Nickel (3,3'-diaminobenzidine) (SK-4100) (Vector Laboratories)

ImmPACT DAB (SK-4105)

## Safety warnings

- ! Use appropriate care when using hydrogen peroxide (reactive, can cause skin/eye damage) and DAB (suspected carcinogen). Collect DAB solution for chemical waste disposal.



## Part I (Day 1)

3h

- 1 Bring tissue to Room temperature in **phosphate-buffered saline (PBS)** on an orbital shaker for **30 minutes**.

00:30:00

30m

- 2 Prepare **Peroxide Solution (3 % H<sub>2</sub>O<sub>2</sub>)** in **dH<sub>2</sub>O**.

*E.g.*, for 10 mL 3% H<sub>2</sub>O<sub>2</sub> use:

- 1000 µL 30% H<sub>2</sub>O<sub>2</sub>
- 9000 µL dH<sub>2</sub>O

5m

- 3 Prepare **Blocking Serum Solution: Normal Goat Serum (NGS)** in **TNB Blocking Buffer**.

*E.g.*, in 10 mL buffer (TNB Blocking Buffer) add:

- 150 µL NGS (or 3 drops of normal serum if using an ABC kit, e.g. Vectastain ABC-HRP Kit Rabbit IgG PK-4001)

5m

- 4 Prepare **Primary Antibody Solution: anti-TH NB300-109** at **1:1000 TNB Blocking Buffer**.

*E.g.*, for 10 mL **Primary Antibody Solution** use:

- 9990 µL TNB Blocking Buffer
- 10 µL anti-TH NB300-109

5m

- 5 **Rinse** in **PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature :

**3 x 3 minutes.** 00:09:00

9m

- 6 Quench endogenous peroxide in **Peroxide Solution (3 % H<sub>2</sub>O<sub>2</sub>)** on a shaker at

Room temperature :

**2 x 10 - 15 minutes.** 00:20:00 - 01:00:00

1h 20m

- 7 **Rinse** in **PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature :

**3 x 3 minutes.** 00:09:00

9m

- 8 Incubate in **Blocking Serum Solution** on a shaker at RT: **1 hour**.

01:00:00

1h

**DO NOT RINSE** after blocking serum.



9 Incubate in **Primary Antibody Solution** on a shaker at 4 °C : **overnight x 3 (60 - 72 hours).**

3d

72:00:00 .

## Part II (Day 2)

4h

10 Bring tissue (in the **Primary Antibody Solution**) to Room temperature on a shaker (**45 minutes**). 00:30:00 - 01:00:00

30m

11 Prepare **ABC Solution** in **PBS with 0.2% Triton X (PBS-Tx)** (at least 30 minutes before use).

5m

00:30:00 .

12 Prepare **Secondary Antibody Solution (1:200)** in **TNB Blocking Buffer**.

5m

In 10 mL TNB Blocking Buffer add:

- 150 µL Normal Goat Serum (NGS) (= 3 drops of normal serum if using an ABC kit, e.g. Vectastain ABC-HRP Kit Rabbit IgG PK-4001).
- 50 µL Secondary biotinylated goat-anti rabbit (= 1 drop of normal goat serum if using an ABC kit, e.g. Vectastain ABC-HRP Kit Rabbit IgG PK-4001).

13 **Rinse** in **PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature : **3 x 3 minutes.** 00:09:00 .

9m

14 Incubate in **Secondary Antibody Solution** on a shaker at Room temperature : **30 minutes.** 00:30:00 .

30m

15 **Rinse** in **PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature : **3 x 3 minutes.** 00:09:00 .

9m

16 Incubate in **ABC Solution** on a shaker at Room temperature : **60 minutes.** 01:00:00 .

1h

17 **Rinse** in **PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature :

9m

**3 x 3 minutes.** 00:09:00 .

18 Prepare **Peroxide Substrate Solution** in dH<sub>2</sub>O.

5m

To use the Vector Labs DAB Peroxidase Substrate Kit (SK-4100):

In 5 mL dH<sub>2</sub>O:

- 2 drops Reagent 1
- 4 drops Reagent 2
- 2 drops Reagent 3
- [optional] 2 drops of Reagent 4 (Nickel) if a black reaction product is desired

**Note: Mix well before use. Use immediately.**

19 Incubate in **Peroxide Substrate Solution** on a shaker at Room temperature :

6m

**3 - 10 minutes.** 00:03:00 - 00:06:00 .

**Note:** Watch the tissue closely to avoid high background staining.

20 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature :

9m

**3 x 3 minutes.** 00:09:00 .

21 Mount tissue on glass slides (subbed or charged) in 1:8 PBS in dH<sub>2</sub>O and let air dry.

22 Rinse slides with dH<sub>2</sub>O and let air dry (preferably in a hood).

23 Coverslip clean and dry slides with Cytoseal 60 (Thermo Fisher #830-16).

## Protocol references

[https://vectorlabs.com/productattachments/protocol/VL\\_SK-4100\\_UserGuide\\_LBL02267.pdf](https://vectorlabs.com/productattachments/protocol/VL_SK-4100_UserGuide_LBL02267.pdf)

<https://vectorlabs.com/products/vectastain-elite-abc-hrp-kit-standard>

[https://my.akoyabio.com/ccrz\\_ProductDetails?sku=FP1012&cclcl=en\\_US](https://my.akoyabio.com/ccrz_ProductDetails?sku=FP1012&cclcl=en_US)