

VERSION 1 AUG 14, 2023

# OPEN BACCESS



DOI:

dx.doi.org/10.17504/protocol s.io.36wgq339klk5/v1

**Protocol Citation:** chuqings un 2023. Efficient recovery of complete gut viral genomes by combined short- and longread sequencing.

#### protocols.io

https://dx.doi.org/10.17504/protocols.io.36wgq339klk5/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Other We attempted this protocol but could not get it to work in our workspace

# € Efficient recovery of complete gut viral genomes by combined short- and long-read sequencing V.1

### chuqingsun<sup>1</sup>

<sup>1</sup>Huazhong University of Science and Techonology



chuqingsun

#### **ABSTRACT**

Current metagenome-assembled human gut phage catalogs contained mostly fragmented genomes. Here, we developed a vigorous gut virome detection procedure involving viral-like particle enrichment from increased amount of feces (~500g) and combined sequencing of short- and long-reads. Applied to 135 fecal samples, we assembled a Chinese Gut Virome Catalog (CHGV) consisting of 21,646 non-redundant phage genomes that were significantly longer than those obtained by short-read sequencing and contained ~35% complete ones, which was ~nine times more than those in the Gut Virome Database (~4%). Interestingly, majority (~60%) of the CHGV genomes were obtained by either long-read or hybrid assemblies, which overlapped little with those assembled from only the short-reads, indicating the necessity of combined sequencing in gut virome discovery. With this dataset, we elucidated the vast diversity of the gut virome from several aspects, including the identification of 32% novel genomes as compared with public gut virome databases, dozens of phages that were more prevalent than the crAssphages and/or Gubaphages, and several viral clades that are more diverse than the two. In the end, we also characterized the functional capacities of the CHGV encoded proteins and constructed a viral-host interaction network to facilitate future research and applications of the gut viruses.

Created: Aug 14, 2023

Last Modified: Aug 14,

2023

#### **PROTOCOL** integer ID:

86442

### **Assembly**

### 1 NGS Assembly

```
#!/bin/bash
#SBATCH --cpus-per-task=16
#SBATCH -o slurm.%N.%j.out
                                  # STDOUT
#SBATCH -e slurm.%N.%j.err
                                  # STDERR
infile=$1
NGS PATH=$2
export
LD LIBRARY PATH=$LD LIBRARY PATH:~/local/lib:/mnt/raid6/sunchuqing/S
oftwares/MCR/v94/runtime/glnxa64:/mnt/raid6/sunchuqing/Softwares/MCR
/v94/sys/os/glnxa64:/mnt/raid6/sunchuging/Softwares/MCR/v94/extern/g
lnxa64:/mnt/raid3/wchen/miniconda2/pkgs/libgcc-7.2.0-h69d50b8 2/lib
cd NGS
mkdir -p 02 trimmed 03 bac cpn60 03 human hg38 04 Stats
mkdir -p 05 Removed 06 Assembly 07 CD-HIT
TRIMMO JAR FILE='/mnt/raid1/tools/ngs tools/Trimmomatic-
0.38/trimmomatic-0.38.jar'
TRIMMO ADAPTOR FILE PE='/mnt/raid1/tools/ngs tools/Trimmomatic-
0.38/adapters/TruSeq3-PE.fa'
R1=`ls ${NGS PATH}/*${infile}*R1*|head -n 1 `
R2=ls ${NGS PATH}/*{infile}*R2*|head -n 1 `
if [ ! -s 02 trimmed/${infile} clean.1.fq ];then
    java -jar $TRIMMO JAR FILE PE -threads 16 ${R1} ${R2}
02 trimmed/${infile} clean.1.fq
02 trimmed/${infile} clean unpaired.1.fq
02 trimmed/${infile} clean.2.fq
02 trimmed/${infile} clean unpaired.2.fq
ILLUMINACLIP: $TRIMMO ADAPTOR FILE PE:2:15:10 LEADING:3 TRAILING:3
SLIDINGWINDOW:15:30 MINLEN:50
fi
#rm human
export
PATH=/home/sunchuqing/bin:/mnt/raid6/sunchuqing/Softwares/miniconda3
/condabin:/mnt/raid6/sunchuqing/Softwares/miniconda3/bin:/mnt/raid8/
```

```
sunchuging/Softwares/bin:/mnt/raid1/puzi/software/metaphlan2:/mnt/ra
id1/puzi/software/metaphlan2/utils:/usr/bin:/usr/local/ncbi/sra-
tools/bin:/usr/local/sbin:/usr/local/bin:/usr/sbin:/usr/bin:/b
in:/usr/games:/usr/local/games:/snap/bin:/mnt/raid6/sunchuging/Softw
ares/miniconda3/bin:/mnt/raid1/sunchuqing/bc/bc/h/parallel-
meta/bin:/mnt/raid6/sunchuqing/Softwares/miniconda3/bin://mnt/raid1/
data/Software/prokka/bin:/mnt/raid6/sunchuqing/Softwares/miniconda3/
bin:/mnt/raid7/sunchuqing/Softwares/bin
bowtie2 -p 16 --un-conc 05 Removed/${infile} %.fastq --no-unal -k 20
-x /mnt/raid5/sunchuqing/Human Gut Phage/ref/hg38 ref -1
02 trimmed/${infile} clean.1.fq -2 02 trimmed/${infile} clean.2.fq
>log
bowtie2 -p 16 --al-conc 05 Removed/${infile} % prophage.fastq --end-
to-end -x /mnt/raid7/sunchuging/Human Gut Phage/db/HGYG prophage -1
05 Removed/${infile} 1.fastq -2 05 Removed/${infile} 2.fastq -S
../04 Abundance/${infile} NGS HGYG prophage.sam
bowtie2 -p 16 --un-conc 05 Removed/${infile} % phage.fastg --end-to-
end -x /mnt/raid7/sunchuqing/Human Gut Phage/db/HGYG -1
05 Removed/${infile} 1.fastq -2 05 Removed/${infile} 2.fastq -S
../04 Abundance/${infile} NGS HGYG.sam >log
cat 05 Removed/${infile} 1 prophage.fastq
05 Removed/${infile} 1 phage.fastq >
05 Removed/${infile} virome bft 1.fastq
cat 05 Removed/${infile} 2 prophage.fastq
05 Removed/${infile} 2 phage.fastq >
05 Removed/${infile} virome bft 2.fastq
R1=`ls 05 Removed/${infile} virome bft 1.fastq`
R2=`ls 05 Removed/${infile} virome bft 2.fastq`
java -jar $TRIMMO JAR FILE PE -threads 16 ${R1} ${R2}
05 Removed/${infile} virome 1.fastq
02 trimmed/${infile} clean unpaired.vir.1.fq
05 Removed/${infile}_virome_2.fastq
02 trimmed/${infile} clean unpaired.vir.2.fq
ILLUMINACLIP: $TRIMMO ADAPTOR FILE PE:2:15:10 LEADING:3 TRAILING:3
SLIDINGWINDOW:15:30 MINLEN:50
export
PYTHONPATH=/mnt/raid6/sunchuqing/Softwares/miniconda3/lib/python3.7/
site-packages:/mnt/raid8/sunchuqing/Softwares/lib/site-packages
# casper 02 trimmed/${infile} clean.1.fq
 02 trimmed/${infile} clean.2.fq
02 trimmed/${infile} clean.merged -t 16
pandaseq -f 02 trimmed/${infile} clean.1.fq -r
02 trimmed/${infile} clean.2.fq -F -w
02 trimmed/${infile} clean.merged.fastq -T 16
export PATH=$PATH:/mnt/raid6/sunchuqing/Softwares/miniconda3/bin
```

```
#16s
/mnt/raid6/sunchuqing/Softwares/ViromeQC/viromeqC/viromeQC.py \
    -i 02 trimmed/${infile} clean.merged.fastq \
    -o 04 Stats/${infile}.viromeqc \
    --bowtie2 threads 16\
    --diamond threads 16
#bac=`tail -n 1 04 Stats/${infile}.viromeqc | awk 'BEGIN {FS="\t"}
{print $4}'`
bacpct=`tail -n 1 04 Stats/${infile}.viromeqc | awk 'BEGIN {FS="\t"}
{print $4}'`
# /mnt/raid5/sunchuqing/Softwares/SortMeRNA/sortmerna-2.1/sortmerna
--ref /mnt/raid5/sunchuqing/Softwares/SortMeRNA/sortmerna-
2.1/rRNA databases/silva-bac-16s-
id90.fasta,/mnt/raid5/sunchuqing/Softwares/SortMeRNA/sortmerna-
2.1/index/silva-bac-16s-db --reads
02 trimmed/${infile} clean.merged.fastq --aligned 03 bac/${infile} -
-blast 1
#bowtie2 -x /mnt/raid6/sunchuqing/Database/Bacteria bowtie/bac -1
02 trimmed/${infile} clean.1.fq -2 02 trimmed/${infile} clean.2.fq
-S 03 bac/${infile}.sam
#bac=`/mnt/raid6/sunchuging/Softwares/miniconda3/bin/samtools
flagstat 03 bac/${infile}.sam | grep 'mapped' | awk 'BEGIN {FS=" "}
{print $1}' | head -n 1`
#cpn60
bowtie2 -x /mnt/raid5/sunchuqing/Human Gut Phage/ref/cpn60 ref -1
02 trimmed/${infile} clean.1.fq -2 02 trimmed/${infile} clean.2.fq
-S 03 bac cpn60/${infile}.sam
cpn=`/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools
flagstat 03 bac cpn60/${infile}.sam | grep 'mapped' | awk 'BEGIN
{FS=" "} {print $1}' | head -n 1`
reads=`/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools
flagstat 03 bac cpn60/${infile}.sam | grep 'total' | awk 'BEGIN
{FS=" "} {print $1}' | head -n 1`
#hq38
bowtie2 -x /mnt/raid5/sunchuqing/Human Gut Phage/ref/hg38 ref -1
02_trimmed/${infile}_clean.1.fq -2 02_trimmed/${infile} clean.2.fq
-S 03 human hg38/${infile}.sam
#/mnt/raid6/sunchuging/Softwares/miniconda3/bin/samtools flagstat
03 human hg38/${infile}.sam
human=`/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools
flagstat 03 human hg38/${infile}.sam | grep 'mapped' | awk 'BEGIN
{FS=" "} {print $1}' | head -n 1`
#bacpct=`awk 'BEGIN {printf "%.10f\n",
("'"${bac}"'"*100/"'"$reads"'")}'`
```

```
#bacpct=${bac}"%"
humanpct=`awk 'BEGIN {printf "%.10f\n",
("'"${human}"'"*100/"'"$reads"'")}'`
cpnpct=`awk 'BEGIN {printf "%.10f\n",
("'"${cpn}"'"*100/"'"$reads"'")}'`
echo
"File name, reads num, 16spct, cpn60 reads, cpn60pct, human reads, humanpc
t">04 Stats/${infile}.csv
echo
"${infile},${reads},${bacpct}%,${cpn},${cpnpct}%,${human},${humanpct}
}%" >>04 Stats/${infile}.csv
#IDBA
/mnt/raid5/sunchuqing/Softwares/idba/bin/fq2fa --merge --filter
05 Removed/${infile} virome 1.fastq
05 Removed/${infile} virome 2.fastq 05 Removed/${infile}.fa
/mnt/raid5/sunchuqing/Softwares/idba/bin/idba ud -r
05 Removed/${infile}.fa --maxk 120 --step 10 -o
O6 Assembly/${infile} --num threads 16 --min contig 1000
cat 06 Assembly/${infile}/contig.fa > 06 Assembly/${infile}.fasta
/mnt/raid6/sunchuqing/Softwares/cdhit-master/cd-hit-est -i
06 Assembly/${infile}.fasta -o 07 CD-HIT/${infile}.fa -c 0.95 -n 5
-T 15 -M 16000 >07 CD-HIT/${infile}.log
cd ..
```

### 2 # TGS assembly

```
#!/bin/bash
#SBATCH --cpus-per-task=16
                                # STDOUT
#SBATCH -o slurm.%N.%j.out
#SBATCH -e slurm.%N.%j.err
                                 # STDERR
infile=$1
NGS PATH=$2
G3 PATH=$3
Software='/mnt/raid6/sunchuging/Softwares'
LD LIBRARY PATH=$LD LIBRARY PATH:~/local/lib:/mnt/raid6/sunchuqing/S
oftwares/MCR/v94/runtime/qlnxa64:/mnt/raid6/sunchuqing/Softwares/MCR
/v94/sys/os/glnxa64:/mnt/raid6/sunchuging/Softwares/MCR/v94/extern/g
lnxa64:/mnt/raid3/wchen/miniconda2/pkgs/libgcc-7.2.0-h69d50b8 2/lib
cd G3
mkdir -p 01 ccs 02 Removed/${infile}
#Run CCS corrction
```

```
CCS=`ls ${G3 PATH}/*${infile}*.subreads.bam `
if [ ! -s 02 Removed/${infile}/${infile}.virome.fastq ];then
 if [ ! -s ${G3 PATH}/CCS/${infile}.subreads.bam ];then
    ${Software}/miniconda3/bin/ccs ${CCS} 01 ccs/${infile}.ccs.fastq
- j 16
 else
    CCS=`ls ${G3 PATH}/CCS/${infile}.subreads.bam`
    bedtools bamtofastq -i ${CCS} -fq 01 ccs/${infile}.ccs.fastq
  fi
 #Remove human genome
  bowtie2 -p 16 --un 02 Removed/${infile}/${infile}.fastq -x
/mnt/raid5/sunchuging/Human Gut Phage/ref/hg38 ref -U
01 ccs/${infile}.ccs.fastq > log
 bowtie2 --end-to-end -x
/mnt/raid7/sunchuqing/Human Gut Phage/db/HGYG_prophage -U
02 Removed/${infile}/${infile}.fastq -S
../04 Abundance/${infile} G3 HGYG prophage.sam -p 16 --al
02 Removed/${infile}/${infile}.prophage.fastq --quiet
  bowtie2 --end-to-end -x
/mnt/raid7/sunchuqing/Human Gut Phage/db/HGYG -U
02 Removed/${infile}/${infile}.fastq -S
../04 Abundance/${infile} G3 HGYG.sam -p 16 --un
02 Removed/${infile}/${infile}.phage.fastq --quiet
  cat 02 Removed/${infile}/${infile}.prophage.fastq
02 Removed/${infile}/${infile}.phage.fastq >
02 Removed/${infile}/${infile}.virome.fastq
  seqtk seq -a 02 Removed/${infile}/${infile}.virome.fastq |
/mnt/raid6/sunchuging/Softwares/miniconda3/bin/segkit rmdup -s -o
 02 Removed/${infile}/${infile}.fa
  rm 02 Removed/${infile}/${infile}.*hage.fastq
02 Removed/${infile}/${infile}.fastq
fi
#Assembly with Canu
rm -rf 03 canu Assembly/${infile}
mkdir -p 03 canu Assembly/${infile}
G3=`pwd`
cd 03 canu Assembly/${infile}
host=`hostname`
if [ "$host" = "bork" ];then
   /mnt/raid7/sunchuging/Softwares/OPERAMS-bork/canu/build/bin/canu
 /
      -p ${infile} \
      -d / genomeSize=20k corOutCoverage=1 \
      -corrected \
      -pacbio ../../02 Removed/${infile}/${infile}.fa \
      useGrid=false
```

```
else
    /mnt/raid6/sunchuqing/Softwares/canu/Linux-amd64/bin/canu \
      -p ${infile} \
      -d ./ genomeSize=20k corOutCoverage=1 \
      -corrected \
      -pacbio ../../02 Removed/${infile}/${infile}.fa \
      useGrid=false
fi
cd ${G3}
#Assembly with Flye
mkdir -p 03 Flye Assembly/${infile}
mkdir -p 05 Assembly/${infile}
zcat 03 canu Assembly/$infile/${infile}.trimmedReads.fasta.gz
>03 canu Assembly/$infile/${infile}.trimmedReads.fasta
flye --pacbio-corr 02 Removed/${infile}/${infile}.fa \
  --meta --genome-size 20k \
  --out-dir 03 Flye Assembly/${infile}/ \
  --threads 16 --min-overlap 1000
cat 03 canu Assembly/${infile}/${infile}.contigs.fasta
03 Flye Assembly/${infile}/assembly.fasta >
05 Assembly/${infile} fc.fa
#Binning with MetaBAT
if [ -s 03 canu Assembly/${infile}/${infile}.unitigs.fasta ];then
 mkdir -p 04 MetaBAT Assembly/${infile}
  cd 04 MetaBAT Assembly/${infile}
 mkdir -p db
 bowtie2-build
../../03 canu Assembly/${infile}/${infile}.unitigs.fasta
db/${infile}
 bowtie2 -x db/${infile} -U ../../01 ccs/${infile}.ccs.fastq -S
${infile}.sam
  /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools view -bS
${infile}.sam -o ${infile}.bam
  /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools sort
${infile}.bam > ${infile}.sort.bam
  ${Software}/berkeleylab-
metabat*/bin/jgi summarize bam contig depths --outputDepth
depth var.txt ${infile}.sort.bam
  ${Software}/berkeleylab-metabat*/bin/metabat -i
../../03 canu Assembly/${infile}/${infile}.unitigs.fasta -a
depth var.txt -o metabat -v
  for file in ./metabat.*.fa
    do
      num=${file//[!0-9]/}
```

```
#echo $num
      sed -e "/^>/ s/$/ ${num}/" metabat.$num.fa >>
metabat binned.concat.fasta
  done
  grep '>' metabat binned.concat.fasta | sed 's/>//g' >
metabat binned.info
  cd ../../05 Assembly/${infile}
  cut -f1 ../../03 canu Assembly/${infile}/${infile}.unitigs.bed |
sort|uniq -c > contig.list
  while read -r contig
  do
    num=`echo ${contig} | awk 'BEGIN {FS=" "} {print $1}'`
    tig=`echo ${contig} | awk 'BEGIN {FS="ctg"} {print $2}'`
    if [ $num -eq 1 ];then
      awk '/^>/ {printf("\n%s\t",$0);next;} {printf("%s",$0);} END
{printf("\n");}' <
../../03 canu Assembly/${infile}/${infile}.contigs.fasta |egrep -v
'^$'|tr "\t" "\n"
>../../03 canu Assembly/${infile}/${infile}.n1.contigs.fasta
      grep "$tig"
../../03 canu Assembly/${infile}/${infile}.n1.contigs.fasta -A 1|sed
's/>tig/>ctg/g' >> ../${infile}.fasta
      echo $tig >> infa.contig.list
      grep "ctg${tig}"
../../03 canu Assembly/${infile}/${infile}.unitigs.bed |cut -f4
>>infa.unitiq.list
    else
      grep "ctg${tig}"
../../03 canu Assembly/${infile}/${infile}.unitigs.bed |cut -f4
>unitig.list
      sed 's/utg/tig/g' unitig.list > uni.list
      binnum=`grep -Ff uni.list
../../04 MetaBAT Assembly/${infile}/metabat binned.info | awk 'BEGIN
{FS=" "} {print $2}' |sort |uniq |wc -l`
      inbin=`grep -Ff uni.list
../../04 MetaBAT Assembly/${infile}/metabat binned.info|wc -l`
      uninum=`cat unitiq.list | wc -l `
      #echo "$binnum,$uninum,${inbin}"
      if [[ $binnum == 1 && $uninum == $inbin ]];then
        grep "$tig"
../../03 canu Assembly/${infile}/${infile}.contigs.fasta -A 100| awk
-v RS='>' 'NR>1{i++}i==1{print ">"$0}' >> ../${infile}.fasta
        echo $tig >> infa.contig.list
        cat unitig.list >> infa.unitig.list
      fi
```

```
fi
```

```
done <"contig.list"</pre>
  cut -f4 ../../03 canu Assembly/${infile}/${infile}.unitigs.bed
  awk '/^>/ {printf("\n%s\t",$0);next;} {printf("%s",$0);} END
{printf("\n");}' <
../../03 canu Assembly/${infile}/${infile}.unitigs.fasta |egrep -v
'^$'|tr "\t" "\n"
>../../03 canu Assembly/${infile}/${infile}.n1.unitigs.fasta
  for unitig in `grep -v -Ff infa.unitig.list unitig.list `
    tig=`echo ${unitig} | awk 'BEGIN {FS="utg"} {print $2}'`
    #echo $tig
    grep "$tig"
../../03 canu Assembly/${infile}/${infile}.n1.unitigs.fasta -A 1 |
sed 's/class=contig/class=unitig/g'|sed 's/>tig/>utg/g' >>
../${infile}.fasta
  done
  cat ../../03 Flye Assembly/${infile}/assembly.fasta
../${infile}.fasta > ../${infile} fc.fa
  cd ../..
fi
mkdir -p 06 CD-HIT/${infile}
/mnt/raid6/sunchuqing/Softwares/cdhit-master/cd-hit-est -i
05 Assembly/${infile} fc.fa -o 06 CD-HIT/${infile}/${infile}.fa -c
0.95 -n 5 -T 16 -M 16000
Software='/mnt/raid6/sunchuqing/Softwares'
PATH="/mnt/raid6/sunchuging/Softwares/miniconda3/bin":$PATH
if [ -s ${NGS PATH}/*${infile}*R1* ];then
  if [ ! -s ../NGS/02 trimmed/*${infile}*clean.1.fq ];then
    TRIMMO JAR FILE='/mnt/raid1/tools/ngs tools/Trimmomatic-
0.38/trimmomatic-0.38.jar'
    TRIMMO ADAPTOR FILE PE='/mnt/raid1/tools/ngs tools/Trimmomatic-
0.38/adapters/TruSeg3-PE.fa'
    R1=`ls ${NGS PATH}/*${infile}*R1*`
    R2=`ls ${NGS PATH}/*${infile}*R2*`
    mkdir -p ../NGS/02 trimmed
    java -jar $TRIMMO JAR FILE PE -threads 4 $R1 $R2
../NGS/02 trimmed/${infile} clean.1.fg
../NGS/02 trimmed/${infile} clean unpaired.1.fq
../NGS/02 trimmed/${infile} clean.2.fq
../NGS/02 trimmed/${infile} clean unpaired.2.fq
ILLUMINACLIP: $TRIMMO ADAPTOR FILE PE:2:15:10 LEADING:3 TRAILING:3
SLIDINGWINDOW: 15:30 MINLEN: 50
  fi
```

```
mkdir -p 07 Pilon/${infile}
  cd 07 Pilon/${infile}
 mkdir -p index
 bwa index -p index/draft ../../06 CD-HIT/${infile}/${infile}.fa
 R1=`ls ../../NGS/02 trimmed/*${infile}*clean.1.fg`
 R2=`ls ../../NGS/02 trimmed/*${infile}*clean.2.fq`
 bwa mem -t 16 index/draft ${R1} ${R2} |
/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools sort -@ 10 -
O bam -o align.bam
 /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools index -@ 10
  /mnt/raid6/sunchuging/Softwares/miniconda3/bin/samtools sort -n
align.bam > align.sort.bam
 /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools fixmate -m
align.sort.bam fixmate.bam
  /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools sort -o
align.som.bam fixmate.bam
  /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools markdup
align.som.bam align markdup.bam
 /mnt/raid6/sunchuging/Softwares/miniconda3/bin/samtools view -@ 10
-q 30 -b align markdup.bam > align filter.bam
  /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/samtools index -@
10 align filter.bam
  java -Xmx5G -jar ${Software}/pilon-1.23.jar --genome ../../06 CD-
HIT/${infile}/${infile}.fa --frags align_filter.bam \
      --fix snps,indels \
      --output ${infile}.pilon
 cd ../..
fi
cd ..
```

### 3 Hybrid Assembly

```
#!/bin/bash
#SBATCH --cpus-per-task=16
#SBATCH -o slurm.%N.%j.out  # STDOUT
#SBATCH -e slurm.%N.%j.err  # STDERR

cd G2G3
infile=$1
NGS_PATH=
G3_PATH=
Software=
```

```
R1=`ls ../NGS/05 Removed/${infile} virome 1.fastq`
R2=`ls ../NGS/05 Removed/${infile} virome 2.fastg`
CCS="../G3/02 Removed/${infile}"
#rm -rf 01 OPERA MS assembly/${infile}
mkdir -p 01 OPERA MS assembly/${infile}
chmod 777 01 OPERA MS assembly/${infile}
perl OPERA-MS.pl \
    --short-read1 $R1 \
    --short-read2 $R2 \
    --long-read ${CCS}/${infile}.virome.fastq \
    --out-dir 01 OPERA MS assembly/${infile} \
    --contig-len-thr 1000 \
    --num-processors 64 \
    --polishing --no-strain-clustering --no-ref-clustering
#metaSPAdes
mkdir -p 02 metaSPAdes/${infile}
/mnt/raid6/sunchuqing/Softwares/SPAdes-3.13.1-
Linux/bin/metaspades.py \
    --pacbio ${CCS}/${infile}.virome.fastq \
    -1 $R1 \
    -2 $R2 \
    -o 02 metaSPAdes/${infile} \
    -t 16 \
    -m 750
mkdir -p 02 CD-HIT/${infile}
cat 02 metaSPAdes/${infile}/contigs.fasta
01 OPERA MS assembly/${infile}/contigs.fasta > 02 CD-
HIT/${infile}/${infile}.all.fa
/mnt/raid6/sunchuging/Softwares/cdhit-master/cd-hit-est -i 02 CD-
HIT/${infile}/${infile}.all.fa -o 02 CD-HIT/${infile}/${infile}.fa
-c 0.95 -n 8 -T 16 -M 16000
cd ..
```

## **Additional analysis**

#### 4 Viral annotation

#!/usr/bin/bash

```
infile=$1
export PATH=/mnt/raid8/sunchuqing/Softwares/blast2.2.26/ncbi-blast-
2.2.26+/bin:/mnt/raid7/sunchuqing/Softwares/bin:$PATH:/home/sunchuqi
ng/.local/bin:/mnt/raid6/sunchuging/Softwares/VirSorter/bi
n:/mnt/raid6/sunchuqing/Softwares/VirSorter/VirSorter:/usr/bin:/mnt/
raid6/sunchuqing/Softwares/miniconda3/envs/virsorter/bin:/mnt/raid6/
sunchuqing/Softwares/PPR-
Meta:/mnt/raid6/sunchuqing/Softwares/miniconda3/envs/virsorter/bin:/
usr/local/bin
export
LD LIBRARY PATH=$LD LIBRARY PATH:/mnt/raid6/sunchuqing/Softwares/MCR
/v94/bin/glnxa64:/mnt/raid6/sunchuging/Softwares/MCR/v94/runtime/gln
xa64:/mnt/raid6/sunchuging/Softwares/MCR/v94/sys/os/glnxa64:/mnt/rai
d6/sunchuqing/Softwares/MCR/v94/extern/glnxa64:/mnt/raid3/wchen/mini
conda2/pkgs/libgcc-7.2.0-h69d50b8 2/lib
echo "Dealing with $infile"
export
LD LIBRARY PATH=/mnt/raid8/sunchuqing/Softwares/lib/perllib:$LD LIBR
ARY PATH
export
PERL5LIB=/mnt/raid8/sunchuqing/Softwares/lib/perllib:$PERL5LIB
fafile="00 CAT/${infile} cdhit.len1.5k.fa"
seq="00_CAT/${infile}_cdhit.len1.5k.fa"
if [ ! -s $seq ];then
    awk '/^>/&&NR>1{print "";}{ printf "%s",/^>/ ? $0" ":$0 }'
 00 CAT/${infile} cdhit.fa |awk 'length($NF)>=1500 {print
1"\n" > 00 CAT/\{infile\} cdhit.len1.5k.fa
fi
n1=`ls 01 Glimmer/${infile}/*.predict|wc -l`
n2=`grep -c ">" 00 CAT/${infile} cdhit.len1.5k.fa `
if [ $n1 -ne $n2 ]; then
    #awk '/^>/&&NR>1{print "";}{ printf "%s",/^>/ ? $0" ":$0 }'
 00 CAT/${infile} cdhit.fa |awk 'length($NF)>=1500 {print
1''n''NF'' > 00 CAT/\{infile\} cdhit.len1.5k.fa
    #awk '/^>/{s=++num}{print >
"01 Glimmer/"file"/"file" "s"";close("01 Glimmer/"file"/"file" "s"")
}' file="${infile}" $fafile
    rm -rf 01 Glimmer/${infile}
    mkdir -p 01 Glimmer/${infile}
    #awk '/^>/{s=++num}{print > "01 Glimmer/"file"/"file" "s""}'
file="${infile}" $fafile
    #awk '/^>/{s=++num}{print >>
"01 Glimmer/"file"/"file" "s"";close("01 Glimmer/"file"/"file" "s"")
}' file="${infile}" $fafile
```

```
awk '/^>/{s=++num}{print > "01 Glimmer/"file"/"file" "s"";}
("01 Glimmer/"file"/"file" "s"")}' file="${infile}" $fafile
    cd 01 Glimmer/${infile}
    ls | grep -v "\." | grep "." |awk 'BEGIN {FS="."} {print
$1}'|sort|unig > ../${infile}.list
    while read -r line
    do
        /mnt/raid5/sunchuqing/Softwares/glimmer3.02/bin/long-orfs -
n -t 1.15 ${line} ${line}.longorfs
        /mnt/raid5/sunchuging/Softwares/glimmer3.02/bin/extract -t
${line} ${line}.longorfs > ${line}.train
        /mnt/raid5/sunchuqing/Softwares/glimmer3.02/bin/build-icm -
r ${line}.icm < ${line}.train
        /mnt/raid5/sunchuqing/Softwares/glimmer3.02/bin/glimmer3 -
o50 -q100 -t30 ${line} ${line}.icm ${line}
        /mnt/raid5/sunchuqing/Softwares/glimmer3.02/bin/extract -t
${line} ${line}.predict > ${line}_predict.fasta
    done < "../${infile}.list"</pre>
    cd ../..
fi
if [ ! -s "01 Glimmer/${infile}/${infile}.faa " ];then
    rm 01 Glimmer/${infile}/${infile}.faa
    while read -r line
    do
            declare -i len=0
            declare -i aorf=0
            declare -i orflen=0
            len=`grep '>' 01 Glimmer/${infile}/${line}.predict |awk
'BEGIN {FS=" "} {print $2}' |awk 'BEGIN {FS=" "} {print $2}'`
            orflen=`grep ">"
01_Glimmer/${infile}/${line}_predict.fasta | awk 'BEGIN {FS="="}
{print $2}'| awk '{sum+=$1}END{print sum}'`
            ((aorf=$orflen*10000))
            if [ $len -gt $aorf ]; then
                echo ${line} >>
 02 Annotate/${infile}/${infile}.no assignmnt under10k
                continue
            fi
            sed "s/^>/>${line} /"
01 Glimmer/${infile}/${line} predict.fasta >>
 01 Glimmer/${infile}/${infile}.faa
    done < "01 Glimmer/${infile}.list"</pre>
fi
#VirSorter
mkdir -p 04 is phage
cd 04_is_phage
```

```
if [ ! -s "04 Positive/${infile}.VirSorter.genome" ];then
    echo "--Start virSorter--"
   #rm -rf ./01 virsorter result/${infile}
    mkdir -p 01 virsorter result
 conda setup="$('/mnt/raid6/sunchuqing/Softwares/miniconda/bin/cond
a' 'shell.bash' 'hook' 2> /dev/null)"
    if [ $? -eq 0 ]; then
        eval "$ conda setup"
    else
        if [ -f
"/mnt/raid6/sunchuqing/Softwares/miniconda/etc/profile.d/conda.sh"
]; then
"/mnt/raid6/sunchuging/Softwares/miniconda/etc/profile.d/conda.sh"
        else
            export
PATH="/mnt/raid6/sunchuging/Softwares/miniconda/bin:$PATH"
        fi
    fi
    unset conda setup
   # # <<< conda initialize <<<
   # source activate virsorter
    conda activate
/mnt/raid6/sunchuqing/Softwares/miniconda3/envs/vs2
    virsorter run \
        -w ./01 virsorter result/${infile} \
        -i ../00 CAT/${infile} cdhit.len1.5k.fa \
        --db-dir /mnt/raid8/sunchuqing/Softwares/Virsorterdb\
        -j 16 --rm-tmpdir\
        --tmpdir ./01 virsorter result/${infile}/temp --min-score
0.7
    conda deactivate
   # mkdir -p 02 vir positive
./01 virsorter result/${infile}/Predicted viral sequences/VIRSorter
cat-1.fasta
./01 virsorter result/${infile}/Predicted viral sequences/VIRSorter
cat-2.fasta > 02 vir positive/${infile}.vir.fasta
    # #VirSprter Positive 基因组 02 vir positive/${infile}.vir.fasta
    mkdir -p 04 Positive
    #grep '>' 02 vir positive/${infile}.vir.fasta |sed
's/>VIRSorter //g'|sed 's/-cat 1//g'|sed 's/-cat 2//g'|sed 's/-
circular//g' > 04 Positive/${infile}.VirSorter.genome
    cat ./01 virsorter result/${infile}/final-viral-score.tsv | awk
'BEGIN {FS="|"} {print $1}'|sort|unig >
04 Positive/${infile}.VirSorter.genome
```

```
echo "--End virSorter--"
fi
#Complete
if [ ! -s "04 Positive/${infile}.cir.genome" ];then
    mkdir -p 03 circle/db 03 circle/${infile}
    /mnt/raid8/sunchuging/Softwares/blast2.2.26/ncbi-blast-
2.2.26+/bin/makeblastdb -in ../00 CAT/${infile} cdhit.len1.5k.fa -
out 03 circle/db/${infile} -dbtype nucl
    blastall -p blastn \
        -i ../00 CAT/${infile} cdhit.len1.5k.fa \
        -d 03 circle/db/${infile} \
        -o 03 circle/${infile}/${infile}.cir.tab \
        -m 8 -e 1e-5 -a 16
    awk '/^>/&&NR>1{print "";}{ printf "%s",/^>/ ? $0" ":$0 }'
../00 CAT/${infile} cdhit.len1.5k.fa |awk '{print
1'', "length(NF)}'|sed 's/>//g' > 03 circle/finfile.len
    awk 'BEGIN {FS="\t"} $1==$2 && $7==1 && $3==100.00 {print $0}'
03 circle/${infile}/${infile}.cir.tab >
03 circle/${infile}/${infile}.cir711.tab
    rm 03 circle/${infile}/${infile}.cir.len.tab
    while read -r line
    do
        contig=`echo $line |awk 'BEGIN {FS=","} {print $1}'`
        sed "s/^$contig\t/$line\t/g"
03 circle/${infile}/${infile}.cir711.tab |grep "$line" >>
03 circle/${infile}/${infile}.cir.len.tab
    done <"03 circle/${infile}/${infile}.len"</pre>
    awk 'BEGIN \{FS="[,\t]"\} ( \$2==\$10 \mid | \$2==\$11 ) && \$2!=\$9  \{print\}
$0}' 03 circle/${infile}/${infile}.cir.len.tab >
03 circle/${infile}/${infile}.circle.tab
    awk 'BEGIN {FS=","} {print $1}'
03 circle/${infile}/${infile}.circle.tab |sed 's/^/>/g' >
03 circle/${infile}/${infile}.complete
    awk '/^>/&&NR>1{print "";}{ printf "%s",/^>/ ? $0"\n":$0 }'
../00 CAT/${infile} cdhit.len1.5k.fa >
03_circle/${infile}/${infile}.fna
    #grep -w -A 1 -Ff 03 circle/${infile}/${infile}.complete
03 circle/${infile}/${infile}.fna |grep -v "-" >
03 circle/${infile}/${infile}.complete.fa
    #mkdir -p 04 complete
    #cp 03 circle/${infile}/${infile}.complete.fa 04 complete/
    sed 's/>//g' 03 circle/${infile}/${infile}.complete >
04 Positive/${infile}.cir.genome
#成环基因组 03_circle/${infile}/${infile}.complete
#02 vir positive/${infile}.vir.fasta
```

```
#pip install numpy
#pip install h5pv
#pip install tensorflow==1.4.1 #CPU version
#pip install keras==2.0.8
if [ ! -s "04 Positive/${infile}.ppr.genome" ];then
    pathh=`pwd`
   #PPR Meta
   # export PATH=/usr/bin:$PATH
   # >>> conda initialize >>>
   # !! Contents within this block are managed by 'conda init' !!
 conda setup="$('/mnt/raid6/sunchuqing/Softwares/miniconda/bin/cond
a' 'shell.bash' 'hook' 2> /dev/null)"
    if [ $? -eq 0 ]; then
        eval "$ conda setup"
    else
        if [ -f
"/mnt/raid6/sunchuging/Softwares/miniconda/etc/profile.d/conda.sh"
]; then
"/mnt/raid6/sunchuqing/Softwares/miniconda/etc/profile.d/conda.sh"
        else
            export
PATH="/mnt/raid6/sunchuqing/Softwares/miniconda/bin:$PATH"
        fi
    fi
    unset conda setup
   # <<< conda initialize <<<
    conda activate
/mnt/raid6/sunchuging/Softwares/miniconda3/envs/tensorflow
    pathh=`pwd`
    mkdir -p 03 PPR META
    cd /mnt/raid6/sunchuqing/Softwares/PPR-Meta
    ./PPR Meta $pathh/../00 CAT/${infile} cdhit.len1.5k.fa
$pathh/03 PPR META/${infile}.csv
    cd $pathh
    awk 'BEGIN {FS=","} $3>0.7 {print $1}'
 03 PPR META/${infile}.csv |awk 'BEGIN {FS=" "} {print $1}' >
04 Positive/${infile}.ppr.genome
   conda deactivate
fi
#VirFinder
#Bork
if [ ! -s "04_Positive/${infile}.virfiner.genome" ];then
   mkdir -p 03 VirFinder
```

```
/usr/bin/Rscript /mnt/raid5/sunchuqing/Buffalo gut/VirFinder.R
../00 CAT/${infile} cdhit.len1.5k.fa 03 VirFinder/${infile}.csv
   #VirFinder输出文件
    awk 'BEGIN {FS=","} $3>0.6 {print $1}'
03 VirFinder/${infile}.csv|awk 'BEGIN {FS=" "} {print $1}' >
04 Positive/${infile}.virfiner.genome
fi
#blast Virus ref
if [ ! -s "04 Positive/${infile}.ref.genome" ];then
    mkdir -p 03 Blast m8/${infile}
    blastall -p blastn \
        -i ../00 CAT/${infile} cdhit.len1.5k.fa \
        -d /mnt/raid6/sunchuging/Database/Virus/virus \
        -o 03 Blast m8/${infile}/${infile}.m8.tab \
        -m 8 -e 1e-10 -a 16
    grep -v -wFf /mnt/raid6/sunchuqing/Database/Virus/not.list
03 Blast m8/${infile}/${infile}.m8.tab >tmp && mv tmp
03 Blast m8/${infile}/${infile}.m8.tab
    #echo "chrom start end" > 03_Blast_m8/${infile}.bed
    awk 'BEGIN {FS="\t"} $3>50 {print $1 "\t" $7 "\t" $8 }'
03 Blast m8/${infile}/${infile}.m8.tab |sort -k 1 -t $'\t' |sed
's/[[:space:]]*$//'|tr -d " "|sort -k1,1 -k2,2n>
03 Blast m8/${infile}/${infile}.bed.1
    /mnt/raid6/sunchuqing/Softwares/miniconda3/bin/bedtools merge -i
03 Blast m8/${infile}/${infile}.bed.1
>03 Blast m8/${infile}/${infile}.bed
    sed 's/,/\t/g' 03 circle/${infile}/${infile}.len >
03 Blast m8/${infile}/${infile}.len
    cut -f 1 03 Blast m8/${infile}/${infile}.len >
03 Blast m8/${infile}/${infile}.list
    /mnt/raid6/sunchuging/Softwares/miniconda3/bin/bedtools
genomecov -i 03 Blast m8/${infile}/${infile}.bed -g
03 Blast m8/${infile}/${infile}.len |awk 'BEGIN {FS="\t"} $2>=1
{print $0}'|grep -v "genome" >
03 Blast m8/${infile}/${infile}.bedtools
    rm 03 Blast m8/${infile}/${infile}.genomecov.csv
   # while read -r line
   # do
          grep "$line" 03 Blast m8/${infile}/${infile}.bedtools |
awk 'BEGIN \{FS="\t"\} \{sum += \$NF\} \{print \$1 "," sum*100\}' | tail -n 1
>> 03 Blast m8/${infile}/${infile}.genomecov.csv
    #
          name=`grep "$line"
03 Blast m8/${infile}/${infile}.bedtools |wc -l`
          if [ $name -eq 0 ]; then
              echo "$line,0" >>
03_Blast_m8/${infile}/${infile}.genomecov.csv
```

```
# done <"03 Blast m8/${infile}/${infile}.list"</pre>
    awk 'BEGIN {FS=","} $NF>0.9 {print $1}'
03 Blast m8/${infile}/${infile}.genomecov.csv >
04 Positive/${infile}.ref.genome
    # awk 'BEGIN {FS=","} $2>90 {print $0}'
03 Blast m8/${infile}/${infile}.genomecov.csv >
03 Blast m8/${infile}/${infile}.genomecov.over90.csv
    # awk 'BEGIN {FS=","} {print $1}'
03 Blast m8/${infile}/${infile}.genomecov.over90.csv >
04 Positive/${infile}.ref.genome
#覆盖度超过90的基因组
03 Blast m8/${infile}/${infile}.genomecov.over90.csv
pathh=`pwd`
#blast pVOGs
if [ ! -s "04 Positive/${infile}.pvog.genome" ];then
    mkdir -p 03 Blast pVOG/${infile}
    rm 03 Blast pVOG/${infile}/${infile} cds num.csv
    cd ../01 Glimmer/${infile}
    ls | grep -v "\." | grep "." |awk 'BEGIN {FS="."} {print
$1}'|sort|uniq > ../${infile}.list
    cd $pathh
    blastall -p blastx -i ../01 Glimmer/${infile}/${infile}.faa -d
/mnt/raid6/sunchuqing/Database/Virus/blastdb/POGseqs \
            -o 03 Blast pVOG/${infile}/${infile}.m8.tab \
            -m 8 -e 1e-10 -a 16
   while read -r line
    do
        file="../01 Glimmer/${infile}/${line}"
        filename=`echo "$file" | awk 'BEGIN {FS="/"} {print $NF}'`
        CDSnum=`grep '>' ${file} predict.fasta |wc -l `
        name=`grep '>' ../01 Glimmer/${infile}/${line} |head -n
1|sed 's/>//q' `
        length=`grep -w "$name"
03 Blast m8/${infile}/${infile}.len`
        hitnum=`grep "${filename} "
03_Blast_pVOG/${infile}/${infile}.m8.tab |awk 'BEGIN {FS="\t"} $3>50
{print $1} ' |sort |uniq |wc -l`
        echo "$length,$CDSnum,$hitnum"|sed 's/\t/,/g'|awk 'BEGIN
{FS=","} $3>3 && $2/5000<$3 && $2/5000<$4 {print $0}' >>
03 Blast pVOG/${infile}/${infile} cds num.csv
    done <"../01 Glimmer/${infile}.list"</pre>
    awk 'BEGIN {FS=","} {print $1}'
 03 Blast pVOG/${infile}/${infile} cds num.csv >
04 Positive/${infile}.pvog.genome
fi
mkdir -p 05 Phage Positive
```

```
awk 'BEGIN {FS=","} {print $1}'
 03 Blast pVOG/${infile}/${infile} cds num.csv >
04 Positive/${infile}.pvog.genome
awk 'BEGIN {FS=","} {print $1}'
03 Blast m8/${infile}/${infile}.genomecov.over90.csv >
04 Positive/${infile}.ref.genome
sed 's/>//g' 03 circle/${infile}/${infile}.complete >
04 Positive/${infile}.cir.genome
awk 'BEGIN {FS=","} $3>0.6 {print $1}'
03 VirFinder/${infile}.csv|awk 'BEGIN {FS=" "} {print $1}' >
04 Positive/${infile}.virfiner.genome
awk 'BEGIN {FS=","} $3>0.7 {print $1}' 03 PPR META/${infile}.csv
|awk 'BEGIN {FS=" "} {print $1}' > 04 Positive/${infile}.ppr.genome
#grep '>' 02 vir positive/${infile}.vir.fasta |sed
's/>VIRSorter //g'|sed 's/-cat 1//g'|sed 's/-cat 2//g'|sed 's/-
circular//g' > 04 Positive/${infile}.VirSorter.genome
#cat 04 Positive/${infile}*|sort |uniq -c|sort -n|awk 'BEGIN {FS="
"} $1>=1 {print $2}'|sed 's/_length/ length/g'|awk 'BEGIN {FS=" "}
{print $1}' >05 Phage Positive/${infile}.phage.genome
cat 04 Positive/${infile}*|sort |uniq -c|sort -n|awk 'BEGIN {FS=" "}
$1>=2 {print $2}'|sed 's/_length/ length/g'|awk 'BEGIN {FS=" "}
{print $1}' >05 Phage Positive/${infile}.phage.genome2
cat 04 Positive/${infile}.cir.genome >>
05 Phage Positive/${infile}.phage.genome2
sed 's/ length/ length/g' 03 circle/${infile}/${infile}.fna >
03 circle/${infile}/${infile}.fna2
#grep -w -A 1 -Ff 05 Phage Positive/${infile}.phage.genome
03 circle/${infile}/${infile}.fna2 |grep -v -e "--" >
05 Phage Positive/${infile}.phage.genome.fa
grep -w -A 1 -Ff 05 Phage Positive/${infile}.phage.genome2
03 circle/${infile}/${infile}.fna2 | grep -v -e "--" | awk
'/^>/&&NR>1{print "";}{ printf "%s",/^>/ ? $0" ":$0 }' |awk '{print
$1","length($NF)}'|sed 's/>//g' |awk 'BEGIN {FS=","} $2>1500 {print
$1}'|sort|uniq > 05_Phage_Positive/${infile} fullfill2
grep -w -A 1 -Ff 05 Phage Positive/${infile}.phage.genome2
03 circle/${infile}/${infile}.fna2 | grep -v -e "--" | awk
'/^>/&&NR>1{print "";}{ printf "%s",/^>/ ? $0" ":$0 }' |awk '{print
$1","length($NF)}'|sed 's/>//g' |awk 'BEGIN {FS=","} $2>1500 {print
$0}'|sort|uniq > 05 Phage Positive/${infile} fullfill2.length.csv
grep -w -A 1 -Ff 05 Phage Positive/${infile} fullfill2
03 circle/${infile}/${infile}.fna2 |grep -v -e "--" >
05 Phage Positive/${infile}.phage.1.5k.genome.fa
```

### 5 UHGG filteration

```
#UHGG filter----
    seq=../CHGV.filtered.fa
    infile=CHGV
    blastall -p blastn \
        -i $seq \
        -d /mnt/raid8/sunchuqing/Database/UHGG/uhgg\
        -o ./${infile}.uhgg.m8.tab \
        -m 8 -e 1e-10 -a 20
    blastall -p blastn \
        -i $seq \
        - d
/mnt/raid7/sunchuqing/Human Gut Phage/HGYG/db/HGYG prophage \
        -o ./${infile}.uhgg.pro.m8.tab \
        -m 8 -e 1e-10 -a 20
    awk 'BEGIN {FS="\t"} $3>90 {print $1 "\t" $7 "\t" $8 }'
./${infile}.uhgg.m8.tab |sort -k 1 -t $'\t' |sed
's/[[:space:]]*$//'|tr -d " "|sort -k1,1 -k2,2n>
${infile}.uhgg.bed.1
    bedtools merge -i ${infile}.uhgg.bed.1 >${infile}.uhgg.bed
    seqtk comp $seq|cut -f 1,2 > ${infile}.uhgg.len
    bedtools genomecov -i ${infile}.uhgg.bed -g ${infile}.uhgg.len
|awk 'BEGIN {FS="\t"} $2>=1 {print $0}'|grep -v "genome" >
${infile}.uhgg.bedtools
    awk 'BEGIN {FS="\t"} $3>90 {print $1 "\t" $7 "\t" $8 }'
./${infile}.uhgg.pro.m8.tab |sort -k 1 -t $'\t' |sed
's/[[:space:]]*$//'|tr -d " "|sort -k1,1 -k2,2n>
${infile}.uhgg.pro.bed.1
    bedtools merge -i ${infile}.uhgg.pro.bed.1
>${infile}.uhgg.pro.bed
    bedtools genomecov -i ${infile}.uhgg.pro.bed -g
finfile.uhgg.len |awk 'BEGIN {FS="\t"} $2>=1 {print $0}'|grep -v
"genome" > ${infile}.uhgg.pro.bedtools
   while read -r lenf
    do
        name=`echo $lenf |awk 'BEGIN {FS=" "} {print $1}'`
        len=`echo $lenf |awk 'BEGIN {FS=" "} {print $2}'`
        line=`grep -w "$name" $infile.uhgg.bedtools`
        if [ `grep -w "$name" $infile.uhgg.bedtools|wc -l` == 1
];then
            pct=`echo $line |awk 'BEGIN {FS=" "} {print $NF}'|awk
```

```
'{printf("%f",$0)} '`
        fi
        pct2=0
        if [ `grep "$name" ${infile}.uhgg.pro.bedtools|wc -l` == 1
];then
            pct2=`grep "$name" ${infile}.uhgg.pro.bedtools |awk
'BEGIN {FS=" "} {print $NF}'|awk '{printf("%f",$0)} '`
        pct3=`echo "$pct-$pct2"|bc`
        if [ \$(echo "\$pct3 < 0"|bc) -eq 1 ]; then
            pct3=0
        fi
        echo "$name,$pct3,$len"
    done < "${infile}.uhgg.len" > ${infile}.np90.csv
    awk 'BEGIN {FS=","wc} $2>0.5' CHGV.np90.csv|wc
    awk 'BEGIN {FS=","} $2>0.5 {print $1}'
UHGG.filtered/CHGV.np90.csv > UHGG.filtered/CHGV.uhgg.txt
    segkit grep -v -f UHGG.filtered/CHGV.uhgg.txt CHGV.filtered.fa
> CHGV.filtered.uhgg.fa
```

### 6 checkv

```
infile=$1
export CHECKVDB=/mnt/raid6/sunchuqing/Softwares/checkv/checkv-db-
v0.6
export PATH=/mnt/raid6/sunchuqing/Softwares/miniconda3/bin:$PATH
awk '/^>/&&NR>1{print "";}{ printf "%s",/^>/ ? $0" ":$0 }'
 00 CAT/${infile} cdhit.fa |awk 'length($NF)>=1500 {print
1"\n" > 00 CAT/\{infile\} cdhit.len1.5k.fa
mkdir -p 05 checkV/${infile} 1.5k
/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/checkv contamination
00 CAT/${infile} cdhit.len1.5k.fa 05 checkV/${infile} 1.5k -t 16
/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/checkv completeness
00 CAT/${infile} cdhit.len1.5k.fa 05 checkV/${infile} 1.5k -t 16
/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/checkv
complete genomes 00 CAT/${infile} cdhit.len1.5k.fa
05 checkV/${infile} 1.5k
/mnt/raid6/sunchuqing/Softwares/miniconda3/bin/checkv
quality summary 00 CAT/${infile} cdhit.len1.5k.fa
05 checkV/${infile} 1.5k
```

### 7 RPKM calculation

```
#prevalence.RPKM
db=CHGV.filtered
seg=CHGV.filtered.fa
# bowtie2-build $seq db/$db
while read -r infile
do
    path="$NGS/"
    R1=`ls $path/$infile*R1*|head -n1`
    R2=`ls $path/$infile*R2*|head -n1`
    if [ -s 04 mapping/${infile} NGS.bam ];then
        echo "Mapped $infile"
        continue
    fi
    $Software/miniconda3/bin/bowtie2 -x db/$db \
        -1 $R1\
        -2 $R2 \
        -S 04 mapping/${infile} NGS.sam -p 8
    samtools view -bS 04 mapping/${infile} NGS.sam >
04 mapping/${infile} NGS.bam
    rm 04 mapping/${infile} NGS.sam
done < "sam.list"</pre>
while read -r infile
do
    if [ ! -s 04 mapping/${infile} NGS.bam ];then
        echo "Unmapped $infile"
        continue
    fi
    if [ ! -s 06 bamst cov/${infile}/chromosomes.report ];then
        echo "Bamdst $infile"
        mkdir -p 06 bamst cov/${infile}
        output=06 bamst cov/${infile}
        samtools sort -@4 04 mapping/${infile} NGS.bam -o
04 mapping/${infile} NGS.sort.bam
        cd $output
        $Software/bamdst/bamdst -p ../../CHGV.bed -o ./
../../04 mapping/${infile} NGS.sort.bam
        cut -f 1,6 chromosomes.report|sed "s/^/${infile}\t/"|awk
'BEGIN {FS="\t"} $3>=50' >> ../Sample.4x.50.cov.tbl
        cut -f 1,6 chromosomes.report|sed "s/^/${infile}\t/"|awk
'BEGIN \{FS="\t"\} $3>=50' \mid cut -f 2 > ./\t"\} $1
        grep -wFf ./$infile.4x.50.cov.list ../../CHGV.bed >
```

```
./$infile.4x.50.cov.bed
        cd ../..
    fi
    if [ ! -s 07 bam2rpkm/${infile}.rpkm.txt ];then
        echo "Bam2rpkm $infile"
            mkdir -p 07 bam2rpkm/
        rm -r 07 bam2rpkm/$infile.rpkm.txt
        bam=04 mapping/${infile} NGS.sort.bam
        samtools index -@4 $bam
        bed=06 bamst cov/${infile}/${infile}.4x.50.cov.bed
        echo $infile
        export total reads=$(samtools idxstats $bam|awk -F '\t'
'{s+=$3}END{print s}')
        echo The number of reads is $total_reads
        bedtools multicov -bams $bam -bed $bed |\
            perl -alne '{$len=$F[2]-$F[1];if($len <1 ){print</pre>
"$.\t$F[3]\t0" }else{$rpkm=(1000000000*$F[3]/($len*
$ENV{total reads}));print "$F[0]\t$F[3]\t$rpkm"}}' |\
            sed "s/^/${infile}\t/" >07_bam2rpkm/${infile}.rpkm.txt
    fi
done < "sam.list"</pre>
```