



FEB 02, 2024

Purification of recombinant Low Density Lipoprotein Receptor Related Protein Associated Protein 1 (LRPAP1, RAP) from Escherichia coli

Patricia Yuste-Checa¹, Silvia Gärtner¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

ABSTRACT

This protocol details how to efficiently purify LDL Receptor Related Protein Associated Protein 1 (LRPAP1 or RAP) from Escherichia coli.

ATTACHMENTS

Purification of
recombinant LRPAP1 from
E. coli.docx

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.rm7vzxbp2gx1/v1

Protocol Citation: Patricia Yuste-Checa, Silvia Gärtner, F Ulrich Hartl 2024. Purification of recombinant Low Density Lipoprotein Receptor Related Protein Associated Protein 1 (LRPAP1, RAP) from Escherichia coli. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzxbp2gx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 25, 2024

Last Modified: Feb 02, 2024

PROTOCOL integer ID: 94598

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's
Grant ID: ASAP-000282

MATERIALS

Buffers

▪ **Lysis buffer:**

A	B
Tris-Cl pH 8.0	50 mM
NaCl	300 mM
Imidazole	10 mM

▪ **High salt buffer:**

A	B
Tris-Cl pH 8.0	50 mM
NaCl	300 mM
Imidazole	250 mM

▪ **Low salt buffer:**


A	B
Tris-Cl pH 8.0	50 mM
NaCl	10 mM


- Size exclusion chromatography (SEC) buffer: 1x Phosphate buffered saline (PBS) pH 7.4

 pQTEV-LRPAP1 **addgene Catalog #31327**

LRPAP1 express4ion


18h 37m 45s


- 1 Thaw RbCl-competent Escherichia coli BI21 cells (DE3)  On ice

2 Add  1 μL of pQTEV-LRPAP1 plasmid without the signal peptide (1-35 amino acids) and incubate

30m



 00:30:00

 On ice



3 Heat shock  00:00:45 at  42 °C

45s


4 Incubate  On ice  00:02:00, then add  850 μL Lysogeny broth (LB) or Super Optimal broth with Catabolite repression (SOC) medium.

2m





5 Shake for  01:00:00 at  37 °C

1h

6 Centrifuge for  3000 x g, 00:05:00 and remove most of the supernatant.


5m




7 Resuspend the pellet with the remaining supernatant and plate the bacteria on LB /Ampicillin agar plates and incubate  Overnight at  37 °C


8h



8 Prepare preculture: Scrap all colonies with the scraper and inoculate 25-50 mL LB/Ampicillin. Shake at  37 °C for 4-6 h.

9 Measure OD₆₀₀ of the preculture and inoculate  6 L of LB media to an OD₆₀₀ = 0.05.

10 Shake flasks at  37 °C until approx. OD₆₀₀ = 0.5-0.8. (2-4 h).

11 Add isopropyl β-D-1-thiogalactopyranoside (IPTG) at final concentration of  1 millimolar (mM).



12 Shake flasks  Overnight at  22 °C.


8h



13 Centrifuge bacterial culture at  4000 rpm, 01:00:00. Discard supernatant.



1h






14 Resuspend each pellet with Lysis buffer (20 mL/1L bacteria) supplemented with Complete EDTA-free protease inhibitor cocktail (Merck). Flash-freeze in liquid nitrogen for storage at  -80 °C.

Lysis



1h 15m 50s

15 Thaw the cell pellets in a water bath at  22 °C and add lysis buffer (final volume 200 mL lysis buffer/ 6L bacteria) supplemented with Complete EDTA-free protease inhibitor cocktail (Merck) and Sm DNase  50 undetermined.

16 Add  1 undetermined lysozyme and incubate gently shaking for  00:30:00 at  4 °C.

30m



17 Sonicate lysate  On ice, 8 cycles  00:00:20 ON,  00:00:30 OFF. 50s

18 Centrifuge lysate at  40000 rpm, 4°C, 00:45:00 . 45m



Ni-NTA chromatography

19 Equilibrate the Ni-NTA column with 10 column volumes (CV, 20 mL) Lysis buffer.


20 Load lysate supernatant to the Ni-NTA column.

21 Wash the Ni-NTA column with 10 CV Lysis buffer.



22 Elute His7-TEV-RAP with 5 CV 100% High salt buffer and collect elution fractions.

23 Analyze eluted fraction by SDS-PAGE and Coomassie blue staining.



Ni-NTA chromatogram and SDS PAGE analysis. P: 20 μ L resuspended pellet + 20 μ L 2x SDS sample buffer, loaded 15 μ L. S: 20 μ L diluted supernatant + 20 μ L 2x SDS sample buffer, loaded 15 μ L. Fractions of interest: 10 μ L + 10 μ L 2x SDS sample buffer; loaded 6 μ L. Green box: Collected elution fractions (Fractions 29-41, 65mL).

Desalting

- 24** In order to reduce the salt concentration, load the eluted protein onto a HiPrep 26/10 desalting column equilibrated with the Low salt buffer.

His-TEV cleavage

8h

- 25** Collect eluted fraction containing protein and add glycerol at final concentration of 10%, DTT at final concentration of [M] 1 millimolar (mM), EDTA at final concentration of [M] 0.25 millimolar (mM) and His-TEV at final concentration of 93U per mg of protein.

- 26** Incubate at 4°C Overnight.

8h




Ni-NTA chromatography (Collect flow through)

- 27** Load the cleavage mixture onto a Ni-NTA column previously equilibrated with Lysis buffer and collect the flow through where the cleaved RAP protein should elute.

28 Wash the column with 5 CV (CV, 20 mL) of Lysis buffer and collect eluted fractions.



29 Analyze flow through fractions by SDS-PAGE and Coomassie blue staining.



Ni-NTA chromatogram and SDS PAGE analysis. TEV digest: 10 μ L sample + 10 μ L 2x SDS sample buffer, loaded 2.0 μ L. Fractions of interest: 10 μ L sample + 10 μ L 2x SDS sample buffer, loaded 5 μ L. Green boxes: Collected flow through (Fractions 20-21, 16 mL), and uncleaved His7-TEV-RAP protein eluted from the Ni-NTA column (Fractions 33-34).



Note

Some uncleaved His7-TEV-RAP protein may be eluted from the Ni-NTA column. Those fractions can be pooled, desalted and TEV digested again.

Size exclusion chromatography

30 Load RAP-containing fractions onto a Superdex-200 column previously equilibrated with SEC buffer.

31 Analyze eluted fractions by SDS-PAGE and Coomassie blue staining.



Size exclusion chromatogram and SDS PAGE analysis. 10 μ L fraction of interest + 10 μ L 2x SDS sample dye, loaded 2.0 μ L. Green box: Collected eluted fractions (72-77, 12mL).

32 Pool fractions containing RAP aliquot and flash-freeze in liquid nitrogen for storage at -80°C .

Note

The protein can be concentrated with a filter device like a VivaSpin20, MWCO 10,000 Da. **Approximate yield:** from 6 L of bacterial culture around 95 mg of pure RAP are obtained.