

Jun 04, 2025

## General western blot protocol

DOI

[dx.doi.org/10.17504/protocols.io.j8nlkomj6v5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkomj6v5r/v1)

Jailyn Izu<sup>1</sup>

<sup>1</sup>NIMH

TIU



Jailyn Izu

NIMH

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.j8nlkomj6v5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkomj6v5r/v1)

**Protocol Citation:** Jailyn Izu 2025. General western blot protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.j8nlkomj6v5r/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** October 23, 2023

**Last Modified:** June 04, 2025

**Protocol Integer ID:** 89749

## Abstract

General western blot protocol



## Materials

- 4-12% Bis-Tris (Invitrogen NuPAGE-SDS) mini gels (1.0 mm x 10 well or 1.5 mm x 15 well) (Invitrogen, #NP0321BOX /#NP0323BOX)
- MOPS-SDS running buffer (Invitrogen, #NP0001)
- 4x LDS buffer (Invitrogen, #NP0007)
- 10x dithiothreitol (DTT) (Thermo Fisher Scientific, #A39255)
- Li-Cor Chameleon Duo Pre-stained Protein Ladder (Li-Cor, P/N: 928-60000)
- Tris Glycine transfer buffer 10X (KD Medical, #RGF-3391)
- HPLC Methanol (Fisher Scientific, #A452SK)
- Li-Cor Intercept (TBS) blocking buffer (Li-Cor, P/N: 927-60001)
- Tris-buffered saline 10X (Corning, #46-012-CM)
- Tween20 (Sigma, #P1379)
- Revert 700 Total Protein Stain for Western Blot Normalization (Li-COR Biosciences, #926-11010)
- Primary antibodies
- Sodium azide (Sigma, #S-8032)
- IRDye 800CW Goat anti-Mouse IgG Secondary Antibody (Li-Cor, P/N: 926-32210), IRDye 680RD Goat anti-Rabbit IgG Secondary Antibody (Li-Cor, P/N: 926-68071), IRDye 800CW Donkey anti-Chicken Secondary Antibody (Li-Cor, P/N: 926-32218)
- Immun-Blot LF PVDF Membrane (BioRad, #1620264)
- Filter papers (Whatman, #1001125)
- XCell II Blot Module (Invitrogen, #EI9051)
- 1.5 mL microcentrifuge tubes
- Transfer sponges
- Aluminum foil
- Li-COR Odyssey M

## Protocol materials

☒ REVERT Total Protein Stain Kit **LI-COR Catalog #926-11010**


☒ Li-Cor Intercept (TBS) blocking buffer (Li-Cor, P/N: 927-60001) **LI-COR Catalog #P/N: 927-60001**






## Sample preparation

- 1 Prepare appropriate cellular extract for analysis (for example, see [Transfecting using FuGENE® 4K Transfection Reagent and cell lysis \(M-PER\) in a 6-well plate](#) protocol).
- 2 Measure the protein content of each sample (for example, see [Bradford Assay](#) protocol).




## Run gel

- 3 Add all to a  1.5 mL microcentrifuge tube:

For loading a 1.0 mm x **10 well**, 4-12% Bis-Tris (Invitrogen NuPAGE-SDS) mini gel (Invitrogen, #NP0321BOX):

-  5 µL 4x LDS buffer (final concentration : 1X)
-  2 µL DTT (final dilution is 1:10)
- Amount of total protein for each sample (calculate from Bradford Assay, 5-20 µg is typical)
- add H<sub>2</sub>O to  20 µL (final volume)

For loading a 1.5 mm x **15 well**, 4-12% Bis-Tris (Invitrogen NuPAGE-SDS) mini gel (Invitrogen, #NP0323BOX):

-  2.5 µL 4x LDS buffer (final concentration : 1X)
-  1 µL DTT (final dilution is 1:10)
- Amount of total protein for each sample (calculate from Bradford Assay, 5-20 µg is typical)
- add H<sub>2</sub>O to  10 µL (final volume)

- 4 Briefly spin the microcentrifuge tubes using a tabletop micro-centrifuge to collect the sample at the bottom of the tube and then load onto the gel. Add 5-10 µl of Li-Cor Chameleon Duo Pre-stained Protein Ladder to Lane 1 for relative molecular weight markers and sample orientation.

*Boil samples for 5 min. at 100°C if needed.*



5 Run gel at 100 volts in 1X MOPS-SDS running buffer for ⌚ 00:50:00

50m

6 Prepare transfer buffer solution at room temperature:

To make 2 L:

- 🧪 200 mL 10X Tris-glycine transblot buffer
- 🧪 400 mL 100% methanol
- ddH<sub>2</sub>O to 🧪 2 L

7 Pre-wet PVDF membrane for ⌚ 00:00:30 with 100% methanol, then soak in transfer buffer for about 1 hour.

30s

8 Soak transfer sponges and Whatman filter paper in transfer buffer for about 1 hour.

## Transfer

9 Assemble the transfer apparatus:

From bottom to top:

- sponge
- Whatman paper
- gel
- PVDF membrane
- Whatman paper
- sponge

10 Transfer proteins to the membrane at 40 volts for ⌚ 02:00:00

2h

## Total Protein Stain

11 After transferring, completely try the membrane by placing it on a piece of Whatman filter paper and incubating in 🌡 37 °C for ⌚ 00:10:00


10m

12 Stain for total protein by following the protocol for the  
🔗 REVERT Total Protein Stain Kit **LI-COR Catalog #926-11010**



- 13 Block the membranes for 1 hr using Li-Cor Intercept (TBS) blocking buffer

1h

 Li-Cor Intercept (TBS) blocking buffer (Li-Cor, P/N: 927-60001) LI-COR Catalog #P/N: 927-60001




(enough to cover the membrane).  01:00:00

## Primary antibodies

- 14 Decant the blocking buffer and add the appropriate primary antibody dilution.
- 15 Incubate the membrane and the primary antibody dilution at 4°C overnight with gentle shaking.

## Secondary antibodies

2h

- 16 Decant the primary antibodies. Wash the membrane with 1X TBS + 0.5% Tween20 + 0.3% sodium azide for 15 minutes with gentle shaking. Repeat the wash 2 more times.  00:45:00
- 17 Rinse the membrane briefly with ddH<sub>2</sub>O. Add secondary antibodies (1:10000) diluted in Licor blocking buffer solution + Tween20 (0.5%) and incubate for 1 hour covered in aluminum foil at room temperature with gentle agitation.  01:00:00
- 18 After the secondary antibody incubation, wash as outlined in Step 16. Add a fourth and final wash with ddH<sub>2</sub>O for 15 minutes.  01:00:00
- 19 Image the membrane using the Li-COR Odyssey M system.

45m

1h

1h