

FEB 08, 2024

# OPEN ACCESS



Protocol Citation: daniel.dautan, Per Svenningsson 2024. Open Field. protocols.io https://protocols.io/view/openfield-c8xuzxnw

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working We use this protocol and it's working

Created: Feb 08, 2024

Last Modified: Feb 08, 2024

PROTOCOL integer ID: 94932

Keywords: ASAPCRN, behavior,

motor

## Open Field

daniel.dautan<sup>1,2</sup>, Per Svenningsson<sup>1,2</sup>

<sup>1</sup>Department of Clinical Neuroscience, Karolinska Institutet, 171 76 Stockholm, Sweden; <sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA



Eileen Ruth Torres
Weill Cornell Medicine

#### **ABSTRACT**

Behavioral assay for general locomotion in mice.

#### **MATERIALS**

- Open field arenas: 40 x 40 dark box
- Recording camera
- Ethovision XT v11



#### **Funders Acknowledgement:**

Aligning Science Across Parkinson's Grant ID: 020608

1 Habituate naive mice to the room for 01:00:00.

2 Transfer each mouse to a custom 40 cm x 40 cm dark box, with a clear floor and a light suspended on the center of the apparatus.

3 Test duration was set at 00:20:00.

20m

4 After each test and between each group, clean apparatus thoroughly with 50% ethanol followed by water.

### 5 Analysis:

- Center of open field is defined by the space located at least 5 cm from any given wall.
- Time in the center, distance travelled, and average speed is obtained by the built-in functions of the tracking software Ethovision XT11.