



Aug 15, 2022

## Expression and purification of mCherry-NDP52

Daniel Bernklau<sup>1</sup>, Dorotea Fracchiolla<sup>1</sup>, Justyna Sawa-Makarska<sup>1</sup><sup>1</sup>Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.5qpvobdr9l4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvobdr9l4o/v1)

Dorotea Fracchiolla

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

### ABSTRACT

This protocol describes how to express in *E. Coli* human NDP52 N-terminally tagged with mCherry and purify it through His affinity purification followed by Size Exclusion Chromatography.

### DOI

[dx.doi.org/10.17504/protocols.io.5qpvobdr9l4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvobdr9l4o/v1)

### PROTOCOL CITATION

Daniel Bernklau, Dorotea Fracchiolla, Justyna Sawa-Makarska 2022. Expression and purification of mCherry-NDP52 . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5qpvobdr9l4o/v1>

### KEYWORDS

mCherry-NDP52, mCh-NDP52, NDP52, E. Coli, Expression , Purification, His Affinity , Size Exclusion Chromatography

### LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Jun 15, 2022

### LAST MODIFIED

Aug 15, 2022

### OWNERSHIP HISTORY

Jun 15, 2022 Justyna Sawa-Makarska

Jul 13, 2022  Dorotea Fracchiolla Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)

### PROTOCOL INTEGER ID

64640

## MATERIALS TEXT

### Protein information:

Molecular weight: 80384.47 Da

Ext. coefficient: 77240 M<sup>-1</sup>cm<sup>-1</sup>

Abs 0.1% (=1 g/l) 0.961, assuming all Cys residues are reduced

### Expression:

pETDuet-1\_6xHis-TEV-mCh-NDP52 (Addgene ID: 187829)

*E. coli* Rosetta pLysS cells

LB medium with antibiotics: 50 µl/ml ampicillin and 34 µl/ml chloramphenicol

IPTG (Isopropyl-β-D-thiogalactopyranosid)

### Lysis Buffer:

50 mM HEPES pH 7.5

300 mM NaCl

2 mM MgCl<sub>2</sub>

2 mM b-Met

Complete inhibitor EDTA free Roche

DNase

### Buffer A:

50mM HEPES pH 7.5

300 mM NaCl

10 mM Imidazole

2 mM Beta-Mercaptoethanol

### Buffer B:

50mM HEPES pH 7.5

300 mM NaCl

300 mM Imidazole

2 mM Beta-Mercaptoethanol

### SEC Buffer:

25 mM HEPES pH 7.5

150 mM NaCl

1 mM DTT

### Columns/Resin:

5-ml His-Trap column (Cytiva)

Superdex 200 increase 16/600 column (Cytiva)

## Expression

- 1 pETDuet-1\_6xHis-TEV-mCh-NDP52 (Addgene ID: 187829) was transformed into *E. coli* Rosetta pLysS cells.
- 2 To express the protein grow *E. coli* Rosetta pLysS cells in 4 L of LB medium (w Amp/Cam) at 37°C until an OD<sub>600 nm</sub> of 0.4. Next, bring the temperature down to 18°C and grow further to an OD<sub>600 nm</sub> of 0.8. Induce protein expression with 0.5 mM IPTG and grow for further 16 h at 18°C.
- 3 Pellet the cells at 4000 rpm ⚡ 4 °C ⌚ 00:15:00 . Re-suspended the cell pellet in a buffer containing 50mM HEPES<sup>15m</sup> pH7.5, 300 mM NaCl, 1 mM MgCl<sub>2</sub>, 10 mM Imidazole, 2 mM Beta-Mercaptoethanol, cOmplete protease inhibitors (Roche), and DNase (Sigma). Snap freeze in liquid nitrogen and store in ⚡ -80 °C until the day of purification.

## Purification

- 4 Open the cells by thawing in RT water bath and sonicating 3 x 30 seconds at 50% power using a Bandelin sonicator.

- 5 Clear the lysate by ultracentrifugation (40,000 rpm for 45 min at 4°C in a Ti45 rotor, Beckman).
- 6 Filter the supernatant with a 0.45 µm syringe filter on ice.
- 7 Apply to a 5-ml His-Trap column (Cytiva) and elute with a stepwise imidazole gradient (50, 75, 100, 150, 200, and 300 mM). Fractions at 75–100 mM imidazole should contain His-TEV-mCh-NDP52. Pool those fractions and subject to TEV protease cleavage over night at 4°C by very gentle rolling.
- 8 After TEV cleavage, concentrate the protein using a 50kDa cut-off Amicon filter to 0.5 ml and inject onto a Superdex 200 increase 16/600 column (Cytiva) pre-equilibrated with a buffer containing 25mM HEPES pH7.5, 150mM NaCl, 1mM DTT at 4°C.
- 9 Pool fractions containing pure protein, concentrate using a 50kDa cut-off Amicon filter, snap freeze in liquid nitrogen, and store at –80°C.