



Jun 22, 2020

# Purification of RNA from the Aqueous Phase Following TRIzol®/Chloroform Extraction using the Monarch® RNA Cleanup Kits

New England Biolabs<sup>1</sup><sup>1</sup>New England Biolabs

1 Works for me

This protocol is published without a DOI.

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com

 Isabel Gautreau  
New England Biolabs

## ABSTRACT

RNA isolation reagents containing guanidine thiocyanate and phenol (e.g., TRIzol, RNAzol®, QIAzol®, etc) combined with chloroform extraction, are often used for sample lysis and RNA purification. The aqueous phase from any guanidinium thiocyanate-phenol-chloroform extraction can be cleaned up using the Monarch RNA Cleanup Kits ([NEB #T2030](#), [T2040](#), [T2050](#)), thereby eliminating the need for tedious RNA precipitation step

## EXTERNAL LINK

<https://neb.com/protocols/2018/06/28/purification-of-rna-from-the-aqueous-phase-following-trizol-chloroform-extraction-using-the-monarch-rna-cleanup-kits>

## PROTOCOL CITATION

New England Biolabs 2020. Purification of RNA from the Aqueous Phase Following TRIzol®/Chloroform Extraction using the Monarch® RNA Cleanup Kits. **protocols.io**  
<https://protocols.io/view/purification-of-rna-from-the-aqueous-phase-followi-7rxhm7n>

## EXTERNAL LINK

<https://neb.com/protocols/2018/06/28/purification-of-rna-from-the-aqueous-phase-following-trizol-chloroform-extraction-using-the-monarch-rna-cleanup-kits>

## LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 27, 2019

## LAST MODIFIED

Jun 23, 2020

## OWNERSHIP HISTORY

Sep 27, 2019  Anita Broellochs protocols.ioJun 20, 2020  New England Biolabs Tech Support New England BiolabsJun 23, 2020  Isabel Gautreau New England Biolabs

## PROTOCOL INTEGER ID

28183

## SAFETY WARNINGS

Please see SDS (Safety Data Sheet) for hazards and safety warnings.

#### BEFORE STARTING

- Add 4 volumes of ethanol ( $\geq 95\%$ ) to the Monarch RNA Wash Buffer before use, as directed on the bottle.
- All centrifugation steps should be carried out at  $16,000 \times g$ . ( $\sim 13K$  RPM in a typical microcentrifuge). This ensures all traces of buffer are eluted at each step.

1 

Following guanidinium-thiocyanate-phenol-chloroform extraction, carefully transfer the upper aqueous phase into an RNase-free tube (not provided).

2 

Add 1 volume of ethanol ( $\geq 95\%$ ).

3 

Mix well by pipetting up and down or flicking the tube. **Do not vortex.**

4 

Insert an RNA cleanup column into a collection tube, load sample onto the column and close the cap.

5 Spin for  **00:01:00**, then discard flow-through.

For diluted samples  $\geq 900 \mu\text{l}$ , load a portion of the sample, spin, and then repeat as necessary.



To save time, spin for 30 seconds, instead of 1 minute.

6 

Re-insert the column into the collection tube. Add  **500  $\mu\text{l}$**  *RNA Cleanup Wash Buffer* and spin for  **00:01:00**.  
Discard the flow-through.




To save time, spin for 30 seconds, instead of 1 minute.

7 

Repeat wash (Step 6).

8 

Transfer the column to an RNase-free 1.5 ml microfuge tube (not provided). Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for 1 minute to ensure traces of salt and ethanol are not carried over to next step.

9 Elute in nuclease-free water according to the table below. The eluted RNA can be used immediately or stored at  **-70 °C**.

Care should be used to ensure the elution buffer is delivered onto the center of the matrix and not the wall of the column to maximize elution efficiency.


KIT	ELUTION VOLUME	INCUBATION TIME	SPIN TIME
T2030	6 – 20 $\mu$ l	N/A	1 minute
T2040	20 – 100 $\mu$ l	N/A	1 minute
T2050**	50 – 100 $\mu$ l	5 minutes (Room temp.)	1 minute

\* When cleaning up large amounts of RNA (> 100  $\mu$ g, NEB #T2050), some precipitation may occur following the addition of the Monarch RNA Cleanup Binding Buffer and ethanol to the sample (Steps 1 – 3). A pellet containing the RNA of interest may form on the side of the column following the first binding spin (Steps 4 and 5). To maximize recovery of this RNA, a second elution is recommended.

\*\* Yield may slightly increase if a larger volume is used, but the RNA will be less concentrated.

To save time, spin for 30 seconds, instead of 1 minute.

Recovery of RNA from Monarch RNA Cleanup Kits with Varying Elution Volumes



rRNA (10, 50 or 500  $\mu$ g, respectively of 16S and 23S Ribosomal Standard from *E. coli*, Sigma) was purified using a Monarch RNA Cleanup Kit (10  $\mu$ g, NEB #T2030) (50  $\mu$ g, NEB #T2040) (500  $\mu$ g, NEB #T2050). Nuclease-free water was used to elute the RNA. The percent recovery of the RNA was calculated from the resulting  $A_{260}$  as measured using a Trinean DropSense 16. ~80% of RNA can be efficiently recovered in 6  $\mu$ l from the Monarch RNA Cleanup Kit (10  $\mu$ g, NEB #T2030), 20  $\mu$ l from the Monarch RNA Cleanup Kit (50  $\mu$ g, NEB #T2040), and 50  $\mu$ l from the Monarch RNA Cleanup Kit (500  $\mu$ g, NEB #T2050).