



MAR 19, 2024

10 X Visium Spatial Gene Expression - Fixed Frozen Tissue Processing with CytAssist

Alberto Pappalardo¹, Kavya Batra¹, Rolando Perez-Lorenzo¹, Angela Christiano¹

¹Columbia University



kavya Batra

Columbia University Irving Medical Center

ABSTRACT

Here we summarize the recommendation released by 10X Genomics for processing fixed-frozen tissue sections while performing 10X Visium Spatial Gene Expression assay using the Visium CytAssist device. Please follow the link for the detailed official protocol on 10X Genomics website. https://cdn.10xgenomics.com/image/upload/v1680564483/support-documents/CG000662_Demonstrated_Protocol_VisiumCytAssist_FixedFrozen_H_E_RevA.pdf

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.q26g71z49gwz/v1

External link:

https://cdn.10xgenomics.com/image/upload/v1680564483/support-documents/CG000662_Demonstrated_Protocol_VisiumCytAssist_FixedFrozen_H_E_RevA.pdf

Protocol Citation: Alberto Pappalardo, Kavya Batra, Rolando Perez-Lorenzo, Angela Christiano 2024. 10 X Visium Spatial Gene Expression - Fixed Frozen Tissue Processing with CytAssist. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g71z49gwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 18, 2024

Last Modified: Mar 19, 2024

PROTOCOL integer ID: 96879

Funders Acknowledgement:
NCI
Grant ID: 5UG3CA275686

MATERIALS

Reagents

Visium FFPE Reagent Kit PN-1000436

Store at -20°C

A	B
Item	Part number
Amp Mix B	2000567
Extension Enzyme	2000389
Extension Buffer	2000409
RNase Enzyme	3000593
RNase Buffer B	2000551
Tissue Removal Enzyme	3000387
Tissue Removal Buffer B	2000543
Tissue Removal Buffer Enhancer	2000557
Decrosslinking Buffer	2000566
TS Primer Mix B	2000537
Block and Stain Buffer	2000554

Visium FFPE Reagent Kit PN-1000436

Other Specific Reagents

A	B	C	D
Item	Alternatives/Options	Vendor	Part Number
Ethanol	Ethyl Alcohol, 200 Proof	Millipore Sigma	E7023
	Ethanol absolute ≥99.5%	VWR	83813.360DP
Eosin	Eosin Y-solution, Alcoholic	Millipore Sigma	HT110116
Hematoxylin	Hematoxylin Solution, Mayer's	Millipore Sigma	MHS16
Bluing Reagent	Bluing Reagent, Dako	Agilent	CS70230-2
1X PBS	Phosphate-Buffered Saline, 1X without calcium and magnesium, PH 7.4	Corning	21-040-CV
Glycerol	Glycerol Solution	Millipore Sigma	49781
0.1 N HCl	Hydrochloric Acid Solution, 0.1 N	Fisher Chemical	SA54-1

Other Materials

Visium CytAssist Slide and Cassettes, 6.5 mm PN-1000519 for 2 runs
Store at ambient temperature

A	B	C
Item	Quantity	Part number
Visium Cassette, 8 port	1	3000811
Visium Tissue Slide Cassette		
Visium CytAssist moveable gasket small (pre-assembled with translator)	2	3000814
Visium CytAssist moveable translator (pre-assembled with gasket)	2	3000816
Visium CytAssist moveable Cassette, frame	2	3000813
Visium CytAssist Slide Seals, 40 pack	1	2000284
Visium CytAssist Spatial Gene Expression Slide v2, 6.5 mm	1	2000549

Visium CytAssist Slide and Cassettes, 6.5 mm
PN-1000519

Equipment

Thermal Cycler
Low Profile Plate Insert PN-3000823
10x Magnetic Separator PN-120250

Rehydration

1 Rehydration

- 1.1
- Place a Low Profile Thermocycler Adapter on a thermal cycler and preheat thermal cycler to 37°C.
- 1.2
- Retrieve the slides holding the tissue from the -80 °C freezer. The tissue sections should be 10 µm thick and placed within the allowable area for the CytAssist device. For Fisherbrand

Superforst slides the allowable area is a rectangle with the length sides distant 5 mm from the slide edge, and with the width sides distant 15 mm from the frosted area and the plus signs.

- 1.3 Place slide on the Low Profile Thermocycler Adapter with the tissue side facing up. Dry for 10 min at 37°C keeping the thermal cycler lid open.
- 1.4 Submerge the slide in a tube containing PBS for 5 min.
- 1.5 Submerge the slide slowly in a new tube containing Milli-Q Water for 3 min.
- 1.6 Submerge the slide slowly in a new tube containing 100% Ethanol for 3 min.
- 1.7 Submerge the slide slowly in a new tube containing 70% Ethanol for 3 min.
- 1.8 Submerge the slide slowly in a new tube containing Milli-Q Water for 20 sec.
- 1.9 Proceed immediately to H&E Staining & Coverslipping.

H&E Staining

2 H & E staining

- 2.1 Prepare 6 water beaker with Milli-Q Water to dip the slide.
- 2.2 Place slide on a flat impermeable staining surface.
- 2.3 Cover the slide with 1 ml of Hematoxylin making sure to uniformly cover the tissue section.
- 2.4 Incubate for 3 min at room temperature.
- 2.5 Remove the Hematoxylin and submerge the slide 5x in a water beaker (I).
- 2.6 Submerge the slide 15x in a new water beaker (II).
- 2.7 Submerge slide 15x in a new water beaker (III).

- 2.8 Place slide on back on the staining surface and add 1 ml of Bluing Buffer making sure to uniformly cover the tissue section. Incubate 1 min at room temperature.
- 2.9 Remove the Bluing Buffer and submerge the slide 10x in a new water beaker (IV).
- 2.10 Place slide on back on the staining surface and add 1 ml of Eosin making sure to uniformly cover the tissue section. Incubate 1 min at room temperature. DO NOT use diluted Eosin.
- 2.11 Remove the Eosin and submerge slide for 30 sec in a new water beaker (V).
- 2.12 Submerge the slide 10x in a new water beaker (VI)

Coverslipping

3 Coverslipping

- 3.1 Remove excess water and add 100 μ l mounting medium to uniformly cover the entire tissue section.
- 3.2 Place the coverslip without introducing bubbles and wait for the mounting medium to settle.

- 3.3** Once the coverslip settles, immediately proceed with imaging or store slide laying flat at 4°C in the dark for up to two weeks. If storing multiple slides avoid any contact between slides. DO NOT exceed two weeks of storage time.

Tissue Imaging

4 Imaging

- 4.1** Image each tissue section individually at the desired magnification using brightfield imaging settings (400x recommended).
- 4.2** After imaging, proceed immediately to Coverslip Removal or store as described above. DO NOT exceed two weeks of storage time.

5 Coverslip Removal

- 5.1** Fill a beaker with 800 ml of Milli-Q water (replace the water after processing 10 slides).
- 5.2** Submerge the slide holding it parallel to the surface and with the coverslipped surface facing sideways.

- 5.3 Hold the slide submerged until the coverslip slowly falls off from the slide
- 5.4 Gently dip the slide 30x in water to remove completely the mounting medium.
- 5.5 Let the slide air dry for a minimum of 5 min until the tissue is mostly dry. DO NOT exceed 20 min.
- 5.6 Incubate slide on the Low Profile Thermocycler Adapter with the thermal cycler lid open for 3 min at 37°C to dry the slide.
- 5.7 Proceed immediately to Decrosslinking.

Decrosslinking

6 Destaining

- 6.1 Place a Low Profile Thermocycler Adapter in the thermal cycler and set the following parameters.
- 6.2 Lid Temperature- 42°C
Reaction Volume - 100 µl 15

Run Time - 15 min

A	B	C
Step	Temperature	Time
Pre-equilibrate	42°C	Hold
Destaining	42°C	00:15:00
Hold	22°C	Hold

Incubation protocol for thermal cycler

- 6.3 Place the slide in the Visium CytAssist Tissue Slide Cassette.
- 6.4 Add 150 µl 0.1 N HCl in a 6.5mm gasket along the side of the wells to uniformly cover the tissue sections (DO NOT introduce bubbles) and ensure uniform coverage.
- 6.5 Remove HCl by careful aspiration.
- 6.6 Add 100 µl 0.1 N HCl along the side of the wells to uniformly cover the tissue sections (DO NOT introduce bubbles) and ensure uniform coverage.
- 6.7 Seal the cassette and place the slide on the Low Profile Thermocycler Adapter at 42°C.
- 6.8 Close the thermal cycler lid and initiate Destaining.

6.9 Remove the slide from the Low Profile Thermocycler Adapter and place on a flat surface. Some color leftover after the destaining is acceptable.

7 Decrosslinking

7.1 Prepare the thermal cycler with the following settings.

7.2 Lid Temperature- 70°C
Reaction Volume - 100 µl
Run Time - 30 min

A	B	C
Step	Temperature	Time
Pre-equilibrate	70°C	Hold
Decrosslinking	70°C	00:30:00
Cooling	22°C	00:10:00
Hold	22°C	Hold

Incubation protocol for thermal cycler

7.3 Remove all the HCl from the well corners.

7.4 Add 150 µl of diluted decrosslinking buffer.

- 7.5** Remove diluted decrosslinking buffer.
- 7.6** Add 100 μ l of diluted decrosslinking buffer.
- 7.7** Re-seal the cassette and place the slide on the Low Profile Thermocycler Adapter at 70°C.
- 7.8** Close the thermal cycler lid and initiate the decrosslinking.
- 7.9** Proceed immediately to Visium CytAssist Spatial Gene Expression User Guide (CG000495).