

VERSION 1  
JUL 27, 2023

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocols.io.8epv5jjq4l1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5jjq4l1b/v1)

**Protocol Citation:** Jesse Delia, Maiah Gaines-Richardson, Sarah C. Ludington, Najva Akbari, Cooper Vasek, Daniel Shaykevich, Lauren A O'Connell 2023. Tissue-specific in vivo transformation of plasmid DNA in Neotropical tadpoles using electroporation.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8epv5jjq4l1b/v1> Version created by Lauren A O'Connell

# Tissue-specific in vivo transformation of plasmid DNA in Neotropical tadpoles using electroporation V.1

Sarah C.

Jesse Delia<sup>1</sup>, Maiah Gaines-Richardson<sup>1</sup>, Ludington<sup>1</sup>,

Cooper

Daniel


Najva Akbari<sup>1</sup>, Vasek<sup>2</sup>,

Shaykevich<sup>1</sup>,

Lauren A O'Connell<sup>1</sup>

<sup>1</sup>Stanford University; <sup>2</sup>De Anza College

O'Connell Lab at Stanford



Lauren A O'Connell  
Stanford University

## ABSTRACT

Electroporation is an increasingly common technique used for exogenous gene expression in live animals, but protocols are largely limited to traditional laboratory organisms. The goal of this protocol is to enable in vivo electroporation techniques in a diverse array of tadpole species. We explore electroporation efficiency in tissue-specific cells of five species from across three families of tropical frogs—poison frogs (Dendrobatidae), forest frogs (Aromobatidae), and glassfrogs (Centrolenidae). These species are well-known for their diverse social behaviors and intriguing physiologies that coordinate chemical defenses, aposematism, and/or transparency. Specifically, we examine the effects of electrical pulse and injection parameters on species- and tissue-specific transfection of plasmid DNA in tadpoles. After electroporation of a plasmid encoding green fluorescent protein (GFP), we found strong GFP fluorescence within brain and muscle cells that increases with the amount of DNA injected and electrical pulse number. We discuss species-related challenges, troubleshooting, and outline ideas for improvement. Extending in vivo electroporation to diverse amphibian species will offer a powerful approach to explore topics in genetics, behavior, and organismal biology.

## ATTACHMENTS

[Delia et al electroporation how to.mp4](#)

## GUIDELINES

These parameters are a suggested starting point. Optimizing the protocol for species-specific applications to maximize transfection efficiency is recommended.

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** May 19, 2023

**Last Modified:** Jul 27, 2023

**PROTOCOL integer ID:**  
82172

**Keywords:** tadpole, brain, muscle, electroporation, plasmid

## MATERIALS

Ethyl 3-Aminobenzoate Methanesulfonate (Millipore Sigma Catalog #10521)  
Sodium Bicarbonate (Millipore Sigma Catalog #S6014)  
Josh's Frogs R/O Rx  
Platinum Foil (Fisher Scientific Catalog #AA11509FF)  
Lead Solder Wire (Amazon Catalog #B075WB98FJ)  
Loctite Fun-Tak Mounting Putty Tabs (Amazon Catalog #1865809-12)  
100 mm X 15 mm Petri Dishes (Fisher Scientific Catalog #FB0875713)  
Serological Pipette (Fisher Scientific Catalog #12-567-600)  
Electrical Tape (Fisher Scientific Catalog #19-047-280)  
Micromanipulator (Sutter Catalog #MM-33)  
Grass Instruments SD9 Square Pulse Stimulator  
Disposable Paired 13 mm Subdermal Needle Electrodes (MFIMedical Catalog #RHL-RLSND121-1-0)  
3.5" Replacement Glass Capillaries (Drummond Scientific Catalog #3-000-203-G/X)  
Sutter Instrument Co P-97 (Sutter Catalog #P-P7)  
Forceps (Fisher Scientific Catalog #12-000-157)  
Mineral Oil (Millipore Sigma Catalog #M8410)  
28 Gauge Metal Hub Blunt Point Needle (Fisher Scientific Catalog #14815616)  
1 mL BD Disposable Syringe (Fisher Scientific Catalog #14-823-30)  
Nanoject II Variable Volume Automatic Injector (Drummond Scientific Catalog #3-000-204)  
pCMV-GFP (Addgene Catalog #11153)  
Fast Green FCF (Millipore Sigma Catalog #F7258)  
Kimwipes (Fisher Scientific Catalog #06-666)  
Standard Disposable Transfer Pipettes (Fisher Scientific Catalog #13-711-7M)  
Micro Detail Paint Brush  
Stereomicroscope with a GFP filter

## SAFETY WARNINGS



MS-222 is a respiratory irritant and the following personal protective equipment should be worn: labcoat, gloves and safety glasses.

## BEFORE START INSTRUCTIONS

Consult with your local animal ethics board prior to experimentation.

## Anesthesia Preparation

5m

- 1 Mix 0.02g ethyl 3-aminobenzoate methanesulfonate (MS-222) and 0.08g sodium bicarbonate with 60 mL tadpole water 5m
- 2 Store at 4C for up to one week

## Electrode Set-Up for Targeting Muscle Fibers

25m

- 3 Remove the tips from two 5 mL serological pipettes using scissors 5m
- 4 Solder two ~ 5 mm X 8 mm pieces of platinum foil to separate electrical lead wires to make an electrode 5m
- 5 Run one electrode wire through each cut serological pipette and secure it with electrical tape 2m
- 6 Construct a platform out of clay evenly spread over the top of a Petri dish 2m
- 7 Embed the anode into the clay with the foil exposed near the center of the Petri dish 5m
- 8 Create a tadpole-sized (~ 1 cm) impression adjacent the anode foil 1m




Figure 1. An anode embedded in a flat clay platform to hold tadpoles during the electroporation procedure.

- 9 Mount the cathode on a micromanipulator and bend the foil parallel to the anode

2m




Figure 2. A cathode mounted on a micromanipulator in position to make contact directly above the anode during the electroporation procedure.

- 10 Set the stimulator parameters to square wave, 1 pps, 0.1 ms delay, 5 ms duration, and 30-50 V

1m

## Electrode Set-Up For Targeting Brain Cells

25m

- 11 Construct a platform out of clay molded into the shape of a hill

5m

- 12 Fix the clay on the top of a Petri dish

1m





Figure 3. A clay platform molded into the shape of a hill to hold tadpoles during the electroporation procedure.

- 13** Remove the tips from two 1 mL serological pipettes using scissors 5m
- 14** Mount the serological pipettes side by side on the micromanipulator and fasten them together with electrical tape 5m
- 15** Run one needle electrode wire through each cut serological pipette 2m
- 16** Position the electrode needles ~ 1 mm apart and parallel to one another 5m




Figure 4. An assembled electrode mounted on a micromanipulator.

- 17 Set the stimulator parameters to 10 pps, 0.1 ms delay, 15 ms duration, and 30 V

1m

#### Note

These settings are a suggested starting point. Optimizing these parameters for each species, as the best pulse shape and settings may vary, is recommended.

## Injection Set-Up

20m

- 18 Pull glass capillaries using a pipette puller

5m

- 19 Break the pipette tip at an angle using forceps to create a beveled tip

2m

- 20 Backfill the pipette with mineral oil using a 28-gauge needle and 1 mL syringe 2m
- 21 Place the pipette onto the injector plunger and tighten the collet 2m
- 22 Select an injection volume of 64.4 nl and set the injection rate to slow 2m
- 23 Empty enough mineral oil to load 2 uL of plasmid DNA solution 2m
- 24 Pipette 2 uL of plasmid DNA solution (0.25–0.27 $\mu$ g/ $\mu$ l) onto a piece of parafilm and mix with 0.2 uL 0.01% Fast Green 2m
- 25 Fill the pipette without introducing air bubbles 1m
- 26 Connect the electrode wires to the stimulator 1m
- 27 Place the platform under a dissection microscope with the electrode on one side and the injector on the other 1m





Figure 5. A dissection microscope focused on a tadpole injection platform that is placed between an electrode mounted on a micromanipulator and a microinjector mounted on a micromanipulator.

## Electroporation

25m

- 28 Anesthetize the tadpole by placing it in a Petri dish of room temperature 0.03% MS-222 for 5 minutes 5m
- 29 Confirm the tadpole is fully sedated by checking for movement in response to stimuli 1m
- 30 Move the tadpole to the platform with a cut transfer pipette 2m

- 31** Adjust the position of the tadpole using a paintbrush 2m
- 31.1** For experiments targeting muscle fibers, place the tadpole flat on its side with its head in the depression and its tail on top of the anode
- 31.2** For experiments targeting brain cells, place the tadpole dorsal side up with its head in the depression
- 32** Orient the platform such that the head of the tadpole is facing toward the injector 1m
- 33** Lower the injector and insert the pipette into the target tissue 2m
- 33.1** For experiments targeting muscle fibers, insert the pipette into a tail myomere
- 33.2** For experiments targeting brain cells, insert the pipette into the brain ventricle
- 34** Inject the plasmid DNA with a 5-10 s interval between each injection 2m

**34.1** For experiments targeting muscle fibers, deliver two injections

**34.2** For experiments targeting brain cells, deliver three injections

**35** Remove the pipette from the tadpole

1m

**36** Orient the platform such that the head of the tadpole is facing toward the electrode

1m

**37** Lower the electrode until it is in full contact with the target tissue

5m

**37.1** For experiments targeting muscle fibers, the tail should lay on top of the anode and the cathode should press on the tail directly above the anode



Figure 6. Muscle electroporation in the tail of a *Ranitomeya imitator* tadpole.

- 37.2** For experiments targeting brain cells, the electrode should be touching the head on either side of the brain




Figure 7. Brain electroporation in a *Ranitomeya imitator* tadpole with the electrode making contact with either side of the brain.

## 38 Deliver the electrical pulses

2m

### Note

The pulse range for targeting brain cells is a suggested starting point. Optimizing the protocol for species-specific applications to maximize transfection efficiency is recommended.

- 38.1** For experiments targeting muscle fibers, deliver 4-8 double pulses with a 1 s interval between each set of pulses



**38.2** For experiments targeting brain cells, deliver 4-10 pulses with a 1 s interval between each pulse

**39** Transfer the tadpole to fresh tadpole water for several hours to recover

4h

**40** Roughly 24 to 48 hours after electroporation, screen tadpoles for plasmid uptake by imaging GFP-positive cells

## In Vivo Screening

20m

**41** Anesthetize the tadpole by placing it in a Petri dish of room temperature 0.03% MS-222 for 5 minutes

5m

**42** Move the tadpole to a new Petri dish with tadpole water and place under a stereomicroscope with a GFP filter

5m

**43** Center the imaging field on the target tissue and capture the fluorescent image

10m

**44** Transfer the tadpole to fresh tadpole water for several hours to recover

4h