

APR 10, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.n92ldpq99l5b/v1

Protocol Citation: Katerina V Djambazova, Martin Dufresne, Angela R.S. Kruse, Jamie Allen, allison.b.esselman, Madeline E. Colley, David Anderson, ali.zahraei, Olof Isberg, Melissa Farrow, Jeff Spraggins 2023. High Spatial Resolution MALDI Imaging Mass Spectrometry Data Acquisition. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldpq99l5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Feb 08, 2023

Last Modified: Apr 10, 2023

PROTOCOL integer ID:
 76683


High Spatial Resolution MALDI Imaging Mass Spectrometry Data Acquisition

Katerina V Djambazova¹, Martin Dufresne¹,
 Angela R.S. Kruse¹, Jamie Allen¹, allison.b.esselman¹,
 Madeline E. Colley¹, David Anderson¹, ali.zahraei¹, Olof Isberg¹,
 Melissa Farrow¹, Jeff Spraggins¹

¹Vanderbilt University

VU Biomolecular Multimodal Imaging Center / Spraggins Research Group

Tech. support email: jeff.spraggins@vanderbilte.du



Katerina V Djambazova
 Vanderbilt University

ABSTRACT

This protocol provides the workflow and instrument parameters for setting up a MALDI imaging mass spectrometry experiment on a Bruker timsTOF Flex instrument. This protocol is intended for 10µm spatial resolution imaging of lipids (m/z 400 - 2000) in positive and negative ionization mode data acquisition in qTOF mode of operation.

MATERIALS

Agilent ESI-L Tune Mix (Agilent Technologies, Santa Clara, CA)

- 1 Place slides in a MTP 2-slide holder. Scan slides using flatbed scanner, allowing for sufficient contrast to visualize the tissue boundary. Ensure the slides have fiducials.

Note

Fiducials are markings on the slide that can be used as reference points for teaching the slides (Step 4).

- 2 Load the 2-slide holder in the instrument. Open timsControl software and select appropriate MALDI IMS data acquisition method. Generate target profile by ensuring "Use Target Profile" is enabled. Check the height at a few different locations on the slide and ensure the difference is <10 μm .

Note

Example Positive/Negative Ionization Mode Parameters for 10 μm MALDI IMS:

m/z 400 - 2000

Laser Power: 35 - 45%, Frequency: 10 000 Hz

Shots: 100 - 150 per pixel

Spot size: 10 μm with Beam scanning function disabled

Source temperature: ~50 °C

- 3 Calibrate the TOF with an external calibration using red phosphorus. Proceed if all calibration points are <1ppm standard deviation error and the calibration score is >99%.

Note

Red phosphorus is commonly used to perform calibration, however Agilent ESI-L Tune Mix can also be used for electrospray ionization.

List of Red Phosphorus Mass Calibration Peaks

A	B
Positive mode	Negative mode
278.7633	278.7644
402.6583	402.6594
464.6058	464.6069
526.5534	526.5545
650.4484	650.4495
774.3435	774.3445
898.2385	898.2396
960.1860	960.1870
1022.1335	1022.1346
1146.0286	1146.0297
1393.8187	1393.8198
1517.7137	1517.7148
1641.6088	1641.6099
1765.5038	1765.5049
1889.3989	1889.4000

- 4 Use the FlexImaging software to load the optical image of the slides. Teach the target position with three teaching points. Check the accuracy of the teaching by moving the stage to a fourth location.
- 5 Select desired tissue measurement region. Start data acquisition in the FlexImaging software.

